

**COOK®****MEDICAL****CE**  
0123**EN**  
3**Biliary Wire Guide (ES)****CS**  
3**Biliární vodící drát (ES)****DA**  
4**Galdevejskateterleder (elektrokirurgisk)****NL**  
5**Biliaire voerdraad (ES)****FR**  
6**Guide biliaire (ES)****DE**  
7**Elektrochirurgischer Gallengangsführungsdraht****EL**  
8**Συρμάτινος οδηγός χοληφόρων (ES)****HU**  
9**Epeúti vezetődrot (ES)****IT**  
9**Guide biliari (elettrochirurgiche)****NO**  
10**Galleledevaier (elektrokirurgisk)****PL**  
11**Prowadnik do dróg żółciowych  
(elektrochirurgiczny)****PT**  
12**Fio guia biliar (eletrocirúrgico)****ES**  
13**Guía biliar (electroquirúrgica)****SV**  
14**Biliär ledare (ES)**


\* 1 0 2 5 0 / 0 6 1 9 \*

1.

2.

3.

**INTENDED USE**

This device is used to assist in cannulation of the biliary and pancreatic ducts and to aid in bridging difficult strictures during ERCP.

**NOTES**

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

**CONTRAINDICATIONS**

Contraindications include those specific to ERCP.

**WARNINGS**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

**POTENTIAL ADVERSE EVENTS**

Potential adverse events associated with ERCP.

Potential adverse events associated with ductal cannulation and bridging of strictures include, but are not limited to: perforation, pancreatitis, bleeding or tissue inflammation.

**PRECAUTIONS**

Use of this wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide.

This wire guide may be left in place when used under normal conditions in conjunction with, and according to the instruction for, a compatible Cook electrosurgical (ES) device.

**INSTRUCTIONS FOR USE**

1. Prior to removing wire guide from holder, flush with 30 cc of sterile water. (See fig. 1)
2. Flush endoscope accessory channel and/or lumen of device with sterile water, then insert wire guide floppy end first. **Angled tip wire guides are not recommended for front loading through these adapters.** (See figs. 2 & 3) Note: For best results, wire guide should be kept wet, if applicable.
3. Fluoroscopically monitor wire guide advancement in ductal system.
4. Proceed with exchange of compatible wire-guided accessories over wire guide.

**Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.**

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**