EN

2

Soehendra® Stent Retriever

CS

3

Extraktor stentu Soehendra®

DA

4

Soehendra® stentudtager

NL

5

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
ENGLISH

INTENDED USE

This device is used to remove stents from the biliary and pancreatic ducts while maintaining wire guide placement.

NOTES

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

WARNINGS

Soehendra Stent Retriever is supplied sterile and intended for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks.

If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization. This device is not intended for dilation, as this could result in pancreatitis, perforation, or tissue inflammation.

PRECAUTIONS

Not compatible with coated wire guides. Use with coated wire guide could result in damage to wire guide coating.

Refer to package label for minimum channel size required for this device.

INSTRUCTIONS FOR USE

1. Introduce and advance threaded end of stent retriever over pre-positioned wire guide until it reaches stent. (See fig. 1)
2. Keeping end of stent retriever that is outside of scope straight, turn stent retriever handle clockwise until threaded end completely connects to stent. (See fig. 2) **Note:** Maximum thread engagement is recommended.
3. Once stent is securely connected to retriever, open elevator of endoscope and slowly pull stent and stent retriever out of channel making sure wire guide is left in place. **Note:** Dislodgement of wire guide from duct can be prevented by simultaneously pushing wire

guide and pulling stent and stent retriever out of endoscope. **Caution:** If stent detaches from stent retriever while in channel, push stent down scope and into elevator with stent retriever. Re-thread stent onto retriever, then slowly pull both out of endoscope.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.