COOK

MEDICAL

CE
0123

EN
2

Biliary Balloon Dilators

CS
4

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
ENGLISH

INTENDED USE

This device is used to dilate strictures of the biliary tree.

NOTES

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP and any procedures to be performed in conjunction with balloon dilation, including, but not limited to: asymptomatic strictures; inability to advance the balloon dilator through the strictured area; known or suspected perforation; severe inflammation or scarring near the dilation site.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP.

Potential adverse events associated with biliary ductal dilation include, but are not limited to: pancreatitis, cholangitis, perforation, and hemorrhage.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Do not inflate the balloon prior to introduction into the scope, as this may cause damage to the scope.

The endoscope must remain as straight as possible with the elevator in the open position (if applicable) when advancing or withdrawing the balloon dilator.

The entire balloon must be extended beyond the tip of the endoscope and be fluoroscopically visualized and positioned before inflation.

This balloon dilator is used in conjunction with a sterile inflation device and a 1:1 mixture of contrast and sterile water/saline. DO NOT use air or gaseous substances to inflate the balloon, as this will result in reduced balloon effectiveness.

During withdrawal of the balloon dilator from the endoscope, negative pressure is mandatory to maintain balloon deflation.

A compromised balloon may prohibit removal from the endoscope accessory channel. Removal of the endoscope along with the compromised balloon may be required.

Do not exceed the recommended balloon inflation pressure as listed on the catheter tag of the balloon dilator.

Do not advance balloon dilator if resistance is encountered. Assess the cause of resistance to determine if dilation should be re-attempted.

WARNINGS

This device is designed for single use only. Attempts to reprocess, re-sterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization

Do not use this device to perform sphincteroplasty, as this may lead to increased incidence of pancreatitis.

SYSTEM PREPARATION

1. Prepare the Inflation device per the manufacturer's instructions.
2. Attach the inflation device to the white inflation hub. (See fig. 1)
3. Create and maintain a vacuum with the inflation device.
4. Remove the protective sheath from the balloon. Apply a water soluble lubricant to the balloon to allow easier passage through the accessory channel.

INSTRUCTIONS FOR USE

1. Introduce the deflated balloon into the accessory channel of the endoscope, and then advance in short increments until the dilator is completely visualized endoscopically.
Note: The balloon may be introduced over a pre-positioned .035 inch or smaller wire guide. Prior to advancing the device over a pre-positioned wire guide, flush the wire guide lumen with sterile water or saline to facilitate advancement.
2. Following cannulation, contrast may be injected to fluoroscopically confirm the position of the device. If a wire guide is in place, remove the wire guide and inject contrast through the wire guide port. (See fig. 1)
3. Using fluoroscopic monitoring, continue advancement into the duct until the desired position is achieved. **Note:** The position may be verified using the metal bands located at either end of the balloon as fluoroscopic markers. If the wire guide was removed for contrast injection, flush the contrast from the catheter with sterile water or saline before reintroducing the wire guide, observing the necessary cautions.
4. Inflate balloon to specified pressure, as listed on the product's labeling. **Note:** The balloon dilator must be inflated with LIQUID ONLY.
5. Monitor the balloon fluoroscopically during the procedure.

6. To deflate the balloon, apply negative pressure to the balloon dilator. **Note:** The indicator on the pressure gauge should be in the vacuum portion during deflation.
7. Maintain negative pressure and observe the balloon for complete deflation. If the balloon does not fully deflate, create negative pressure again as indicated in step 6. Repeat this step as necessary until the balloon is completely deflated. Caution: The balloon must be thoroughly deflated, with all fluid removed, before withdrawal.
8. To remove the deflated balloon from the accessory channel, straighten the endoscope tip and open the elevator (if applicable). Maintain negative pressure while withdrawing the balloon using a continuous twisting motion. The previously placed wire guide may be left in position, in order to facilitate introduction of other wire guided devices.

Upon completion of this procedure, dispose of this device(s) per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.