EN 2	Bronchoscope Cytology Brush Instructions for Use
DA 4	Cytologibørste til bronkoskop Brugsanvisning
DE 7	Bronchoskop-Zytologiebürste Gebrauchsanweisung
EL 10	Ψήκτρα κυτταρολογίας βρογχοσκοπίου Οδηγίες χρήσης
ES 14	Cepillo de citología broncoscópica Instrucciones de uso
FI 17	Bronkoskoopin sytologiaharja Käyttöohjeet
FR 19	Brosse de cytologie pour bronchoscope Mode d'emploi
IT 22	Spazzolino citologico per broncoscopio Istruzioni per l'uso
NL 25	Cytologieborstel voor bronchoscoop Gebruiksaanwijzing
PT 28	Escova de citologia para broncoscópico Instruções de utilização
RO 31	Perie pentru citologie a bronhoscopului Instrucțiuni de utilizare
SV 34	Cytologiborste för bronkoskop Bruksanvisning

Bronchoscope Cytology Brush

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Performance Characteristics

The Bronchial Cytology Brush is comprised of a nylon brush attached to a stainless-steel hard wire and a polycarbonate handle. The device also has a PTFE sheath in which the brush can advance and retract. When the polycarbonate handle is manipulated, the brush is advanced and retracted. Refer to the product label for the device sheath length, sheath size, and minimum accessory channel.

Device Compatibility

Refer to package label for minimum channel size required for this device. The device can be injected with contrast.

Patient Population

This device is indicated for adult use only.

Intended User

This device is restricted to a trained healthcare professional.

Contact with Body Tissue

This device has contact with body tissue in line with intended use.

Operating Principle

The device is advanced through the endoscope to the desired sampling site. The brush of the device is advanced through the sheath by pushing forward on the handle. Once the brush is out of the sheath, the sample site is brushed several times to collect cellular material. After collecting the sample, the brush is pulled back into the sheath and the device is withdrawn from the endoscope.

INTENDED USE

This device is used to collect cells from the bronchi.

INDICATIONS FOR USE

This device is indicated for patients that require cell collection from the bronchial system such as for pulmonary/lung lesions/cancer or suspected pulmonary/lung disease (i.e., chronic obstructive pulmonary disease – COPD).

CLINICAL BENEFITS

Cook Bronchial Cytology Brushes are used to collect cells from the bronchi to provide samples for diagnostic testing or analysis.

WARNINGS

Do not use this device for any purpose other than stated intended use.

This device contains nickel which may cause an allergic reaction and/or systemic complications related to nickel allergy.

Sterile and Single Use

This device is supplied sterile and designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device contamination with biological or chemical agents and/or mechanical integrity failure of device.

CONTRAINDICATIONS

Those specific to the primary endoscopic procedure to be performed in gaining access to desired cytology site.

PRECAUTIONS

This device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL COMPLICATIONS

Those associated with pulmonary endoscopy include: airway injury that induces bronchospasm, epistaxis, laryngospasm, laryngeal edema, and/or laryngeal injury • aspiration • bacteremia • cardiac arrhythmia or arrest • death • fever • hemorrhage • hypotension • infection • nausea • perforation • pneumonia • pneumothorax • respiratory depression or arrest • respiratory failure • vomiting.

HOW SUPPLIED

This device is supplied ethylene oxide (EtO) sterile in a peel-open pouch.

STORAGE

Store in a dry location, away from temperature extremes.

INSPECTION OF DEVICE

Do not use the device if the sterile packaging is damaged or unintentionally opened before use. Visually inspect with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

INSTRUCTIONS FOR USE

1. Uncoil the catheter and prepare the device for use.
2. Remove the protective sheath with label from the brush (if included) and extend/retract brush.

3. Endoscopically inspect the area to be sampled.
4. With the cytology brush retracted, insert the device into the accessory channel of the endoscope.
5. Advance in short increments until the tip of the device is visualized exiting the endoscope.
6. Advance the device to the desired sampling site. Endoscopically and/or fluoroscopically monitor device advancement.
7. Advance the brush out of the sheath by pushing forward on the handle.
8. Brush the area to be sampled several times to obtain adequate cellular material.
9. Pull the brush back into the sheath and withdraw the device from the endoscope.
10. Prepare the specimen for examination per institutional guidelines.

DISPOSAL OF DEVICE

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, the competent authority, and/or regulatory authority of the country where the device was used.

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The symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class 1B). However, the device contains cobalt as a stainless steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

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If symbol appears on product label, X = minimum accessory channel

DA Hvis symbolet vises på produktets mærkning, er X = minimum for tilbehørskanal • **DE** Sofern das Symbol auf dem Produktetikett aufgeführt ist: X = Minstdurchmesser des Arbeitskanals • **EL** Εάν εμφανίζεται το σύμβολο στην επισήμανση του προϊόντος, X = ελάχιστο κανάλι εργασίας • **ES** Si el símbolo aparece en la etiqueta del producto, X = diámetro mínimo del canal de accesorios • **FI** Jos symboli esiintyy tuotemerkinnössä, X = pienin työskentelykanava • **FR** Si le symbole est visible sur l'étiquette du produit, X = canal opérateur minimal • **IT** Se questo simbolo compare sull'etichetta del prodotto, X = canale operativo minimo • **NL** Als dit symbool op het productetiket staat: X = minimale diameter van werkkanaal • **PT** Se o símbolo aparecer no rótulo do produto, X = canal acessório mínimo • **RO** Dacă simbolul apare pe eticheta produsului, X = canalul pentru accesorii minim • **SV** Om symbolen finns på produktens etikett, X = minsta arbetskanal

If symbol appears on product label, X = maximum insertion portion width

DA Hvis symbolet vises på produktets mærkning, er X = maximum for tilbehørskanal • **DE** Sofern das Symbol auf dem Produktetikett aufgeführt ist: X = maximale Breite des einzuführenden Teils • **EL** Εάν εμφανίζεται το σύμβολο στην επισήμανση του προϊόντος, X = μέγιστο πλάτος τμήματος εισαγωγής • **ES** Si el símbolo aparece en la etiqueta del producto, X = anchura máxima de la parte que se va a insertar • **FI** Jos symboli esiintyy tuotemerkinnössä, X = sisäänvientiosan suurin leveys • **FR** Si le symbole est visible sur l'étiquette du produit, X = largeur maximale de la section d'insertion • **IT** Se questo simbolo compare sull'etichetta del prodotto, X = larghezza massima della sezione inseribile • **NL** Als dit symbool op het productetiket staat: X = maximale breedte van inbrenggedeelte • **PT** Se o símbolo aparecer no rótulo do produto, X = largura máxima da parte introduzida • **RO** Dacă simbolul apare pe eticheta produsului, X = lățimea maximă a porțiunii de introducere • **SV** Om symbolen finns på produktens etikett, X = maximal bredd på införsdelen

EU Importer

DA EU-importør • **DE** EU-Importeur • **EL** Εισαγωγέας στην ΕΕ • **ES** Importador en la UE • **FI** EU-maahantuoja • **FR** Importateur au sein de l'UE • **IT** Importatore UE • **NL** Importeur in de EU • **PT** Importador para a UE • **RO** Importator în UE • **SV** EU-importör

Units per package

DA Enheder pr. pakke • **DE** Einheiten pro Packung • **EL** Μονάδες ανά συσκευασία • **ES** Unidades por envase • **FI** Yksikköä pakkausta kohden • **FR** Unités par emballage • **IT** Unità per confezione • **NL** Eenheden per verpakking • **PT** Unidades por embalagem • **RO** Unități per ambalaj • **SV** Enheter per förpackning

Australian Sponsor

Australian Sponsor

DA Australsk sponsor • **DE** Australischer Sponsor • **EL** Χορηγός για την Αυστραλία • **ES** Patrocinador australiano • **FI** Australialainen rahoittaja • **FR** Partenaire australien • **IT** Sponsor australiano • **NL** Australische opdrachtgever • **PT** Promotor australiano • **RO** Sponsor australian • **SV** Australisk sponsor

Rx only

This device is intended
for single use only.

cookmedical.com/symbol-glossary

Australian Sponsor

William A. Cook Australia Pty. Ltd.

95 Brandl Street,
Eight Mile Plains QLD 4113
Australia

Wilson-Cook Medical, Inc.

4900 Bethania Station Road
Winston-Salem, North Carolina
27105
USA

Cook Medical Europe LTD

O'Halloran Road
National Technology Park
Limerick
Ireland

2025-09
10512_0925
© COOK 2025
cookmedical.com