
EN 2	Non-lithotripsy Extraction Basket
CS 4	Extrakční košík, není kompatibilní s litotrypsí
DA 6	Udtrækningskurv, ikke til litotripsi
NL 8	Extractiebasket, niet compatibel met lithotripsie
FR 11	Panier d'extraction non compatible avec la lithotripsie
DE 13	Extraktionskörbchen (nicht für Lithotripsien)
EL 16	Καλάθι εξαγωγής λίθων χωρίς λιθοτριψία
HU 18	Nem lithotripsiás extraháló kosár
IT 20	Cestello estrattore non per litotripsia
NO 23	Kurv for uttrekking uten litotripsi
PL 25	Koszyk do ekstrakcji bez litotrypsji
PT 27	Cesto de extração não compatível com litotripsia
ES 30	Cesta de extracción no compatible con litotricia
SV 32	Korg för extraktion utan litotripsi

- 1 Connector Cap
 Víčko spojky
 Konnektorhætte
 Connectorodop
 Capuchon de connecteur
 Anschlusskappe
 Πώμα συνδέσμου
 Csatlakozósapka
 Cappuccio del connettore
 Tilkoblingshette
 Nasadka łącznika
 Tampa do conector
 Tapa del conector
 Kopplingsskydd

ENGLISH

INTENDED USE

This device is used for the endoscopic removal of stones in the biliary system and foreign bodies.

NOTES

Do not use this device for any purpose other than the stated intended use.

Store in a dry location. Recommended long term storage temperature is 22°C (71.6°F). Device may be exposed to temperatures up to 50°C (122°F) during transportation.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in conjunction with the removal of the foreign body.

Perforation.

If this device is to be used for biliary stone extraction, contraindications include an ampullary opening inadequate to allow for unimpeded passage of the stone and basket.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy and ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction, respiratory depression or arrest, cardiac arrhythmia or arrest.

Others potential complications associated with biliary stone or foreign body removal include, but are not limited to: impaction of the object, aspiration of the foreign body, localized inflammation, pressure necrosis, perforation and ductal trauma.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

If the device is compatible with a wire guide, the wire guide diameter and the inner lumen of the wire guided device must be compatible.

If this device is to be used for the removal of biliary stones, assessment of the stone size and ampullary orifice must be made to determine the necessity of sphincterotomy. In the event sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

WARNINGS

This device is **not** compatible with the Soehendra® Lithotripter or any other mechanical lithotripter. Use of this device with a mechanical lithotripter may result in device failure and/or basket impaction. Surgical intervention may be required if impaction occurs.

When retrieving a foreign body from the upper GI tract, extreme caution must be exercised to avoid occlusion of the airway.

Baskets should not be used in the pancreatic duct. Use of baskets in the pancreatic duct may result in pancreatitis and/or tearing, rupturing or laceration of the ductal system.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

INSTRUCTIONS FOR USE

Note, if applicable: Flush the wire guide port with sterile water to facilitate easy passage of a wire guide.

1. With the basket fully retracted into sheath, insert the device into the endoscope accessory channel. Advance the device through the channel, in short increments, until the basket sheath exits the endoscope.

2. Endoscopically or fluoroscopically visualize the extraction target and advance the basket to the appropriate position. If applicable, fluoroscopic visualization may be enhanced by injecting contrast via the side arm fitting. **Note, if applicable:** Tighten the connector cap on the handle before injection of contrast. (See fig. 1)
3. Place the sheath of the basket beyond the target. Confirm the desired position of the basket sheath relative to the target.
4. Advance the basket out of the sheath by pushing forward on the handle. **Caution:** Pulling on the sheath while advancing or retracting the basket may damage the device, rendering it inoperable.
5. Place the basket around the extraction target or foreign body and slowly retract the handle until entrapment is achieved. **Note:** It may be necessary to move the basket back and forth along the target in order to capture the target.

For stone removal: Withdraw the basket into the channel until the stone is captured. Pull back on the sheath to remove the basket and stone from the duct. For multiple stones, repeat the process until all stones have been removed. **Caution:** If withdrawal of the basket from the biliary duct is restricted, surgical intervention may be necessary.

For foreign body removal: Withdraw the basket into the channel while maintaining entrapment. Support the foreign body against the endoscope tip and withdraw the endoscope from the patient. Maintain endoscopic visualization to ensure continued entrapment of the object.

6. Prepare or dispose of the extracted specimen per institutional guidelines.

Upon completion of the procedure, dispose of the device(s) per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.