en 6	Fusion® Sphincterotome Instructions for Use
bg 9	Сфинктеротом Fusion® Инструкции за употреба
cs 14	Sfinkterotom Fusion® Návod k použití
da 18	Fusion® sfinkterotom Brugsanvisning
de 22	Fusion® Sphinkterotom Gebrauchsanweisung
el 27	Σφιγκτηροτόμος Fusion® Οδηγίες χρήσης
es 31	Esfinterótomo Fusion® Instrucciones de uso
et 36	Fusion® Sfinkterotoom Kasutusjuhend
fi 40	Fusion®-sfinkterotoomi Käyttööhjeet
fr 44	Sphinctérotome Fusion® Mode d'emploi
hr 48	Sfinkterotom Fusion® Upute za uporabu
hu 52	Fusion® sphinkterotóm Használati utasítás
it 56	Sfinterotomo Fusion® Istruzioni per l'uso
lt 60	„Fusion®“ sfinkterotomas Naudojimo instrukcija
lv 64	Fusion® sfinkterotoms Lietošanas instrukcija
nl 68	Fusion®-sfincterotoom Gebruiksaanwijzing
no 72	Fusion®-sfinkterotom Bruksanvisning
pl 76	Sfinkterotom Fusion® Instrukcja użycia

pt
81

Esfincterótomo Fusion®

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ro
85

Sfincterotom Fusion®

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sk
89

Sfinkterotóm Fusion®

Návod na použitie

sl
93

Sfinkterotom Fusion®

Navodila za uporabo

sv
97


Fusion® sfinkterotom

Bruksanvisning

tr
101

Fusion® Sfinkterotom

Kullanma Talimatı

1.

1 Active Cord Connector (Metal Hot Pin)

bg Активен кабелен конектор (метален гъвкав щифт) • **cs** Konektor aktivního vodiče (kovový aktivní pin) • **da** Stik med strømførende ledning (aktivt metalstikben) • **de** HF-Kabel-Anschluss (elektrisch aktiver Metallstift) • **el** Σύνδεσμος ενεργού καλωδίου (μεταλλική θερμή καρφίδα) • **es** Conector del cable de conexión (pin activo) • **et** Aktiivjuhtme konnektor (metallist aktiivtihvt) • **fi** Aktivointijohdon liitin (metallinen jännitteinen piikki) • **fr** Connecteur de cordon actif (broche métallique sous tension) • **hr** Priključak za aktivni kabel (metalna igla pod naponom) • **hu** Aktiv vezeték csatlakozója (fém aktiv tű) • **it** Connettore del cavo diatermico (pin metallico attivo) • **lt** Aktyviojo laidų jungtis (metalinis elektrinis kištukas) • **lv** Aktivā vada savienotājs (metāla karstā tapa) • **nl** Connector actief snoer (metalen hot pin) • **no** Kontakt for strømførende ledning (strømførende metallpinne) • **pl** Złącze przewodu zasilania (metalowy styk aktywny) • **pt** Conector do cabo ativo (pino metálico aquecido) • **ro** Conector cablu activ (pin metalic încălzit) • **sk** Konektor aktívnej šnúry (kovový aktívny kolík) • **sl** Priključek aktivnega kabla (kovinski aktivni priključek) • **sv** Kontakt till strömförande kabel (kontaktstift av metall) • **tr** Aktif Kablo Konektörü (Aktif Metal Pim)

2 Wire Stop

bg Стопер за теления водач • **cs** Zarážka drátu • **da** Kateterlederstop • **de** Drahtstopper • **el** Стоп σύρματος • **es** Tope de la guía • **et** Traadi piirik • **fi** Johdinsalpa • **fr** Butée du guide • **hr** Graničnik žice • **hu** Dróthatároló • **it** Meccanismo di arresto della guida • **lt** Vielos stabdiklis • **lv** Stīgas aizturis • **nl** Draadstop • **no** Væierstopp • **pl** Przrząd hamujący ruch przewodnika • **pt** Batente do fio guia • **ro** Opritor fir • **sk** Zarážka drôtu • **sl** Zaustavitveni nastavek žičnega vodila • **sv** Ledarspärr • **tr** Tel Durdurucu

3 Proximal Wire Port

bg Проксимален вход за теления водач • **cs** Proximální vodičí port • **da** Proksimal kateterlederport • **de** Proximaler Drahtport • **el** Εγγύς θύρα ούριατος (PWP) • **es** Acceso proximal de la guía • **et** Proksimaalne traadiport • **fi** Proksimaalinen johdinaukko • **fr** Orifice pour guide proximal • **hr** Proksimalni otvor za žicu • **hu** Proximális drótnyílás • **it** Ingresso prossimale della guida • **lt** Proksimalinis vielos prievadas • **lv** Proksimālais stigas ievadports • **nl** Proximale voerdraadpoort • **no** Proksimal vaierport • **pl** Proksymalny port dla prowadnika • **pt** Orifício proximal do fio • **ro** Port pentru fir proximal • **sk** Port na proximálny drôt • **sl** Proksimalni vhod za žično vodilo • **sv** Proximal ledarport • **tr** Proksimal Tel Portu

4 Retaining Clip

bg Задържащ клипс • **cs** Přidrzná svorka • **da** Holdeclips • **de** Halteklammer • **el** Κλιπ συγκράτησης • **es** Clip de retención • **et** Kinnitusklamber • **fi** Kiinnityspidike • **fr** Pince de retenue • **hr** Pričvrtna kopča • **hu** Rögzítőkapocs • **it** Clip di fissaggio • **lt** Fiksuojažiamasis spaustukas • **lv** Fiksācijas skava • **nl** Borgklem • **no** Holdeklips • **pl** Zacisk mocujący • **pt** Grampo de fixação • **ro** Clemă de retenție • **sk** Upevňovacia spona • **sl** Zadrževalna sponka • **sv** Clipshällare • **tr** Tutma Klipsi

5 Injection Port

bg Вход за инжектиране • **cs** Injekční port • **da** Injektionsport • **de** Injektionsport • **el** Θύρα έγχυσης • **es** Orificio de inyección • **et** Süstimisport • **fi** Ruisikutusaukko • **fr** Orifice d'injection • **hr** Otvor za ubrizgavanje • **hu** Befecskendezési nyílás • **it** Raccordo di iniezione • **lt** Injekcijos prievadas • **lv** Injekcijas pieslēgvieta • **nl** Injectiepoort • **no** Injeksjonsport • **pl** Port do wstrzykiwania • **pt** Orifício de injeção • **ro** Port de injectare • **sk** Vstrekovací port • **sl** Vhod za injiciranje • **sv** Injektionsport • **tr** Enjeksiyon portu

6 IDE Port

bg IDE вход • **cs** Port IDE • **da** IDE-port • **de** Intraduktaler Wechsellport (IDE) • **el** Θύρα IDE • **es** Acceso de intercambio intraductal (Intra Ductal Exchange, IDE) • **et** IDE-porti • **fr** Orifice IDE • **hr** IDE otvor • **hu** Intraduktális cserenyílás (IDE-nyílás) • **it** Accesso IDE • **lt** IDE prievadas • **lv** Intraduktālais nomainas (IDE) ievadports • **nl** IDE-poort • **no** IDE-port • **pl** Port wymiany wewnątrzprowadowej (IDE) • **pt** Orifício IDE • **ro** Port IDE • **sk** Port IDE • **sl** Vhod IDE • **sv** IDE-port • **tr** IDE Portu

7* Mid Point of Cutting Wire

bg Средна точка на режещата тел • **cs** Střed disekčního drátu • **da** Midtpunkt på skæretråd • **de** Mittelpunkt des Schneidedrahts • **el** Μεσαίο σημείο της χορδής τομής • **es** Punto medio del alambre de corte • **et** Lõikeetraadi keskpunkt • **fi** Leikkauslangan keskikohta • **fr** Point central du fil de coupe • **hr** Srednja točka rezne žice • **hu** A vágódrót középpontja • **it** Punto centrale del filo di taglio • **lt** Pjovimo vielos vidurio taškas • **lv** Griešanas stigas viduspunkts • **nl** Midden van snijdraad • **no** Skjæretrådens midtpunkt • **pl** Środek cięciwy tnącej • **pt** Ponto central do fio de corte • **ro** Punctul median al firului de tăiere • **sk** Stredový bod rezacieho drôtu • **sl** Srednja točka rezalne žice • **sv** Mittpunkt på skärtråden • **tr** Kesme Telinin Orta Noktası

8 Radiopaque Band


bg Рентгеноположителна лента • **cs** Rentgenkontrastní proužek • **da** Røntgenfast markering • **de** Röntgendichte Markierung • **el** Ακτινοοκιερή ταινία • **es** Banda radiopaca • **et** Röntgenkontrastne riba • **fi** Röntgenpositiivinen juova • **fr** Bande radio-opaque • **hr** Rendgenski vidljiva traka • **hu** Sugárfogó sáv • **it** Banda radiopaca • **lt** Rentgenokonstrastinė juostelė • **lv** Starojuma necaurlaidīga josla • **nl** Radiopake band • **no** Radiopakt bånd • **pl** Cieniodajny pasek • **pt** Banda radiopaca • **ro** Bandă radioopacă • **sk** Rádioopaký krúžok • **sl** Radioneprepustni obroček • **sv** Röntgentätt band • **tr** Radyopak Bant

9* Catheter

bg Катетър • **cs** Katetr • **da** Kateter • **de** Katheter • **el** Καθετήρας • **es** Catéter • **et** Kateteer • **fi** Katetri • **fr** Cathéter • **hr** Kateter • **hu** Katéter • **it** Catetere • **lt** Kateteris • **lv** Katetrs • **nl** Katheter • **no** Kateter • **pl** Cewnik • **pt** Cateter • **ro** Cateter • **sk** Katéter • **sl** Kateter • **sv** Kateter • **tr** Kateter

* Patient Contact

bg Контакт с пациента • **cs** Kontakt s pacientem • **da** Patientkontakt • **de** Patientenkontakt • **el** Τμήμα που έρχεται σε επαφή με τον ασθενή • **es** En contacto con el paciente • **et** Patsiendi kontakt • **fi** Potilasta koskettava osa • **fr** Contact patient • **hr** Kontakt s pacijentom • **hu** Beteggel érintkező • **it** Contatto con il paziente • **lt** Kontaktas su pacientu • **lv** Saskare ar pacientu • **nl** Patiëntcontact • **no** Pasientkontakt • **pl** Element wchodzący w kontakt z ciałem pacjenta • **pt** Contacto com o doente • **ro** Contact pacient • **sk** Kontakt s pacientom • **sl** Stik s pacientom • **sv** Patientkontakt • **tr** Hasta Teması

2.

1 Active Cord Connector (Metal Hot Pin)

bg Активен кабелен конектор (метален гъвкав щифт) • **cs** Konektor aktivního vodiče (kovový aktivní pin) • **da** Stik med strømførendende ledning (aktivt metalstikben) • **de** HF-Kabel-Anschluss (elektrisch aktiver Metallstift) • **el** Σύνδεσμος ενεργού καλωδίου (μεταλλική θερμή καρφίδα) • **es** Conector del cable de conexión (pin activo) • **et** Aktiivjuhtme konnektor (metallist aktiivtihvt) • **fi** Aktivointijohdon liitin (metallinen jännitteinen piikki) • **fr** Connecteur de cordon actif (broche métallique sous tension) • **hr** Priključak za aktivni kabel (metalna igla pod naponom) • **hu** Aktiv vezeték csatlakozója (fém aktiv tű) • **it** Connettore del cavo diatermico (pin metallico attivo) • **lt** Aktyviojo laidų jungtis (metalinis elektrinis kištukas) • **lv** Aktivā vada savienotājs (metāla karstā tapa) • **nl** Connector actief snoer (metalen hot pin) • **no** Kontakt for strømførende ledning (strømførende metallpinne) • **pl** Złącze przewodu zasilania (metalowy styk aktywny) • **pt** Conector do cabo ativo (pino metálico aquecido) • **ro** Conector cablu activ (pin metalic încălzit) • **sk** Konektor aktívnej žnúry (kovový aktivný kolík) • **sl** Priključek aktivnega kabla (kovinski aktivni priključek) • **sv** Kontakt till strömförande kabel (kontaktstift av metall) • **tr** Aktif Kablo Konektörü (Aktif Metal Pim)

2 Proximal Wire Port

bg Проксимален вход за теления водач • **cs** Proximální vodičí port • **da** Proximal kateterlederport • **de** Proximaler Drahtport • **el** Εγγύς θύρα σύρματος (PWP) • **es** Acceso proximal de la guía • **et** Proksimaalne traadiport • **fi** Proksimaalinen johdinaukko • **fr** Orifice pour guide proximal • **hr** Proksimalni otvor za žicu • **hu** Proximális drótnyílás • **it** Ingresso prossimale della guida • **lt** Proksimalinis vielos prievadas • **lv** Proksimālais stigas ievadports • **nl** Proximale voerdraadpoort • **no** Proksimal vaierport • **pl** Proksymalny port dla prowadnika • **pt** Orifício proximal do fio • **ro** Port pentru fir proximal • **sk** Port na proximálny drôt • **sl** Proksimalni vhod za žično vodilo • **sv** Proximal ledarport • **tr** Proximal Tel Portu

3 Injection Port

bg Вход за инжектиране • **cs** Injekční port • **da** Injektionsport • **de** Injektionsport • **el** Θύρα έγχυσης • **es** Orificio de inyección • **et** Süstimisport • **fi** Ruiskutusaukko • **fr** Orifice d'injection • **hr** Otvor za ubrizgavanje • **hu** Befecskendezési nyílás • **it** Raccordo di iniezione • **lt** Injekcijos prievadas • **lv** Injekcijas pieslēgvietā • **nl** Injectiepoort • **no** Injekjonsport • **pl** Port do wstrzykiwania • **pt** Orifício de injeção • **ro** Port de injectare • **sk** Vstrekovací port • **sl** Vhod za injiciranje • **sv** Injektionsport • **tr** Enjeksiyon portu

4 Slide Cover

bg Плъзгащ се капак • **cs** Posuvný kryt • **da** Glidehylster • **de** Schiebeabdeckung • **el** Συρόμενο κάλυμμα • **es** Cubierta deslizable • **et** Liuguri kate • **fi** Liukusuojuus • **fr** Protection coulissante • **hr** Klizni poklopac • **hu** Elcsúsztható fedél • **it** Copertura scorrevole • **lt** Slankusis dangtelis • **lv** Slidņa pārsegis • **nl** Schuifhoesje • **no** Glidehylster • **pl** Ruchoma pokrywa • **pt** Cobertura deslizante • **ro** Capac glisor • **sk** Posúvací kryt • **sl** Drsní pokrovček • **sv** Glidhölje • **tr** Sürgülü Kapak

5 Wire Control Port

bg Вход за контрол на тения водач • **cs** Port pro ovládání drátu • **da** Kateterlederkontrolport • **de** Drahtkontroll-Port • **el** Θύρα ελέγχου σύρματος • **es** Acceso de control de la guía • **et** Traadi juhtport • **fi** Johtimen ohjausaukko • **fr** Orifice pour contrôle du guide • **hr** Kontrolni otvor žice • **hu** Drótvézellő nyílás • **it** Apertura di controllo della guida • **lt** Vielos valdymo prievadas • **lv** Stīgas vadības pieslēgvietā • **nl** Draadcontrolepoort • **no** Vaierkontrollport • **pl** Port kontrolny przewodnika • **pt** Orifício de controlo do fio • **ro** Portul de control al firului • **sk** Port na ovládanie drôtu • **sl** Vhod za nadzor žičnega vodila • **sv** Ledarkontrollport • **tr** Tel Kontrol Portu

6* Mid Point of Cutting Wire

bg Средна точка на режещата тел • **cs** Střed disekčního drátu • **da** Midtpunkt på skæretråd • **de** Mittelpunkt des Schneidedrahts • **el** Μεσοίο σημείο της χορδής τομής • **es** Punto medio del alambre de corte • **et** Lõiketraadi keskpunkt • **fi** Leikkauslangan keskikohta • **fr** Point central du fil de coupe • **hr** Srednja točka rezne žice • **hu** A vágódrótt középpontja • **it** Punto centrale del filo di taglio • **lt** Pjovimo vielos vidurio taškas • **lv** Griešanas stīgas viduspunkts • **nl** Midden van snijdraad • **no** Skjæretrådens midtpunkt • **pl** Środek cięciwy tnącej • **pt** Ponto central do fio de corte • **ro** Punctul median al firului de tăiere • **sk** Stredový bod rezacieho drôtu • **sl** Srednja točka rezalne žice • **sv** Mittpunkt på skärtråden • **tr** Kesme Telinin Orta Noktası

7 Radiopaque Band

bg Рентгеноположителна лента • **cs** Rentgenkontrastní proužek • **da** Røntgenfast markering • **de** Röntgendichte Markierung • **el** Ακτινοσκοπική ταινία • **es** Banda radiopaca • **et** Röntgenkontrastne riba • **fi** Röntgenpositiivinen juova • **fr** Bande radio-opaque • **hr** Rendgenski vidljiva traka • **hu** Sugárfogó sáv • **it** Banda radiopaca • **lt** Rentgenokonstrastinė juostelė • **lv** Starojuma necaurīdīga josla • **nl** Radiopake band • **no** Radiopakt bånd • **pl** Cieniodajny pasek • **pt** Banda radiopaca • **ro** Bandă radioopacă • **sk** Rádioopaký krúžok • **sl** Radioneprepustni obroček • **sv** Röntgentätt band • **tr** Radyopak Bant

8* Catheter

bg Катетър • **cs** Katetr • **da** Kateter • **de** Katheter • **el** Καθετήρας • **es** Catéter • **et** Kateeter • **fi** Katetri • **fr** Cathéter • **hr** Kateter • **hu** Katéter • **it** Catetere • **lt** Kateteris • **lv** Katetrs • **nl** Katheter • **no** Kateter • **pl** Cewnik • **pt** Cateter • **ro** Cateter • **sk** Katéter • **sl** Kateter • **sv** Kateter • **tr** Kateter

* Patient Contact

bg Контакт с пациента • **cs** Kontakt s pacientem • **da** Patientkontakt • **de** Patientenkontakt • **el** Τίμημα που έρχεται σε επαφή με τον ασθενή • **es** En contacto con el paciente • **et** Patsiendi kontakt • **fi** Potilasta koskettava osa • **fr** Contact patient • **hr** Kontakt s pacijentom • **hu** Beteggel érintkező • **it** Contatto con il paziente • **lt** Kontaktas su pacientu • **lv** Saskare ar pacientu • **nl** Patiëntcontact • **no** Pasientkontakt • **pl** Element wchodzący w kontakt z ciałem pacjenta • **pt** Contacto com o doente • **ro** Contact pacient • **sk** Kontakt s pacientom • **sl** Stik s pacientom • **sv** Patientkontakt • **tr** Hasta Teması

Fusion® OMNI-TOME®**Fusion® OMNI-TOME® Sphincterotome****Fusion® Pre-Loaded OMNI-TOME®****Fusion® OMNI-TOME® Pre-Loaded Sphincterotome****Fusion® OMNI Pre-Loaded with Acrobat® 2 Wire Guide****Fusion® Triple Lumen Sphincterotome**

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION**Performance Characteristics**

Sphincterotomes consist of a three-ring handle with a metal hot pin for connection to an electrosurgical generator via an active cord. The hot pin delivers current to the cutting wire that is exposed at the distal end of the device. A multi-lumen catheter allows for injection of contrast or passage over a prepositioned wire guide using ports on the proximal end of the device. Devices contain an intraductal exchange (IDE) port for use with a short-wire technique.

Maximum rated input voltage for this device is 1.5 kVp-p (750 Vp).

Device Compatibility

Refer to the package label for device dimensions, minimum endoscope compatibility, and wire guide compatibility.

This device is recommended for use with Cook electrosurgical compatible wire guides.

This device must only be used with an active cord compatible with a 3 mm diameter connector.

This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

Electrosurgical generator - type BF or CF (generators that provide protection from electric shock in compliance with IEC60601-1).

Patient Population

This device is indicated for adult use only.

Intended Users

The sphincterotome should only be used by or under the supervision of physicians who have been trained in ERCP procedures.

Contact with Body Tissue

This device has contact with body tissue in line with the intended use.

Operating Principle

In a sphincterotomy, the sphincterotome is passed through a duodenoscope and positioned adjacent to the major papilla. The papilla is cannulated using the tip of the sphincterotome. The papilla is then incised using electrosurgical current applied to the cutting wire by an electrosurgical generator.

INTENDED USE

This device is used for cannulation of the ductal system and for sphincterotomy. If preloaded, also aids in bridging difficult strictures during ERCP.

INDICATIONS FOR USE

Sphincterotomes are used for selective deep cannulation of the ductal system and sphincterotomy to facilitate the completion of diagnostic and therapeutic endoscopic retrograde cholangiopancreatography.

CLINICAL BENEFITS

Sphincterotomes facilitate the completion of diagnostic procedures and interventions during endoscopic retrograde cholangiopancreatography (ERCP).

WARNINGS

- Do not use this device for any purpose other than stated intended use.
- Any electrosurgical device constitutes a potential electrical hazard to patient and operator. Possible adverse effects include: burns • cardiac arrhythmia • fulguration • nerve and/or muscle stimulation.
- The electrosurgical effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.
- Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper selection, placement, and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.
- For patients with electrically conductive implants, monopolar diathermy or electrosurgical cautery can result in localized thermal injury at the implant site.
- Do not use this device with an output higher than the rated voltage 1.5kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, sphincterotome and/or active cord.
- Do not use this device with an active cord which has a maximum voltage rating less than 1.5kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant.

- This device must be used with an endoscope compatible with high frequency equipment.
- If cutting is not observed at normal operating settings, this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.
- The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.
- Skin-to-skin contact should be avoided (for example between the patient's arms and body) by using dry gauze.
- Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.
- Avoid contact between the active cords and the patient's body, or any other electrodes.
- This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electric shock in compliance with IEC60601-1). Use of generator other than type BF or CF could cause electric shock to patient, operator, or assistant.
- For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.
- This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere, or in the presence of explosive gases. Doing so may result in explosion.
- In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electrosurgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the sphincterotome in these patients.
- This device must be used with an endoscope that is intended for use with high frequency (HF) electrosurgical devices. If the device is used with an endoscope not compatible with HF electrosurgery, a burn to the patient or the operator may occur.
- To avoid neuromuscular stimulation and minimize electric arc, activate the cutting wire only when it is in contact with the tissue to be cauterized.
- This device contains nickel which may cause an allergic reaction and other potential systemic complications related to nickel allergy.
- Additional warnings may be listed in procedural steps below.

Sterile and Single Use

This device is supplied sterile and designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device contamination with biological or chemical agents, electrical failure, and/or mechanical integrity failure of device.

CONTRAINDICATIONS

Contraindications of biliary sphincterotomy include general contraindications of ERCP, i.e., unstable or uncooperative patients.

Contraindications specific to biliary sphincterotomy include: newly created gastrointestinal anastomosis • the inability to orient the cutting wire of the sphincterotome towards the axis of the CBD • the patient's inability to provide informed consent • uncorrected coagulopathy.

PRECAUTIONS

Switch the electrosurgical unit to the "off" position when it is not in use.

Do not over flex or bow the tip beyond 90°, as this may damage or cause the cutting wire to break.

The elevator should remain open/down when advancing or retracting the sphincterotome.

If preloaded, use of a wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide.

This device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

TECHNICAL DESCRIPTION

See fig.1 for the Fusion Triple Lumen Sphincterotome and fig. 2 for the Fusion OMNI Sphincterotome. Components marked with an asterisk (*) denote applied parts (the part of the device that comes into contact with the patient, as defined by IEC60601-1), see figs. 1. and 2. Prior to starting the procedure in which the device will be used, become familiarized with distal mark location with respect to the cutting wire. Handle reference markings may be used for approximation of tip deflection.

This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6°F).

WARNING

No modification of this equipment is allowed.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • pancreatitis • respiratory depression or arrest • perforation • sepsis.

Potential Complications associated with ductal cannulation and bridging of strictures include: bleeding • pancreatitis • perforation • tissue inflammation.

HOW SUPPLIED

This device is supplied ethylene oxide (EtO) sterile in a peel-open pouch.

STORAGE

Store in a dark, dry location. Recommended storage temperature 22°C (71.6°F). Device may be exposed to temperatures up to 50°C (122°F) during transportation.

INSPECTION OF DEVICE

Do not use the device if the sterile packaging is damaged or unintentionally opened before use. Visually inspect with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

DEVICE PREPARATION

PRECAUTIONS

Do not apply manual pressure to the tip or the cutting wire of the sphincterotome in an attempt to influence orientation, as this may result in damage to the device.

Do not exercise the handle while the device is coiled or the precurved stylet is in place, as this may cause damage to the sphincterotome and render it inoperable.

1. Upon removing the device from the package, uncoil and straighten the sphincterotome. Carefully remove the precurved stylet wire from the cannulating tip.

WARNINGS

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

Before using this device, follow the recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout procedure.

2. With the electrosurgical unit off, prepare the equipment.
3. Attach the Wire Guide Locking Device to the endoscope accessory port.
4. Upon removing the device from the package, uncoil and straighten sphincterotome. Carefully remove the precurved stylet from the cannulating tip.
5. Prior to use, flush the catheter with sterile water or saline.

Note: If using the Fusion OMNI Sphincterotome (see fig. 2). If not preloaded with a wire guide, the Slide Cover should be slid down to cover the Wire Control Port prior to flushing to prevent leakage. If flushing over a wire guide that is preloaded through the Wire Control Port, the Slide Cover should be slid down to partially cover the Wire Control Port prior to flushing to minimize leakage. After flushing over a preloaded wire guide, the Slide Cover should be moved back proximal to the Wire Control Port to enable the wire guide to be moved freely.

INSTRUCTIONS FOR USE

WARNING

The sphincterotome should never be connected to the active cord prior to inserting the sphincterotome through the

endoscope. Doing so may injure the patient and/or damage the equipment due to improper electrical circuit grounding

1. Insert the catheter into the endoscope.

Note: If using the Intra Ductal Exchange (IDE) port

and short wire (See fig. 1) remove the Wire Stop from its retaining clip. Ensure that the Wire Stop is disengaged from the Proximal Wire Port (PWP) hub and that the distal tip of the Wire Stop is proximal to the IDE port. Insert the distal tip of the wire guide into the IDE port and advance until it is flush with the distal tip of the sphincterotome. Advance the Wire Stop handle until it reaches the PWP hub and connect the Luer lock securely to the PWP hub. Advance the tip of the sphincterotome through the cap of the Wire Guide Locking Device and continue advancing until it is endoscopically visible. Disengage the Wire Stop by releasing the handle from the PWP hub and remove the Wire Stop from the catheter.

Note: If using the Proximal Wire Port (PWP) and a

pre-positioned long wire guide (See fig. 1) remove the Wire Stop (if applicable). Advance the sphincterotome over the pre-positioned wire guide ensuring that the wire exits the catheter at the PWP. Continue advancing the device until it is endoscopically visible. For best results, wire guide should be kept wet with sterile water or saline.

Note: If using the Wire Control Port (WCP) and a

short wire or long wire guide (See fig. 2) Advance the tip of the sphincterotome through the cap of the Wire Guide Locking Device and continue advancing until it is endoscopically visible. If not preloaded, a wire guide may be introduced into the WCP and separated from the catheter.

CANNULATION

WARNING

Failure to confirm the cutting wire has fully exited endoscope may result in contact between the cutting wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

1. Using endoscopic visualization, confirm the cutting wire has exited the endoscope.
2. Cannulate using the sphincterotome. **Note:** A wire guide may be utilized to facilitate cannulation.
3. Following cannulation, confirm position fluoroscopically via contrast injection through the injection port or visualization of the radiopaque marker.

SPHINCTEROTOMY

WARNING

To avoid unintended disconnection during procedure and damage to the device, device should only be used with compatible active cords (Reference Device Compatibility for compatibility).

1. Attach the active cord to the connector on the sphincterotome handle. Active cord fittings should fit snugly into both the device handle and the electrosurgical unit.

2. If using a wire guide that is not compatible with electrosurgical procedures, remove the wire guide from the sphincterotome before activating the electrosurgical current. If the wire is not removed, there is a risk of thermal injury to the patient or operator.

WARNINGS

Do not use a wire guide (or sphincterotome) that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

The sphincterotome should not be energized prior to performing sphincterotomy. Energizing the cutting wire prior to use will cause premature cutting wire fatigue and will compromise the cutting wire's integrity.

Verify that the cutting wire has exited the endoscope by visualizing it on the endoscope monitor. Failure to do so may result in contact between the cutting wire and the endoscope while electrical current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire and/or damage to the endoscope.

User must maintain constant contact with tissue when applying electrocautery current. Failure to do so may result in increased current resulting in a broken cutting wire, damage to the endoscope and/or patient injury.

When applying current, ensure cutting wire is completely out of endoscope. Contact of cutting wire with endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope. If a wire guide that is not compatible with electrosurgical procedures is used in sphincterotome, it must be removed prior to applying electrosurgical current.

Cook wire guides packaged with this device may be left in place during sphincterotomy.

3. Following the electrosurgical unit manufacturer's instructions, verify desired settings and proceed with the sphincterotomy.
4. Upon completion of the sphincterotomy, turn the electrosurgical unit off.
5. Disconnect the active cord from the device handle and from the electrosurgical unit. Clean and store the active cord per the active cord manufacturer's instructions.

DEVICE REMOVAL FROM ENDOSCOPE

1. A previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices. If the wire guide is to remain in place while the device is withdrawn, utilize the following steps:

Note: If using the Intra Ductal Exchange (IDE) port and short wire (See fig. 1), prior to removing the device, utilize the reference marks on the catheter to ensure that the IDE port is within the ductal system. Fluoroscopically visualize the radiopaque band at the IDE port. Retract the wire guide until the radiopaque distal tip of the wire guide passes the band. A disengagement from the wire

guide lumen will occur. Advance the disengaged wire guide to maintain ductal access. Lock the wire guide into the Wire Guide Locking Device and remove the sphincterotome from the endoscope accessory channel.

Note: If using the Proximal Wire Port (PWP) and a pre-positioned long wire guide (See fig. 1), remove the device using standard long wire exchange technique.

Note: If using the Wire Control Port (WCP) and a short wire or long wire guide. To withdraw the sphincterotome from the endoscope, ensure the wire guide is locked in place (monitor the position of the wire guide at the Wire Guide Locking Device and distal end to ensure the wire guide is locked in place for the duration of the exchange), and then pull back on the catheter, allowing the wire guide lumen to separate from the wire guide until the metal band is visible at the Wire Guide Locking Device and resistance is felt. Unlock the wire guide from the Wire Guide Locking Device, completely remove the sphincterotome from the wire guide, and re-lock the wire guide.

DISPOSAL OF DEVICE

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, the competent authority, and/or regulatory authority of the country where the device was used.

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© COOK 2026
Wilson-Cook Medical, Inc.
Cook Medical Europe LTD
William A. Cook Australia Pty. Ltd.