COOK®

MEDICAL

CE
0123

EN
2

Classic ERCP Catheters

CS
3

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DA
4

**Klassiske katetre til endoskopisk retrograd
kolangiopancreatikografi (ERCP)**

NL
5

Klassieke ERCP-katheters

FR
6

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DE
7

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EL
8

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* 1 0 6 8 2 / 0 6 1 9 *

INTENDED USE

This device is used for endoscopic cannulation of the ductal system.

This device is indicated for adult use only.

NOTES

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP.

Those associated with ductal cannulation include, but are not limited to: bacteremia, bleeding, pancreatitis, tissue inflammation, and the inability to successfully cannulate.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device and to determine wire guide compatibility.

WARNINGS

This device is supplied sterile and intended for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

INSTRUCTIONS FOR USE

1. Gently advance the stylet completely into the catheter. Lock the stylet cap into the catheter fitting.
2. Introduce the catheter into the endoscope and advance in short increments until the device is endoscopically visualized exiting the scope.
3. Following cannulation, contrast may be injected to fluoroscopically confirm the position of the device by attaching a syringe to the side port of the catheter fitting.
4. To insert a compatible wire guide, remove the stylet from the catheter and advance the wire guide prior to removing the catheter from the endoscope.
Note: Before each introduction of a wire guide into a device, wet the wire guide in a sterile water or saline bath. A contrast-filled catheter may make the wire guide difficult to advance. Flush the contrast from the catheter with sterile water or saline before each introduction of the wire guide, observing necessary cautions if the device is positioned in the pancreatic duct.

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.