COOK®

MEDICAL

EN

2

Savary-Gilliard® Wire Guide

FR

3

Guide Savary-Gilliard®

PT

4

Fio guia Savary-Gilliard®


ES

5

Guía Savary-Gilliard®

* 1 0 6 8 6 / 0 1 2 0 *

ENGLISH

INTENDED USE

This device is used to introduce the Savary-Gilliard Dilator.

NOTES

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to upper GI endoscopy.

Contraindications to dilation include, but are not limited to: uncooperative patient; asymptomatic strictures; inability to advance the dilator through the strictured area; coagulopathy; known or suspected perforation; severe inflammation or scarring near the dilation site, recent myocardial infarction, active ulcer and severe cervical arthritis.

POTENTIAL COMPLICATIONS

Potential complications associated with upper gastrointestinal endoscopy and esophageal dilation include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

This device is supplied non-sterile, and therefore should not be used in a sterile field.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

INSTRUCTIONS FOR USE

1. Perform screening endoscopy and identify the strictured area.
2. Introduce the device, floppy tip first, into the accessory channel of the endoscope and advance until it is endoscopically visualized well beyond the tip of the scope. Example: Esophageal strictures in an unoperated GI tract would require wire guide placement in the distal body or antrum of the stomach.
3. When the wire guide is in position well beyond the strictured area, slowly begin to withdraw the endoscope in 5-10 cm increments while simultaneously advancing the wire guide in 5-10 cm increments to ensure the wire guide remains in position. Caution: Continuous fluoroscopic monitoring of the wire guide is essential in order to ensure it remains in the proper position.

4. When the endoscope is removed completely, fluoroscopically confirm that the wire guide has not been displaced. Note: Markings on the distal portion of the wire guide may also be utilized to determine the wire guide position. Markings occur in 20 cm increments from 40 to 140 cm. (See fig. 1)
5. Generously lubricate the dilator with water-soluble lubricant and advance it, under fluoroscopic monitoring, over, pre-positioned wire guide through the strictured area. Note: When multiple dilator sizes are used, repeat this step for each.
6. Upon completion of the procedure, remove the dilator and wire guide from the patient.

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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A symbol glossary can be found at cookmedical.com/symbol-glossary

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X mm

If symbol appears on product label, X = maximum insertion portion width

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X mm

If symbol appears on product label, X = minimum accessory channel

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X mm

If symbol appears on product label, X = wire guide compatibility

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If symbol appears on product label, X = quantity per box

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**This device is intended
for single use only.**

RxOnly

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA