**COOK®**


**MEDICAL**

**CE**  
**0123**

EN 3	<b>D.A.S.H.® Sphincterotome</b>
CS 5	<b>D.A.S.H.® sfinkterotom</b>
DA 7	<b>D.A.S.H.® sfinkterotom</b>
NL 10	<b>D.A.S.H.®-sfincterotoom</b>
FR 13	<b>Sphinctérotome D.A.S.H.®</b>
DE 15	<b>D.A.S.H.® Sphinkterotom</b>
EL 18	<b>Σφιγκτηροτόμος D.A.S.H.®</b>
HU 21	<b>D.A.S.H.® sphincterotóm</b>
IT 24	<b>Sfinterotomo D.A.S.H.®</b>
NO 27	<b>D.A.S.H.®-sfinkterotom</b>
PL 29	<b>Sfinkterotom D.A.S.H.®</b>
PT 32	<b>Esfincterótomo D.A.S.H.®</b>
ES 34	<b>Esfinterótomo D.A.S.H.®</b>
SV 37	<b>D.A.S.H.® sfinkterotom</b>

\* 1 0 7 5 6 / 0 8 1 9 \*

\* Patient Contact  
 Kontakt s pacientem  
 Patientkontakt  
 Contact met patiënt  
 Contact patient  
 Patientenkontakt  
 Τμήμα που έρχεται σε επαφή με τον ασθενή  
 A beteggel érintkező rész  
 A contatto con il paziente  
 Pasientkontakt  
 Element mający kontakt z ciałem pacjenta  
 Contacto com o doente  
 En contacto con el paciente  
 Patientkontakt

1 Tuohy-Borst adapter  
 Adaptér Tuohy-Borst  
 Tuohy-Borst-adapter  
 Tuohy-Borst-adapter  
 Adaptateur Tuohy-Borst  
 Tuohy-Borst-Adapter  
 Προσαρμογέας Tuohy-Borst  
 Tuohy-Borst-adapter  
 Adattatore Tuohy-Borst  
 Tuohy-Borst-adapter  
 Łącznik Tuohy-Borst  
 Adaptador Tuohy-Borst  
 Adaptador Tuohy-Borst  
 Tuohy-Borst-adapter

2 Wire Guide Port  
 Port pro vodič drát  
 Kateterlederport  
 Voerdraadpoort  
 Orifice pour guide  
 Führungsdraht-Port  
 Θύρα συμπίπτουσα οδηγού  
 Vezetődrótnyílás  
 Raccordo per la guida  
 Ledevierport  
 Port dla prowadnika  
 Orifício do fio guia  
 Orificio para la guía  
 Ledarport

3 Injection Port  
 Injekční port  
 Injektionsport  
 Injectiepoort  
 Orifice d'injection  
 Injektionsport  
 Θύρα έγχυσης  
 Befecskendezési nyílás  
 Raccordo di iniezione  
 Injektionsport  
 Port do wstrzykiwań  
 Orifício de injeção  
 Orificio de inyección  
 Injektionsport

4 Reference Markings  
 Referenční značky  
 Referencemarker  
 Referentiemarkeringen  
 Repères de référence  
 Bezugsmarkierungen  
 Σημάνσεις αναφοράς  
 Referenciajelzések  
 Marker di riferimento  
 Referansemarkører  
 Oznaczenia odnośne  
 Marcas de referência  
 Marcas de referencia  
 Referensmarkeringar

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## ENGLISH

### INTENDED USE

This device is used for cannulation of the ductal system and for sphincterotomy. If preloaded, also aids in bridging difficult strictures during ERCP.

### CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with sphincterotomy.

Contraindications to sphincterotomy include, but are not limited to: coagulopathy and inability to properly position the sphincterotome cutting wire.

### POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

### USER QUALIFICATIONS

The sphincterotome should only be used by or under the supervision of physicians who have been trained in ERCP procedures.

### WARNINGS

Do not use this device for any purpose other than stated intended use.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Any electrological device constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

The electrological effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.

Before using this device, follow recommendations provided by the electrological unit manufacturer to ensure patient safety through proper selection, placement, and utilization of patient return electrode. Ensure a proper path from patient return electrode to electrological unit is maintained throughout procedure.

For patients with electrically conductive implants, monopolar diathermy or electrological cautery can result in localized thermal injury at the implant site.

Do not use this device with an output higher than the rated voltage 1.5kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, sphincterotome and/or active cord.

Do not use this device with an active cord which has a maximum voltage rating less than 1.5kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant. If cutting is not observed at normal operating settings this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.

The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.

Skin-to-skin contact should be avoided (for example between the arms and body of the patient) by using dry gauze.

Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.

Avoid contact between active cords and the patient's body, or any other electrodes.

This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electric shock in compliance with IEC60601-1). Use of generator other than type BF or CF could cause electric shock to patient, operator, or assistant.

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere or in the presence of explosive gases. Doing so may result in explosion.

In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electrological cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the sphincterotome in these patients.

This device must be used with an endoscope that is intended for use with high frequency (HF) electrological devices. If the device is used with an endoscope not compatible with HF electrological surgery, a burn to the patient or the operator may occur.

To avoid neuromuscular stimulation and minimize electric arc, activate the cutting wire only when the cutting wire is in contact with the tissue to be cauterized.

Additional warnings may be listed in procedural steps below.

### PRECAUTIONS

Switch the electrological unit to "off" position when it is not in use.

Do not over flex or bow the tip beyond 90°, as this may damage or cause the cutting wire to break.

The elevator should remain open/down when advancing or retracting the sphincterotome.

If preloaded, use of a wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide.

### TECHNICAL DESCRIPTION

See fig. 1 for DASH Sphincterotomes. Components marked with an asterisk (\*) denote applied parts (the part of the device that comes into contact with the patient, as defined by IEC60601-1), see fig. 1. Prior to starting the procedure in which, the device will be used, become familiarized with distal mark location with respect to the cutting wire. Handle reference markings may be used for approximation of tip deflection.

Refer to the package label for device dimensions, minimum endoscope channel compatibility, and wire guide compatibility.

This device is recommended for use with Cook electrosurgical compatible wire guides.

This device must only be used with an active cord compatible with a 3 mm diameter connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

Maximum rated input voltage for this device is 1.5 kVp-p (750 Vp).

This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6°F).

Store in a dry location. Recommended storage temperature 22°C (71.6° F). Device may be exposed to temperatures up to 50°C (122°F) during transportation.

### WARNING

No modification of this equipment is allowed.

### SYSTEM PREPARATION

#### PRECAUTIONS

Do not apply manual pressure to the tip or the cutting wire of the sphincterotome in an attempt to influence orientation, as this may result in damage to the device.

Do not exercise the handle while the device is coiled or the precurved stylet is in place, as this may cause damage to the sphincterotome and render it inoperable.

1. Upon removing the device from the package, uncoil and straighten the sphincterotome. Carefully remove the precurved stylet wire from the cannulating tip.

### WARNINGS

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

2. With the electrosurgical unit off, prepare the equipment.

3. Verify the wire guide position within the sphincterotome. The distal tip of the wire guide should be even with the tip of the catheter.
4. Once the distal tips of the wire guide and the device are aligned, secure the wire guide in position with the Tuohy-Borst adapter to avoid risk of wire movement.
5. Flush the injection port with sterile water or saline. For best results, the wire guide should be kept wet. **Note:** Make sure the Tuohy-Borst adapter is securely tightened around the wire guide prior to injecting.

### INSTRUCTIONS FOR USE

#### WARNING

The sphincterotome should never be connected to the active cord prior to inserting the sphincterotome through the endoscope. Doing so may injure the patient and/or damage the equipment because of improper electrical circuit grounding.

1. Advance the tip of the sphincterotome into the endoscope accessory channel and continue to advance in short increments until the device is endoscopically visible.

### CANNULATION

#### WARNING

Failure to confirm the cutting wire has fully exited endoscope may result in contact between the cutting wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

1. Using endoscopic visualization, confirm cutting wire has exited the endoscope.
2. Cannulate with the sphincterotome. **Note:** A wire guide may be utilized to facilitate cannulation, if desired.
3. Fluoroscopically monitor the wire guide advancement in the ductal system. The radiopaque tip of the wire guide may be used to fluoroscopically verify position of the wire guide.
4. Following cannulation, contrast may be injected through the injection port to fluoroscopically confirm the position of the device. **Note:** The Tuohy-Borst adapter must be tightened when injecting to prevent backflow. When injecting, use a 10 cc syringe (or smaller) to assure proper pressure.

### SPHINCTEROTOMY

#### WARNING

To avoid unintended disconnection during procedure and damage to the device, device should only be used with compatible active cords (Reference Technical Description for compatibility).

1. Attach the active cord to the connector on the sphincterotome handle. Active cord fittings should fit snugly into both the device handle and the electrosurgical unit.
2. If using a wire guide that is not compatible with electrosurgical procedures, remove the wire guide from the sphincterotome before activating the electrosurgical

current. If the wire is not removed, there is a risk of thermal injury to the patient or operator.

## WARNINGS

Do not use a wire guide (or sphincterotome) that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

The sphincterotome should not be energized prior to performing sphincterotomy. Energizing the cutting wire prior to use will cause premature cutting wire fatigue and will compromise the cutting wire's integrity.

Verify that the cutting wire has exited the endoscope by visualizing it on the endoscope monitor. Failure to do so may result in contact between the cutting wire and the endoscope while electrical current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire and/or damage to the endoscope.

User must maintain constant contact with tissue when applying electrocautery current. Failure to do so may result in increased current resulting in broken cutting wire, damage to the endoscope and/or patient injury.

When applying current, ensure cutting wire is completely out of endoscope. Contact of cutting wire with endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope.

If a wire guide that is not compatible with electrosurgical procedures is used in sphincterotome, it must be removed prior to applying electrosurgical current.

Cook wire guides packaged with this device may be left in place during sphincterotomy.

3. Following the electrosurgical unit manufacturer's instructions, verify desired settings and proceed with the sphincterotomy.
4. Upon completion of the sphincterotomy, turn the electrosurgical unit off.
5. Disconnect the active cord from the device handle and from the electrosurgical unit. Clean and store the active cord per the active cord manufacturer's instructions.

## DEVICE REMOVAL FROM ENDOSCOPE

1. A previously placed wire guide may be left in position, in order to facilitate introduction of compatible devices (if applicable).

## DEVICE DISPOSAL

1. Upon completion of the procedure, dispose of the device(s) per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**