**COOK®**

**MEDICAL**



**CE**  
**0123**

**CCPT, CT, UTS**



<b>EN</b> <b>3</b>	<b>PreCurved Double Lumen Sphincterotome</b>
<b>CS</b> <b>5</b>	<b>Dvojlumenný preformovaný sfinkterotom</b>
<b>DA</b> <b>8</b>	<b>Buede sfinkterotomer med dobbelt lumen</b>
<b>NL</b> <b>10</b>	<b>Voorgevormde sfincterotoom met dubbel lumen</b>
<b>FR</b> <b>13</b>	<b>Sphinctérotome à double lumière précourbé</b>
<b>DE</b> <b>15</b>	<b>Vorgebogenes doppelumiges Sphinkterotom</b>
<b>EL</b> <b>18</b>	<b>Προκεκαμμένος σφιγκτηροτόμος διπλού αυλού</b>
<b>HU</b> <b>21</b>	<b>Előhajlított kettős lumenű sphincterotom</b>
<b>IT</b> <b>23</b>	<b>Sfinterotomo precurvato a doppio lume</b>
<b>NO</b> <b>26</b>	<b>Forhåndsbøyd sfinkterotom med dobbeltlumen</b>
<b>PL</b> <b>28</b>	<b>Zakrzywiony sfinkterotom dwukanałowy</b>
<b>PT</b> <b>31</b>	<b>Esfinterótomo pré-curvado de duplo lúmen</b>
<b>ES</b> <b>33</b>	<b>Esfinterótomo de doble luz precurvado</b>
<b>SV</b> <b>36</b>	<b>Förböjda sfinkterotomer med dubbellumen</b>


\* 1 0 7 5 8 / 0 8 1 9 \*

- |  |  |  |  |  |
|--|--|--|--|--|
| <p><b>1</b> Handle Reference Markers<br/>Referenční značky rukojeti<br/>Håndtagets referencemarkører<br/>Referentiemarkeringen op handgriep<br/>Repères de référence sur la poignée<br/>Bezugsmarkierungen am Griff<br/>Δείκτες αναφοράς λαβής<br/>Fogantyú referencialelei<br/>Marker di riferimento dell'impugnatura<br/>Referansemarkører på håndtak<br/>Marcadores de referência do punho<br/>Marcadores de referencia del mango<br/>Handtagets referensmarkörer</p> | <p><b>2</b> Wire Guide Port<br/>Port pro vodiči drát<br/>Kateterlederport<br/>Voerdraadpoort<br/>Orifice pour guide<br/>Führungsdraht-Port<br/>Θύρα συμμάτνιου οδηγού<br/>Vezetődrótnyílás<br/>Raccordo per la guida<br/>Ledevalierport<br/>Port prowadnika<br/>Orificio do fio guia<br/>Orificio para la guía<br/>Ledarport</p> | <p><b>3</b> Injection Port<br/>Injekční port<br/>Injektionsport<br/>Injctiepoort<br/>Orifice d'injection<br/>Injektionsport<br/>Θύρα έγχυσης<br/>Befecskendezési nyílás<br/>Raccordo di iniezione<br/>Injektsjonsport<br/>Port do wstrzykiwania<br/>Orificio de injeção<br/>Orificio de inyección<br/>Injektionsport</p> | <p><b>4*</b> Mid Point of Cutting Wire<br/>Střed disekčního drátu<br/>Midtpunkt på skærewire<br/>Midden van snijdraad<br/>Point central du fil de coupe<br/>Mittelpunkt des Schneidedrahts<br/>Μεσοίο σημείο της χορδής τομής<br/>A vágódrót középontja<br/>Punto centrale del filo di taglio<br/>Skjæretrådets midtpunkt<br/>Port do wstrzykiwania<br/>Ponto central do fio de corte<br/>Punto medio del alambre de corte<br/>Mittpunkt på skärtråden</p> | <p><b>5*</b> Catheter*<br/>Katetr<br/>Kateter*<br/>Katheter*<br/>Cathéter<br/>Katheter*<br/>Καθετήρας*<br/>Katéter*<br/>Catetere*<br/>Kateter*<br/>Cewnik*<br/>Cateter<br/>Catéter*<br/>Kateter*</p> |
|--|--|--|--|--|

- |  |   |   |   |
|--|---|---|---|
| <p><b>1</b> Injection Port / Wire Guide Port<br/>Injekční port / Port pro vodiči drát<br/>Injektionsport/kateterlederkontrolport<br/>Injctiepoort / Voerdraadpoort<br/>Orifice d'injection/Orifice pour guide<br/>Injektionsport / Führungsdraht-Port<br/>Θύρα έγχυσης / Θύρα συμμάτνιου οδηγού<br/>Befecskendezési nyílás / Vezetődrótnyílás<br/>Raccordo di iniezione / raccordo per la guida<br/>Injektsjonsport/ledevalierport<br/>Port do wstrzykiwania / Port prowadnika<br/>Orificio de injeção/Orificio do fio guia<br/>Orificio de inyección/Orificio para la guía<br/>Injektionsport/Ledarport</p> | <p><b>2*</b> Catheter<br/>Katetr<br/>Kateter<br/>Katheter<br/>Cathéter<br/>Katheter<br/>Καθετήρας<br/>Katéter<br/>Catetere<br/>Kateter<br/>Cewnik<br/>Cateter<br/>Catéter<br/>Kateter</p> | <p><b>3*</b> Cutting Wire<br/>Disekční drát<br/>Skærertråd<br/>Snijdraad<br/>Fil de coupe<br/>Schneidedraht<br/>Χορδή τομής<br/>Vágódrót<br/>Filo di taglio<br/>Skjæretråd<br/>Cięciwa tnąca<br/>Fio de corte<br/>Alambre de corte<br/>Skärtråd</p> | <p><b>*</b> Patient Contact<br/>Kontakt s pacientem<br/>Patientkontakt<br/>Contact met patiënt<br/>Contact patient<br/>Patientenkontakt<br/>Τμήμα που έρχεται σε επαφή με τον ασθενή<br/>A beteggel érintkező rész<br/>A contatto con il paziente<br/>Pasientkontakt<br/>Element wchodzący w kontakt z ciałem pacjenta<br/>Contacto com o doente<br/>En contacto con el paciente<br/>Patientkontakt</p> |
|--|---|---|---|

1	Slide	2	On	3	Off
	Jezdec		Zapnuto		Vypnuto
	Glidehyster		Til		Fra
	Schuif		Aan		Uit
	Glissière		Marche		Arrêt
	Schieber		Ein		Aus
	Ολοθητήρας		Ανοικτός		Κλειστός
	Retesz		Be		Ki
	Elemento scorrevole		Attivazione		Disattivazione
	Glideknapp		På		Av
	Suwak		Włączony		Wyłączony
	Peça deslizante		Ligado		Desligado
	Interruptor deslizante		Encendido		Apagado
	Glidknapp		Till		Från

## ENGLISH

### INTENDED USE

This device is used for cannulation of the ductal system and for sphincterotomy.

### CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with sphincterotomy.

Contraindications to sphincterotomy include, but are not limited to: coagulopathy and inability to properly position sphincterotome cutting wire.

### POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

### USER QUALIFICATIONS

The sphincterotome should only be used by or under the supervision of physicians who have been trained in ERCP procedures.

### WARNINGS

Do not use this device for any purpose other than stated intended use.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of diseases.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Any electrocautery device constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

The electrocautery effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.

Before using this device, follow recommendations provided by the electrocautery unit manufacturer to ensure patient safety through proper selection, placement, and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrocautery unit is maintained throughout the procedure.

For patients with electrically conductive implants, monopolar diathermy or electrocautery can result in localized thermal injury at the implant site.

Do not use this device with an output higher than the rated voltage 1.5kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, sphincterotome and/or active cord.

Do not use this device with an active cord which has a maximum voltage rating less than 1.5kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant.

If cutting is not observed at normal operating settings this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.

The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.

Skin-to-skin contact should be avoided (for example between the arms and body of the patient) by using dry gauze.

Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.

Avoid contact between the active cords and the patient's body, or any other electrodes.

This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electric shock in compliance with IEC60601-1). Use of generator other than type BF or CF could cause electric shock to patient, operator, or assistant.

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere or in the presence of explosive gases. Doing so may result in explosion.

In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electro-surgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the sphincterotome in these patients.

This device must be used with an endoscope that is intended for use with high frequency (HF) electro-surgical devices. If the device is used with an endoscope not compatible with HF electro-surgery, a burn to the patient or the operator may occur.

To avoid neuromuscular stimulation and minimize electric arc, activate the cutting wire only when it is in contact with the tissue to be cauterized.

Additional warnings may be listed in procedural steps below.

## PRECAUTIONS

Switch the electro-surgical unit to "off" position when it is not in use.

Do not over flex or bow the tip beyond 90°, as this may damage or cause the cutting wire to break.

The elevator should remain open/down when advancing or retracting the sphincterotome.

## TECHNICAL DESCRIPTION

See fig. 1 for CT devices, see fig. 2 for CCPT and UTS devices, see fig. 3 for CT-25M-P slide. Components marked with an asterisk (\*) denote applied parts (the part of the device that comes into contact with the patient, as defined by IEC60601-1), see figs.1 and 2. Prior to starting the procedure in which the device will be used, become familiarized with distal mark location with respect to the cutting wire. Handle reference markings may be used for approximation of tip deflection.

Refer to the package label for device dimensions, minimum endoscope channel compatibility, and wire guide compatibility.

This device is recommended for use with Cook electro-surgical compatible wire guides.

This device must only be used with an active cord compatible with a 3 mm diameter connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

Maximum rated input voltage for this device is 1.5 kVp-p (750 Vp).

This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6° F)

Store in a dry location. Recommended storage temperature 22°C (71.6°F). Device may be exposed to temperatures up to 50°C (122°F) during transportation.

## WARNING

No modification of this equipment is allowed.

## SYSTEM PREPARATION

### PRECAUTIONS

Do not apply manual pressure to tip or cutting wire of sphincterotome in an attempt to influence orientation, as this may result in damage to device.

Do not exercise the handle while the device is coiled or the precurved stylet is in place, as this may cause damage to the sphincterotome and render it inoperable.

1. Upon removing device from package, uncoil and straighten the sphincterotome. Carefully remove the precurved stylet from the cannulating tip.

## WARNINGS

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

Before using this device, follow recommendations provided by the electro-surgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electro-surgical unit is maintained throughout the procedure.

2. With the electro-surgical unit off, prepare the equipment.

## INSTRUCTIONS FOR USE

### WARNING

The sphincterotome should never be connected to the active cord prior to inserting the sphincterotome through the endoscope to avoid injury to patient and/or equipment resulting from improper electrical circuit grounding.

1. Advance the tip of sphincterotome into the endoscope accessory channel and continue to advance in short increments until the device is endoscopically visible. Note: This sphincterotome may be placed over a pre-positioned wire guide. Prior to advancing the device over a pre-positioned wire guide, flush the wire guide lumen with sterile water or saline.

## CANNULATION

### WARNING

Failure to confirm the cutting wire has fully exited the endoscope may result in contact between the cutting wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

1. Using endoscopic visualization, confirm cutting wire has exited the endoscope.
2. Cannulate with the sphincterotome. **Note:** A wire guide may be utilized to facilitate cannulation, if desired.
3. Following cannulation, contrast may be injected to fluoroscopically confirm position of device. **Note:** If applicable, move the slide to "off" position prior to injecting over a pre-positioned wire guide. To manipulate the wire guide, move the slide to "on" position (See fig. 3).

## SPHINCTEROTOMY

### WARNING

To avoid unintended disconnection during procedure and damage to the device, device should only be used with compatible active cords (Reference Technical Description for compatibility).

1. Attach the active cord to the connector on the sphincterotome handle. Active cord fittings should fit snugly into both the device handle and the electrosurgical unit.
2. If using a wire guide that is not compatible with electrosurgical procedures, remove the wire guide from the sphincterotome before activating the electrosurgical current. If the wire is not removed, there is a risk of thermal injury to the patient or operator.

### WARNINGS

Do not use a wire guide (or sphincterotome) that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

The sphincterotome should not be energized prior to performing sphincterotomy. Energizing the cutting wire prior to use will cause premature cutting wire fatigue and will compromise the cutting wire's integrity.

Verify that the cutting wire has exited the endoscope by visualizing it on the endoscope monitor. Failure to do so may result in contact between the cutting wire and the endoscope while electrical current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire and/or damage to the endoscope.

User must maintain constant contact with tissue when applying electrocautery current. Failure to do so may result in increased current resulting in a broken cutting wire, damage to the endoscope and/or patient injury

When applying current, ensure cutting wire is completely out of endoscope. Contact of cutting wire with endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope.

If a wire guide that is not compatible with electrosurgical procedures is used in sphincterotome, it must be removed prior to applying electrosurgical current.

3. Following the electrosurgical unit manufacturer's instructions, verify the desired settings and proceed with the sphincterotomy.
4. Upon completion of the sphincterotomy, turn the electrosurgical unit off.
5. Disconnect the active cord from the device handle and from the electrosurgical unit. Clean and store the active cord per the active cord manufacturer's instructions.

## DEVICE REMOVAL FROM ENDOSCOPE

1. A previously placed wire guide may be left in position, in order to facilitate introduction of compatible devices (if applicable).

## DEVICE DISPOSAL

1. Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**