<b>EN</b> 2	<b>Classic ERCP Catheter</b>
<b>CS</b> 2	<b>Katetr Classic pro ERCP</b>
<b>DA</b> 3	<b>Classic-kateter til endoskopisk retrograd kolangiopancreatikografi</b>
<b>NL</b> 4	<b>Classic ERCP-katheter</b>
<b>FR</b> 5	<b>Cathéter de CPRE Classic</b>
<b>DE</b> 6	<b>Classic ERCP-Katheter</b>
<b>EL</b> 7	<b>Καθετήρας Classic ERCP</b>
<b>HU</b> 8	<b>Classic ERCP-katéter</b>
<b>IT</b> 9	<b>Catetere per ERCP Classic</b>
<b>NO</b> 10	<b>Classic ERCP-kateter</b>
<b>PL</b> 11	<b>Cewnik do ECPW Classic</b>
<b>PT</b> 12	<b>Cateter de CPRE Classic</b>
<b>ES</b> 13	<b>Catéter de ERCP Classic</b>
<b>SV</b> 14	<b>Classic ERCP-kateter</b>

**INTENDED USE**

This device is used for endoscopic cannulation of the ductal system.

This device is indicated for adult use only.

**NOTES**

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device is restricted to a trained healthcare professional.

**CONTRAINDICATIONS**

Those specific to ERCP.

**POTENTIAL COMPLICATIONS**

Those associated with ERCP.

Those associated with ductal cannulation include, but are not limited to: bacteremia, bleeding, pancreatitis, tissue inflammation, and the inability to successfully cannulate.

**PRECAUTIONS**

Refer to the package label for the minimum channel size required for this device.

This device is not wire guide compatible.

**WARNINGS**

This device is supplied sterile and intended for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

**INSTRUCTIONS FOR USE**

1. Gently advance the stylet completely into the catheter. Lock the stylet cap into the catheter fitting.
2. Introduce the catheter into the endoscope and advance in short increments until the device is endoscopically visualized exiting the scope.
3. Following cannulation, contrast may be injected to fluoroscopically confirm the position of the device by attaching a syringe to the side port of the catheter fitting.

**Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.**

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**