EN
5

AcuSnare® Polypectomy Snare
Instructions for Use

CS
5

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

AcuSnare® Polypectomy Snare

AcuSnare® Polypectomy Snare Soft

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Performance Characteristics

The maximum rated input voltage for this device is 2.0 kVp-p (1000 Vp) for cut mode and 5 kVp-p (2500 Vp) for coagulation mode.

Patient Population

The device is indicated for adult use only.

Intended User

This device should only be used by or under the supervision of physicians who have been trained in endoscopic procedures.

Contact with Body Tissue

Tissue contact for this device is in line with the intended use.

Operating Principle

The snare is inserted into the endoscopic accessory channel and advanced to the polyp target site. The snare head is deployed and looped around the polyp. Once snared, the polyp may be removed en bloc or in a piecemeal fashion through a combination of loop closure pressure and electrosurgical cutting and/or cauterization. The excised polyp may be left to pass naturally or may be captured with a polyp retrieval basket and removed.

INTENDED USE

This device is used endoscopically in the removal and cauterization of sessile polyps and pedunculated polyps from within the gastrointestinal tract.

INDICATIONS FOR USE

The Polypectomy devices are used by trained professionals in conjunction with a compatible endoscope and electrosurgical unit for resection of sessile and/or pedunculated polyps in the gastrointestinal tract.

CLINICAL BENEFITS

The Polypectomy devices provide the clinical benefit of polyp removal during an endoscopic procedure.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the polypectomy site.

Contraindications to polypectomy include: coagulopathy • insufficiently prepped bowel.

WARNINGS

- Please read all instructions before using this device.
- Do not use this device for any purpose other than the stated intended use.
- Any electrosurgical device constitutes a potential electrical hazard to the patient and operator. Possible adverse effects include: burns • cardiac arrhythmia • fulguration • nerve and/or muscle stimulation.
- The electrosurgical effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.
- Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through the proper selection, placement and utilization of the patient return electrode. Ensure that a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.
- For patients with electrically conductive implants, monopolar diathermy or electrosurgical cautery can result in localized thermal injury at the implant site.
- Do not use this device with an output higher than the rated voltage 2.0 kVp-p (1000 Vp) for cut mode and 5 kVp-p (2500 Vp) for coagulation mode. This could cause thermal injury to the patient, operator or assistant and could also damage the endoscope, snare, and/or active cord.
- Do not use this device with an active cord which has a maximum voltage rating less than 5kVp-p (2500 Vp). This could cause thermal injury to the patient, operator, or assistant.
- This device must be used with an endoscope compatible with high frequency equipment.
- If cutting is not observed at normal operating settings, this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.
- The patient should not come into contact with metal parts that are earthed or that have an appreciable capacitance to earth.
- Skin-to-skin contact should be avoided (for example between the patient's arms and body) by using dry gauze.
- Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.
- Avoid contact between the active cords and the patient's body or any other electrodes.
- This device must only be used in conjunction with a type BF or CF generator (generators that provide protection

from electric shock in compliance with IEC60601-1). Use of a generator other than type BF or CF could cause electric shock to patient, operator, or assistant.

- For recommendations on current settings for polypectomy, refer to the monopolar generator's recommended power settings for polypectomy.
- This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere, or in the presence of explosive gases. Doing so may result in an explosion.
- In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electrosurgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the snare in these patients.
- This device must be used with an endoscope that is intended for use with high frequency (HF) electrosurgical devices. If the device is used with an endoscope not compatible with HF electrosurgery, a burn to the patient or the operator may occur.
- To avoid neuromuscular stimulation and minimize electric arc, activate the snare wire only when it is in contact with the tissue to be cauterized.
- This device contains nickel which may cause an allergic reaction and other potential systemic complications related to nickel allergy.
- Additional warnings may be listed in procedural steps below.

• **Sterile and Single Use**

This device is supplied sterile and designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device contamination with biological or chemical agents, electrical failure, and/or mechanical integrity failure of device.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

Polypectomy should be performed under direct endoscopic visualization.

This device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • fever • hemorrhage •

hypotension • infection • perforation • respiratory depression or arrest.

Complications which can occur with endoscopic polypectomy include: explosion • thermal injury to the patient • transmural burns.

TECHNICAL DESCRIPTION

See fig 1. Components highlighted with an asterisk (*) denote applied parts (the part of the device that comes into contact with the patient). Handle reference markings may be used for approximation of snare extension.

See fig 2. Snares with oval, hexagonal, or duckbill snare head are available for use based on physician preference. Snares ranging in size (width x length) from 1 cm x 1.5 cm to 3 cm x 6 cm are available for use based on physician preference related to performing polypectomy on target anatomy. This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6° F).

WARNING

No modification of this equipment is allowed.

Device Compatibility

Electrosurgical generator - type BF or CF (generators that provide protection from electric shock in compliance with IEC60601-1)

Refer to the package label for device dimensions and minimum endoscope compatibility.

This device must only be used with an active cord compatible with a 3 mm diameter connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

HOW SUPPLIED

This device is supplied ethylene oxide (EtO) sterile in a peel-open pouch. There are 10 pouches to a box. Refer to Storage and Technical Description for storage and handling conditions.

STORAGE

Store in a dry location. Recommended storage temperature 22°C (71.6° F). This device may be exposed to temperatures up to 50°C (122° F) during transportation.

INSPECTION OF DEVICE

Do not use the device if the sterile packaging is damaged or unintentionally opened before use. Visually inspect with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

DEVICE PREPARATION

1. Inspect the active cord. The cord must be free of kinks, bends, breaks, and exposed wires to allow for the accurate transfer of current. If an abnormality is noted, do not use the active cord.

WARNINGS

For recommendations on current settings for polypectomy, refer to the monopolar generator's recommended power settings for polypectomy.

Before using this device, follow the recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through the proper placement and utilization of the patient return electrode. Ensure that a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

2. With the electrosurgical unit off, prepare the equipment. Following the instructions from the electrosurgical unit manufacturer, position the patient return electrode and connect it to the electrosurgical unit.
3. Fully retract and extend the snare to confirm smooth operation of the device.

INSTRUCTIONS FOR USE

WARNINGS

The snare should never be connected to the active cord prior to inserting the snare through the endoscope. Doing so may injure the patient and/or damage the equipment due to improper electrical circuit grounding.

To ensure patient safety, power to the electrosurgical unit should remain off until the snare is properly positioned around the polyp.

1. When the polyp is in endoscopic view, introduce the sheath and retracted snare into the endoscope accessory channel.
2. Advance the device, in small increments, until it is endoscopically viewed exiting endoscope.
3. Advance the snare wire out of the sheath and position it around the polyp to be removed.

WARNINGS

Do not use a snare that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

The snare should not be energized prior to performing polypectomy. Energizing the snare prior to use will cause premature wire fatigue and will compromise the wire's integrity.

Using the endoscope monitor, verify that the snare wire has exited the endoscope by visually ensuring that the distal end of the catheter has exited the endoscope. Failure to confirm the snare wire has fully exited the endoscope may result in contact between the snare wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken snare wire, and/or damage to the endoscope.

When applying current, tissue must be isolated from the surrounding mucosa. Failure to isolate tissue may cause fulguration of normal mucosa and/or perforation.

4. Following the electrosurgical unit manufacturer's instructions for settings, verify the desired settings and activate the electrosurgical unit. **Note:** The maximum rated input voltage for this device is 2 kVp-p for cut mode and 5 kVp-p for coagulation mode.

WARNINGS

To avoid unintended disconnection during the procedure and damage to the device, the device should only be used with compatible active cords, (Reference Technical Description for Device Compatibility.)

5. Securely connect the active cord to the device handle and electrosurgical unit. The active cord fittings should fit snugly into both the device handle and the electrosurgical unit.
6. Proceed with the polypectomy.
7. Upon completion of the polypectomy, turn the electrosurgical unit off.
8. Disconnect the active cord from the device handle and from the electrosurgical unit. Clean and store the active cord per the active cord manufacturer's instructions.
9. Retract the snare into the sheath and remove the device from the endoscope.
10. Retrieve the polyp and prepare the specimen per institutional guidelines.

DISPOSAL OF DEVICE

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, the competent authority, and/or regulatory authority of the country where the device was used.

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The symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

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If symbol appears on product label, X = maximum insertion portion width

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If symbol appears on product label, X = minimum accessory channel

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RxOnly

This device is intended for single use only.

cookmedical.com/symbol-glossary

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina
27105
USA

Cook Medical Europe LTD
O'Halloran Road
National Technology Park
Limerick
Ireland

Australian Sponsor

William A. Cook Australia Pty. Ltd.
95 Brandt Street,
Eight Mile Plains QLD 4113
Australia

2025-04
10850_0425
© COOK 2025
cookmedical.com