EN 4 Instinct Plus® Endoscopic Clipping Device

Instructions for Use

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


2.

3.

Instinct Plus® Endoscopic Clipping Device

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Performance Characteristics

The Instinct Plus Endoscopic Clipping Device consists of a single-use clip pre-loaded on a flexible delivery system. The handle of the device allows the clip to be opened, closed, rotated, and deployed at the target site.

Device Compatibility

This device is compatible with a side viewing or end-viewing endoscope.

Refer to the package label for minimum channel size required for this device.

Patient Population

This device is intended for adult use only.

Intended Users

Use of this device is restricted to a trained healthcare professional.

Contact with Body Tissue

Tissue contact with this device is in line with the intended use.

Operating Principle

The clip is closed onto the tissue by applying pressure on the handle spool. To permanently deploy the clip, the handle spool is pulled toward the handle thumb ring until the clip detaches.

Qualitative and Quantitative Information

Table 1

Implant Material Qualitative Information		Quantitative Information Weight (grams)
455 Stainless Steel	Jaws	Up to 0.024
Nitinol	Strip	Up to 0.006
304 Stainless Steel	Spring Wire	Up to 0.035
304 Stainless Steel	Driver	Up to 0.013
303 Stainless Steel	Catheter Attachment	Up to 0.030
301 Stainless Steel, Delrin®	Housing	Up to 0.098

INTENDED USE

This device is endoscopically placed within the gastrointestinal tract to apply mechanical pressure to tissue.

INDICATIONS FOR USE

This device is used for endoscopic clip placement within the gastrointestinal tract for the purpose of:

1. Endoscopic marking,
2. Hemostasis for
 - Mucosal/submucosal defects less than 3 cm,
 - Bleeding ulcers,
 - Arteries less than 2 mm,
 - Polyps less than 1.5 cm in diameter,
 - Diverticula in the colon, and
 - Prophylactic clipping to reduce the risk of delayed bleeding post lesion resection.
3. Anchoring to affix jejunal feeding tubes to the wall of the small bowel,
4. As a supplementary method for closure of GI tract luminal perforations less than 20 mm that can be treated conservatively,
5. Anchoring to affix fully covered esophageal self-expanding metal stents to the wall of the esophagus in patients with fistulas, leaks, perforations, or disunion.

CLINICAL BENEFITS

Table 2

Indication	Intended Clinical Benefit
Endoscopic Marking	• Identify GI lesion prior to subsequent procedure
Hemostasis	• Stop gastrointestinal bleeding • Relief of symptoms associated with GI bleeding
Prophylactic clipping to reduce the risk of delayed bleeding post lesion resection	• Prevent gastrointestinal bleeding • Prevent symptoms associated with GI bleeding
Anchoring to affix jejunal tube to the wall of the small bowel	• Prevent migration of feeding tube
As a supplementary method for closure of GI tract luminal perforations less than 20 mm that can be treated conservatively	• Relief of symptoms associated with GI perforations
Anchoring to affix fully covered esophageal self-expanding metal stents to the wall of the esophagus in patients with fistulas, leaks, perforations, or disunion	• Relief of symptoms associated with esophageal fistulas, leaks, perforations, or disunion • Prevent migration of esophageal stent

WARNINGS

This device has not been evaluated for anchoring to affix fully covered esophageal self-expanding metal stents to the wall of the esophagus in patients with esophageal strictures or malignant obstructions.

This device contains nickel which may cause an allergic reaction and other potential systemic complications related to nickel allergy.

Additional warnings may be listed in procedural steps below. Do not use this device for any purpose other than stated intended use.

Sterile and Single Use

This device is supplied sterile and designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device contamination with biological or chemical agents and/or mechanical integrity failure of device.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Those specific to endoscopic hemostasis include: coagulopathy • cricopharyngeal or esophageal narrowing or stricture • tortuous esophagus • uncooperative patient.

PRECAUTIONS

This device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that this Instinct Plus Endoscopic Clipping Device is MR Conditional. A patient with this device can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 3.0 or 1.5 T
- Maximum spatial field gradient of 1900 Gauss/cm (19T/m)
- Maximum MRI system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg

Under the scan conditions defined above, non-clinical testing results indicated the Instinct Plus Endoscopic Clipping Device is expected to produce a maximum temperature rise of 4°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 82 mm from the Instinct Plus Endoscopic Clipping Device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

PRECAUTION

It is recommended that healthcare providers distribute patient MR information cards with the name of the clip and date it was placed. **Note:** A patient MR information card is

provided with each device. Please complete and provide to patient to ensure safety for future MRI procedures.

WARNING

Failure to follow the recommended MR conditional labeling may result in a deployed GI hemoclip dislodging from tissue or heating of tissue at the GI hemoclip location. A GI hemoclip dislodgement may result in rebleeding requiring additional intervention or surgery, serious injury, or death.

WARNING

This clip contains ferromagnetic material. Follow your institutional protocols to determine whether or not an x-ray should be performed prior to an MRI exam. There may be a small potential risk of clip dislodgment and rebleeding if the clip is used in friable or healing tissues due to magnetic forces acting on the clip when in or near an MRI scanner.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy and endoscopic hemostasis include: allergic reaction to medication • anatomy disruption • aspiration • aspiration pneumonia • cardiac arrhythmia or arrest • fever • hematemesis • hemorrhage • hypotension • infection • migration of clip into the bile duct • minimal acute inflammatory tissue reaction • perforation • respiratory depression or arrest • transient dysphagia • transitory local irritation • wound dehiscence.

HOW SUPPLIED

This device is supplied ethylene oxide (EtO) sterile in a peel open pouch. The pouches are packaged as multiples in a box.

STORAGE REQUIREMENTS

Store in a dark, dry place. Avoid extended exposure to light and extremes of temperature.

INSPECTION OF DEVICE

Do not use the device if the sterile packaging is damaged or unintentionally opened before use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

DEVICE PREPARATION

1. Uncoil the device. Verify smooth handle operation and clip action. Open the clip by gently moving the handle spool distally (away from the handle thumb ring). Once the clip is fully open, do not continue advancing the handle spool as the clip may prematurely detach from the catheter. Close the clip by moving the handle spool proximally until the clip is fully closed (See Fig. 3).

PRECAUTION

Do not continue to pull the handle spool beyond tactile resistance as this may prematurely deploy the clip.

INSTRUCTIONS FOR USE

1. Visually determine the desired tissue site.
2. With the clip closed, hold the red handle cap (See Fig. 1), and advance the device in small increments into the accessory channel of the gastroscope, duodenoscope, or colonoscope. **Note:** If difficult to advance the device, relax the elevator and/or straighten the endoscope.

PRECAUTIONS

- Holding the handle spool or the handle in a way that prevents handle spool mobility during clip advancement may prematurely deploy the clip (See Fig. 1).
- Failure to keep the endoscope as straight as possible when inserting the device can result in a difficult passage.
- Failure to keep the endoscope elevator as relaxed as possible when inserting the device can result in a difficult passage.
3. Endoscopically confirm that the device has exited the endoscope.
 4. Position the distal tip of the device toward the targeted site and open the clip. The clip can be rotated by turning the handle until the clip is in proper position (See Fig. 2).
 5. Advance the device into contact with the targeted site. **Note:** Ensure the clip is not advanced into tissue with excessive force.

PRECAUTION

- Advancing the device into tissue with excessive force may cause the device to lacerate tissue during approximation.
6. When satisfied with the clip position, close the clip onto the tissue by using slight pressure on the handle spool, until tactile resistance is felt. **Note:** Avoid clipping onto hard or severely fibrotic tissue. **Note:** Avoid clipping onto previously deployed clips.

PRECAUTIONS

- Clipping hard or severely fibrotic lesions may result in difficulty approximating tissue or premature clip deployment.
- Clipping onto a previously deployed device may result in difficulty approximating tissue or premature clip deployment.
7. Clip position may now be assessed prior to deployment.
 8. If the clip is not in the desired position, the clip may be reopened and repositioned up to 5 times.

PRECAUTIONS

- Reopening and repositioning the clip greater than 5 times can result in difficulty opening and closing clip.
- Reopening and repositioning the clip greater than 5 times can result in difficulty positioning the device or approximating tissue.
9. To permanently deploy the clip, pull the handle spool toward the handle thumb ring until the clip detaches (See Fig. 3). **Note:** If separation of the clip is not immediate, gently move the catheter back and forth or use other endoscopic maneuvers to separate the catheter from the clip.

10. After clip deployment, continue to apply slight pressure on the handle spool as the device is removed from the endoscope (See Fig. 3). **Note:** If it is difficult to remove the device, relax the elevator and/or straighten the endoscope.

WARNING

Failure to apply slight pressure to the handle during removal can result in tissue damage from the device.

PRECAUTIONS

- Failure to apply slight pressure to the handle during removal can result in difficulty removing the device from the endoscope.
- Failure to keep the endoscope as straight as possible when removing the device can result in difficulty removing the device from the endoscope.
- Failure to keep the endoscope elevator as relaxed as possible when removing the device can result in difficulty removing the device from the endoscope.

DISPOSAL OF DEVICE

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

Summary of Safety and Clinical Performance (SSCP)

When available, you can use the basic unique device identifier for this device (0827002WCM00000000020ME) to find the Summary of Safety and Clinical Performance (SSCP) on the Eudamed website: www.ec.europa.eu/tools/eudamed

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, the competent authority, and/or regulatory authority of the country where the device was used.

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If symbol appears on product label, X = wire guide compatibility

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If symbol appears on product label, X = minimum accessory channel

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If symbol appears on product label, X = minimum accessory channel for side viewing scope

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The symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a cobalt alloy (or stainless steel alloy containing cobalt), which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence

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
Hold here when advancing into accessory channel

CS Při posouvání do aksesorního kanálu držte zde • **DA** Hold her ved fremføring ind i tilbehørskanal • **DE** Beim Verschieben in den Arbeitskanal hier halten • **EL** Κρατήστε εδώ κατά την προώθηση στο κανάλι εργασίας • **ES** Sujete aquí al hacer avanzar el dispositivo al interior del canal de accesorios • **ET** Hoidke siin lisakanalisse liikumisel • **FI** Pidä tästä kiinni, kun työnnät työskentelykanavaan • **FR** Tenir ici lors de la progression dans le canal opérateur • **HR** Držite ovdje prilikom uvođenja u pomoćni kanal • **HU** Tartsa itt, amikor a munkacsatornába tolja • **IT** Afferrare qui durante l'avanzamento nel canale operativo • **NL** Hier vasthouden bij het opvoeren in het werkkanaal • **NO** Hold her når du fører anordningen inn i arbeidskanalen • **PL** Przytrzymać tutaj podczas wsuwania do kanału roboczego • **PT** Segure aqui quando avançar para o canal acessório • **RO** Țineți aici atunci când avansați în canalul pentru accesorii • **SL** Pri potiskanju v kanal za dodatke držite tukaj • **SV** Håll här när du för in anordningen i arbetskannalen

Rx only

cookmedical.com/symbol-glossary

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina
27105
USA

Cook Medical Europe LTD
O'Halloran Road
National Technology Park
Limerick
Ireland

William A. Cook Australia Pty. Ltd.
95 Brandl Street,
Eight Mile Plains
QLD 4113
Australia

2025-10
10894_1025
© COOK 2025
cookmedical.com