COOK®

MEDICAL

- EN**
2
Biliary Wire Guide
Instructions for Use
- ES**
3
Guía biliar
Instrucciones de uso
- FR**
4
Guide biliaire
Mode d'emploi
- PT**
5
Fio guia biliar
Instruções de utilização

* 1 0 9 7 2 / 1 1 2 3 *

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

INTENDED USE

This device is used to assist in cannulation of the biliary and pancreatic ducts and to aid in bridging difficult strictures during ERCP.

NOTES

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL COMPLICATIONS

Potential adverse events associated with ERCP include, but are not limited to: perforation, allergic reaction to contrast or medication.

PRECAUTIONS

Use of this wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide (except THSF wire guides).

WARNINGS

Remove this wire guide before applying electrical current.


This device is designed for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use.

Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

INSTRUCTIONS FOR USE

1. Prior to removing wire guide from holder, flush with 30 cc of sterile water. (See fig.)
2. Flush endoscope accessory channel and/or lumen of device with sterile water, **then insert wire guide floppy end first.** **Note:** For best results, wire guide should be kept wet.

3. Fluoroscopically monitor wire guide advancement in ductal system.
4. Proceed with exchange of compatible wire-guided accessories over wire guide.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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A symbol glossary can be found at cookmedical.com/symbol-glossary

ES En cookmedical.com/symbol-glossary, puede consultarse un glosario de símbolos • **FR** Pour un glossaire des symboles, consulter le site Web cookmedical.com/symbol-glossary • **PT** Existe um glossário de símbolos em cookmedical.com/symbol-glossary

If symbol appears on product label, X = maximum insertion portion width

ES Si el símbolo aparece en la etiqueta del producto, X = anchura máxima de la parte que se va a insertar • **FR** Si le symbole est visible sur l'étiquette du produit, X = largeur maximale de la section d'insertion • **PT** Se o símbolo aparecer no rótulo do produto, X = largura máxima da parte introduzida

If symbol appears on product label, X = minimum accessory channel

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If symbol appears on product label, X = wire guide compatibility

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If symbol appears on product label, X = quantity per box

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RxOnly

**This device is intended
for single use only.**

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina
27105
USA

Australian Sponsor

William A. Cook Australia Pty. Ltd.
95 Brandl Street,
Eight Mile Plains
QLD 4113
Australia

2023-11
cookmedical.com
© COOK 2023