FS-OMNI, FS

EN
4

Fusion® Sphincterotome

Instructions for Use

ES
6

Esfinterótomo Fusion®

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FR
10

Sphinctérotome Fusion®

Mode d'emploi


PT
13

Esfinterótomo Fusion®

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* 1 0 9 7 4 / 1 1 2 3 *

1 Wire Stop
 ES Tope de la guía
 FR Butée du guide
 PT Batente do fio guia

2 Proximal Wire Port
 ES Acceso proximal de la guía
 FR Orifice pour guide proximal
 PT Orifício proximal do fio

3 Retaining Clip
 ES Clip de retención
 FR Pince de retenue
 PT Grampo de fixação

4 Injection Port
 ES Orificio de inyección
 FR Orifice d'injection
 PT Orifício de injeção

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


CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

INTENDED USE

This device is used for cannulation of the ductal system and for sphincterotomy. If preloaded, also aids in bridging difficult strictures during ERCP.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with sphincterotomy.

Contraindications to sphincterotomy include but are not limited to: coagulopathy and inability to properly position the sphincterotome cutting wire.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

USER QUALIFICATIONS

The sphincterotome should only be used by or under the supervision of physicians who have been trained in ERCP procedures.

WARNINGS

Do not use this device for any purpose other than stated intended use.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Any electro-surgical device constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

The electro-surgical effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.

Before using this device, follow recommendations provided by the electro-surgical unit manufacturer to ensure patient safety through proper selection, placement, and utilization

of the patient return electrode. Ensure a proper path from the patient return electrode to the electro-surgical unit is maintained throughout the procedure.

For patients with electrically conductive implants, monopolar diathermy or electro-surgical cautery can result in localized thermal injury at the implant site.

Do not use this device with an output higher than the rated voltage 1.5 kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, sphincterotome and/or active cord.

Do not use this device with an active cord which has a maximum voltage rating less than 1.5 kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant. If cutting is not observed at normal operating settings, this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.

The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.

Skin-to-skin contact should be avoided (for example between the patient's arms and body) by using dry gauze.

Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.

Avoid contact between the active cords and the patient's body, or any other electrodes.

This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electric shock in compliance with IEC 60601-1). Use of generator other than type BF or CF could cause electric shock to patient, operator, or assistant.

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere, or in the presence of explosive gases. Doing so may result in explosion.

In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electro-surgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the sphincterotome in these patients.

This device must be used with an endoscope that is intended for use with high frequency (HF) electro-surgical devices. If the device is used with an endoscope not compatible with HF electro-surgery, a burn to the patient or the operator may occur.

To avoid neuromuscular stimulation and minimize electric arc, activate the cutting wire only when it is in contact with the tissue to be cauterized.

Additional warnings may be listed in procedural steps below.

PRECAUTIONS

Switch the electro-surgical unit to the "off" position when it is not in use.

Do not over flex or bow the tip beyond 90°, as this may damage or cause the cutting wire to break.

The elevator should remain open/down when advancing or retracting the sphincterotome.

If preloaded, use of a wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide.

TECHNICAL DESCRIPTION

See fig. 1 for Fusion IDE Tome and fig. 2 for Fusion OMNI Sphincterotome. Components marked with an asterisk (*) denote applied parts (the part of the device that comes into contact with the patient, as defined by IEC 60601-1), see figs. 1 and 2. Prior to starting the procedure in which the device will be used, become familiarized with distal mark location with respect to the cutting wire. Handle reference markings may be used for approximation of tip deflection. Refer to the package label for device dimensions, minimum endoscope compatibility, and wire guide compatibility.

This device is recommended for use with Cook electro-surgical compatible wire guides.

This device must only be used with an active cord compatible with a 3 mm diameter connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

Maximum rated input voltage for this device is 1.5 kVp-p (750 Vp).

This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6°F).

Store in a dry location. Recommended storage temperature 22°C (71.6°F). Device may be exposed to temperatures up to 50°C (122°F) during transportation.

WARNING

No modification of this equipment is allowed.

SYSTEM PREPARATION

PRECAUTIONS

Do not apply manual pressure to the tip or the cutting wire of the sphincterotome in an attempt to influence orientation, as this may result in damage to the device.

Do not exercise the handle while the device is coiled or the precurved stylet is in place, as this may cause damage to the sphincterotome and render it inoperable.

WARNINGS

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

Before using this device, follow the recommendations

provided by the electro-surgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electro-surgical unit is maintained throughout procedure.

1. With the electro-surgical unit off, prepare the equipment.
2. Attach the Wire Guide Locking Device to the endoscope accessory port.
3. Upon removing the device from the package, uncoil and straighten sphincterotome. Carefully remove the precurved stylet from the cannulating tip.
4. Prior to use, flush the catheter with sterile water or saline.

Note: If using the Fusion OMNI Sphincterotome (see fig. 2). If not preloaded with a wire guide, the Slide Cover should be slid down to cover the Wire Control Port prior to flushing to prevent leakage. If flushing over a wire guide that is preloaded through the Wire Control Port, the Slide Cover should be slid down to partially cover the Wire Control Port prior to flushing to minimize leakage. After flushing over a preloaded wire guide, the Slide Cover should be moved back proximal to the Wire Control Port to enable the wire guide to be moved freely.

INSTRUCTIONS FOR USE

WARNING

The sphincterotome should never be connected to the active cord prior to inserting the sphincterotome through the endoscope. Doing so may injure the patient and/or damage the equipment due to improper electrical circuit grounding.

1. Insert the catheter into the endoscope.

Note: If using the Intra Ductal Exchange (IDE) port and short wire (See fig. 1) remove the Wire Stop from its retaining clip. Ensure that the Wire Stop is disengaged from the Proximal Wire Port (PWP) hub and that the distal tip of the Wire Stop is proximal to the IDE port. Insert the distal tip of the wire guide into the IDE port and advance until it is flush with the distal tip of the sphincterotome. Advance the Wire Stop handle until it reaches the PWP hub and connect the Luer lock securely to the PWP hub. Advance the tip of the sphincterotome through the cap of the Wire Guide Locking Device and continue advancing until it is endoscopically visible. Disengage the Wire Stop by releasing the handle from the PWP hub and remove the Wire Stop from the catheter.

Note: If using the Proximal Wire Port (PWP) and a pre-positioned long wire guide (See fig. 1) remove the Wire Stop (if applicable). Advance the sphincterotome over the pre-positioned wire guide ensuring that the wire exits the catheter at the PWP. Continue advancing the device until it is endoscopically visible. For best results, wire guide should be kept wet with sterile water or saline.

Note: If using the Wire Control Port (WCP) and a short wire or long wire guide (See fig. 2) Advance the tip of the sphincterotome through the cap of the Wire Guide Locking

Device and continue advancing until it is endoscopically visible. If not preloaded, a wire guide may be introduced into the WCP and separated from the catheter (See Figs. 3 and 4).

CANNULATION

WARNING

Failure to confirm the cutting wire has fully exited endoscope may result in contact between the cutting wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

1. Using endoscopic visualization, confirm the cutting wire has exited the endoscope.
2. Cannulate using the sphincterotome. **Note:** A wire guide may be utilized to facilitate cannulation.
3. Following cannulation, confirm position fluoroscopically via contrast injection through the injection port or visualization of the radiopaque marker.

SPHINCTEROTOMY

WARNING

To avoid unintended disconnection during procedure and damage to the device, device should only be used with compatible active cords (Reference Technical Description for compatibility).

1. Attach the active cord to the connector on the sphincterotome handle. Active cord fittings should fit snugly into both the device handle and the electrosurgical unit.
2. If using a wire guide that is not compatible with electrosurgical procedures, remove the wire guide from the sphincterotome before activating the electrosurgical current. If the wire is not removed, there is a risk of thermal injury to the patient or operator.

WARNINGS

Do not use a wire guide (or sphincterotome) that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

The sphincterotome should not be energized prior to performing sphincterotomy. Energizing the cutting wire prior to use will cause premature cutting wire fatigue and will compromise the cutting wire's integrity.

Verify that the cutting wire has exited the endoscope by visualizing it on the endoscope monitor. Failure to do so may result in contact between the cutting wire and the endoscope while electrical current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire and/or damage to the endoscope.

User must maintain constant contact with tissue when applying electrocautery current. Failure to do so may result in increased current resulting in a broken cutting wire, damage to the endoscope and/or patient injury.

When applying current, ensure cutting wire is completely out of endoscope. Contact of cutting wire with endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope. If a wire guide that is not compatible with electrosurgical procedures is used in sphincterotome, it must be removed prior to applying electrosurgical current.

Cook wire guides packaged with this device may be left in place during sphincterotomy.

3. Following the electrosurgical unit manufacturer's instructions, verify desired settings and proceed with the sphincterotomy.
4. Upon completion of the sphincterotomy, turn the electrosurgical unit off.
5. Disconnect the active cord from the device handle and from the electrosurgical unit. Clean and store the active cord per the active cord manufacturer's instructions.

DEVICE REMOVAL FROM ENDOSCOPE

1. A previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices. If the wire guide is to remain in place while the device is withdrawn, utilize the following steps:

Note: If using the Intra Ductal Exchange (IDE) port and short wire (See fig. 1), prior to removing the device, utilize the reference marks on the catheter to ensure that the IDE port is within the ductal system. Fluoroscopically visualize the radiopaque band at the IDE port. Retract the wire guide until the radiopaque distal tip of the wire guide passes the band. A disengagement from the wire guide lumen will occur. Advance the disengaged wire guide to maintain ductal access. Lock the wire guide into the Wire Guide Locking Device and remove the sphincterotome from the endoscope accessory channel.

Note: If using the Proximal Wire Port (PWP) and a pre-positioned long wire guide (See fig. 1) remove the device using standard long wire exchange technique.

Note: If using the Wire Control Port (WCP) and a short wire or long wire guide. To withdraw the sphincterotome from the endoscope, ensure the wire guide is locked in place, and then pull back on the catheter, allowing the wire guide lumen to separate from the wire guide until the metal band is visible at the Wire Guide Locking Device and resistance is felt. (See fig. 5) Unlock the wire guide from the Wire Guide Locking Device, completely remove the sphincterotome from the wire guide, and re-lock the wire guide.

DEVICE DISPOSAL

1. Upon completion of the procedure, dispose of the device(s) per institutional guidelines for biohazardous medical waste.

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A symbol glossary can be found at cookmedical.com/symbol-glossary

ES En cookmedical.com/symbol-glossary, puede consultarse un glosario de símbolos • FR Pour un glossaire des symboles, consulter le site Web cookmedical.com/symbol-glossary • PT Existe um glossário de símbolos em cookmedical.com/symbol-glossary

If symbol appears on product label, X = maximum insertion portion width

ES Si el símbolo aparece en la etiqueta del producto, X = anchura máxima de la parte que se va a insertar • FR Si le symbole est visible sur l'étiquette du produit, X = largeur maximale de la section d'insertion • PT Se o símbolo aparecer no rótulo do produto, X = largura máxima da parte introduzida

If symbol appears on product label, X = minimum accessory channel

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If symbol appears on product label, X = wire guide compatibility

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If symbol appears on product label, X = quantity per box

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This device is intended
for single use only.

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina
27105
USA

Australian Sponsor

William A. Cook Australia Pty. Ltd.
95 Brandl Street,
Eight Mile Plains
QLD 4113
Australia

2023-11
cookmedical.com
© COOK 2023