
EN
2

SureTome™ SW
Sphincterotome with DomeTip®
Instructions for Use

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SureTome™ SW

Sphincterotome with DomeTip®

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Performance Characteristics

Sphincterotomes consist of a three-ring handle with a metal hot pin for connection to an electrosurgical generator via an active cord. The hot pin delivers current to the cutting wire that is exposed at the distal end of the device. A multi-lumen catheter allows for injection of contrast or passage over a prepositioned wire guide using ports on the proximal end of the device. This device is optimized for use with short-wire (SW) technique.

Maximum rated input voltage for this device is 1.5 kVp-p (750 Vp).

Device Compatibility

Refer to the package label for device dimensions and minimum endoscope channel size.

This device is recommended for use with electrosurgical compatible wire guides. Refer to package label for compatible wire guide size.

This device is compatible with Cook wire guide locking devices to lock the wire guide(s) in place.

This device must only be used with an active cord compatible with a 3 mm diameter connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

Electrosurgical generator - type BF or CF (generators that provide protection from electric shock in compliance with IEC60601-1)

The active length information provided is intended to be used with the guidance from the Electrosurgical Generator manufacturer to determine if the overall length of an Electrosurgical Accessory remains within the maximum allowable length. The overall active length of this Electrosurgical Accessory is 225 cm. The overall active length includes the cutting wire, catheter, and handle. The length of other connecting cords must be considered when used with the Electrosurgical Accessory.

Patient Population

This device is indicated for adult use only.

Intended User

The sphincterotome should only be used by or under the supervision of physicians who have been trained in ERCP procedures.

Contact with Body Tissue

This device has contact with body tissue in line with the intended use.

Operating Principle

In a sphincterotomy, the sphincterotome is passed through a duodenoscope and positioned adjacent to the major papilla. The papilla is cannulated using the tip of the sphincterotome. The papilla is then incised using electrosurgical current applied to the cutting wire by an electrosurgical generator.

INTENDED USE

This device is used for cannulation of the ductal system and for sphincterotomy. If preloaded, also aids in bridging difficult strictures during ERCP.

INDICATIONS FOR USE

Sphincterotomes are used for selective deep cannulation of the ductal system and sphincterotomy to facilitate the completion of diagnostic and therapeutic endoscopic retrograde cholangiopancreatography.

CLINICAL BENEFITS

Sphincterotomes facilitate the completion of diagnostic procedures and interventions during endoscopic retrograde cholangiopancreatography (ERCP).

WARNINGS

- Do not use this device for any purpose other than stated intended use.
- Any electrosurgical device constitutes a potential electrical hazard to patient and operator. Possible adverse effects include: burns • cardiac arrhythmia • fulguration • nerve and/or muscle stimulation.
- The electrosurgical effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.
- Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper selection, placement, and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.
- For patients with electrically conductive implants, monopolar diathermy, or electrosurgical cautery can result in localized thermal injury at the implant site.
- Do not use this device with an output higher than the rated voltage 1.5 kVp-p (750 Vp). This could cause thermal injury

to the patient, operator, or assistant and could also damage the endoscope, sphincterotome, and/or active cord.

- Do not use this device with an active cord which has a maximum voltage rating less than 1.5 kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or assistant.
- This device must be used with an endoscope compatible with high frequency equipment.
- If cutting is not observed at normal operating settings, this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.
- The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.
- Skin-to-skin contact should be avoided (for example between the patient's arms and body) by using dry gauze.
- Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.
- Avoid contact between the active cords and the patient's body, or any other electrodes.
- This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electrical shock in compliance with IEC60601-1). Use of generator other than type BF or CF could cause electric shock to patient, operator, or assistant.
- For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.
- This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere, or in the presence of explosive gases. Doing so may result in explosion.
- In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electrosurgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the sphincterotome in these patients.
- This device must be used with an endoscope that is intended for use with high frequency (HF) electrosurgical devices. If the device is used with an endoscope not compatible with HF electrosurgery, a burn to the patient or the operator may occur.
- To avoid neuromuscular stimulation and minimize electric arc, activate the cutting wire only when it is in contact with the tissue to be cauterized.
- This device contains nickel which may cause an allergic reaction and other potential systemic complications related to nickel allergy.
- Attempts to cannulate through a metal stent may result in damage to the device requiring endoscopic retrieval.

- Cutting wire orientation and the direction/length of the cut should be assessed before applying current and continually monitored during delivery of current.
- Additional warnings may be listed in procedural steps below.

Sterile and Single Use

This device is supplied sterile and designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device contamination with biological or chemical agents, electrical failure, and/or mechanical integrity failure of device.

CONTRAINDICATIONS

Contraindications of biliary sphincterotomy include general contraindications of ERCP, i.e., unstable or uncooperative patients. Contraindications specific to biliary sphincterotomy include: the inability to orient the cutting wire of the sphincterotome towards the axis of the CBD • the patient's inability to provide informed consent • newly created gastrointestinal anastomosis • uncorrected coagulopathy.

PRECAUTIONS

Switch the electrosurgical unit to the "off" position when it is not in use.

Do not over flex or bow the tip beyond 90°, as this may damage or cause the cutting wire to break.

The elevator should remain open/down when advancing or retracting the sphincterotome.

If preloaded, use of a wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide.

TECHNICAL DESCRIPTION

Components marked with an asterisk (*) denote applied parts (the part of the device that comes into contact with the patient, as defined by IEC60601-1), see Fig. 1. Prior to starting the procedure in which the device will be used, become familiarized with distal mark location with respect to the cutting wire. Handle reference markings may be used for approximation of tip deflection.

This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6°F).

WARNING

No modification of this equipment is allowed.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include: allergic reaction to contrast or medication • aspiration • cardiac arrhythmia or arrest • cholangitis • hemorrhage • hypotension • infection • pancreatitis • perforation • respiratory depression or arrest • sepsis.

Potential Complications associated with ductal cannulation and bridging of strictures include: bleeding • pancreatitis • perforation • tissue inflammation.

HOW SUPPLIED

This device is supplied ethylene oxide (EtO) sterile in a peel-open pouch.

STORAGE

Store in a dark, dry place. Recommended storage temperature 22°C (71.6°F). Device may be exposed to temperatures up to 50°C (122°F) during transportation.

INSPECTION OF DEVICE

Do not use the device if the sterile packaging is damaged or unintentionally opened before use. Visually inspect with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

DEVICE PREPARATION

PRECAUTIONS

Do not apply manual pressure to the tip or the cutting wire of the sphincterotome in an attempt to influence orientation, as this may result in damage to the device.

Do not exercise the handle while the device is coiled or the precurved stylet is in place, as this may cause damage to the sphincterotome and render it inoperable.

1. Upon removing the device from the package, uncoil and straighten the sphincterotome. Carefully remove the precurved stylet wire from the cannulating tip.

WARNINGS

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

Before using this device, follow the recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout procedure.

2. With the electrosurgical unit off, prepare the equipment.
3. Attach the wire guide locking device to the endoscope accessory port.

INSTRUCTIONS FOR USE

WARNING

The sphincterotome should never be connected to the active cord prior to inserting the sphincterotome through the endoscope. Doing so may injure the patient and/or damage the equipment due to improper electrical circuit grounding.

1. Insert the catheter into the endoscope.

Note: If introducing a wire guide, using a preloaded sphincterotome, advance the tip of the sphincterotome through the cap of the wire guide locking device and continue advancing until it is endoscopically visible. If not preloaded, a wire guide may be introduced into the Wire Port (see Fig.1).

Note: If using a pre-positioned wire guide (.025" and .035" wires recommended), advance the sphincterotome over the pre-positioned wire guide. Slightly bend catheter at the 5 cm mark (see Fig. 1) and advance proximal end of wire through bend. Continue advancing the device until it is endoscopically visible. For best results, wire guide should be kept wet with sterile water or saline.

CANNULATION

WARNING

Failure to confirm the cutting wire has fully exited endoscope may result in contact between the cutting wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

1. Using endoscopic visualization, confirm the cutting wire has exited the endoscope.

Note: Optimize the tip orientation for cannulation by grasping the proximal end of the catheter marked 'Rotate' and rotating the catheter handle.

2. Separate the wire guide lumen from the wire guide up to the wire locking device to allow for physician wire control.
3. Cannulate using the sphincterotome. **Note:** A wire guide may be utilized to facilitate cannulation.
4. Following cannulation, inject contrast through the injection port to fluoroscopically confirm position.

SPHINCTEROTOMY

WARNING

To avoid unintended disconnection during the procedure and damage to the device, the device should only be used with compatible active cords (Reference Technical Description for compatibility).

1. Attach the active cord to the connector on the sphincterotome handle. Active cord fittings should fit snugly into both the device handle and the electrosurgical unit.
2. If using a wire guide that is not compatible with electrosurgical procedures, remove the wire guide from the sphincterotome before activating the electrosurgical current. If the wire is not removed, there is a risk of thermal injury to the patient or operator.

WARNINGS

Do not use a wire guide (or sphincterotome) that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

The sphincterotome should not be energized prior to performing sphincterotomy. Energizing the cutting wire prior to use will cause premature cutting wire fatigue and will compromise the cutting wire's integrity.

Verify that the cutting wire, including the Protector™ insulated sheath, has fully exited the endoscope by visualizing the proximal end of the cutting wire on the

endoscope monitor. Failure to do so may result in contact between the cutting wire and the endoscope while electrical current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

User must maintain constant contact with tissue when applying electrocautery current. Failure to do so may result in increased current resulting in a broken cutting wire, damage to the endoscope and/or patient injury.

When applying current, ensure cutting wire is completely out of endoscope. Contact of cutting wire with endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope.

If a wire guide that is not compatible with electrosurgical procedures is used in sphincterotome, it must be removed prior to applying electrosurgical current.

Wire guides packaged with this device may be left in place during sphincterotomy.

3. Following the electrosurgical unit manufacturer's instructions, verify desired settings and proceed with the sphincterotomy.
4. Upon completion of the sphincterotomy, turn the electrosurgical unit off.
5. Disconnect the active cord from the device handle.
6. Device may be rotated to selectively cannulate as needed.

DEVICE REMOVAL FROM ENDOSCOPE

A wire guide may be left in position to facilitate introduction of other wire-guided devices.

1. To withdraw the sphincterotome from the endoscope, ensure the wire guide is locked in place, and then pull back on the catheter, allowing the wire guide lumen to separate from the wire guide to the 5 cm mark (see Fig.1).
2. To complete removal from the endoscope do one of the following:
 - Completely separate the wire guide through the tip, or
 - unlock the wire guide from the wire guide locking device, completely remove the sphincterotome from the wire guide, and re-lock the wire guide.

DISPOSAL OF DEVICE

After completion of the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contradictions, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, the competent authority, and/or regulatory authority of the country where the device was used.

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© COOK 2025
Wilson-Cook Medical, Inc.
William A. Cook Australia Pty. Ltd.