


EN	Instinct™ Endoscopic Hemoclip
CS	Endoskopická svorka Instinct™ Hemoclip
DA	Instinct™ endoskopisk hæmostaseclips
NL	Instinct™ endoscopische hemostatische clip
FR	Hémoclip endoscopique Instinct™
DE	Endoskopischer Hämoclip Instinct™
EL	Ενδοσκοπικό κλιπ απολίνωσης αγγείων Instinct™
HU	Instinct™ endoszkópos vérzéscsillapító kapocs
IT	Hemoclip endoscopica Instinct™
NO	Instinct™ endoskopisk hemoklips
PL	Endoskopowy klips hemostatyczny Instinct™
PT	Clipe hemostático para endoscopia Instinct™
ES	Clip hemostático endoscópico Instinct™
SV	Instinct™ endoskopiskt hemoclips

1.

2.

3.

ENGLISH

INTENDED USE

This device is used for endoscopic clip placement within the gastrointestinal tract for the purpose of endoscopic marking, hemostasis for mucosal/submucosal defects less than 3 cm in the upper GI tract, bleeding ulcers, arteries less than 2 mm, and polyps less than 1.5 cm in diameter in the GI tract. This device is not intended for the repair of GI tract luminal perforations.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

The number of clips required for hemostasis may vary depending upon anatomical site, histology, lesion type, and patient condition and history. A sufficient quantity of clips should be prepared in consideration of all of these factors prior to procedure.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Those specific to endoscopic hemostasis include, but are not limited to: uncooperative patient, coagulopathy, cricopharyngeal or esophageal narrowing or stricture, and tortuous esophagus.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy and endoscopic hemostasis include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, hematemesis, transient dysphagia, aspiration pneumonia, wound dehiscence, minimal acute inflammatory tissue reaction, transitory local irritation, migration of clip into the bile duct, and anatomy disruption.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Endoscope must remain as straight as possible when inserting or withdrawing device.

Clip must remain closed during introduction into, advancement through and removal from endoscope. If clip is open, damage to clip and endoscope may occur.

Exercising handle while clip is coiled may result in damage to clip.

Endoscopic hemostasis should be performed under direct endoscopic visualization.

Clipping hard or severely fibrotic lesions to achieve hemostasis may be more difficult.

If clip deployment device is used with endoscope in a torqued or retroflexed position, clip deployment difficulties can occur.

Removing undeployed device through a retroflexed scope can cause detachment of the clip.

Lesions located in esophagus and lesser curvature of stomach may be difficult to treat with a forward-viewing endoscope.

Treatment of esophageal varices may require clipping in combination with a sclerosing agent.

Use of clips in presence of bacterial contamination may potentiate or prolong infection.

MRI INFORMATION

MR Conditional

Non-clinical testing has demonstrated that this Instinct Endoscopic Hemoclip is **MR Conditional** according to ASTM F2503. A patient with this clip can be scanned safely immediately after placement under the following conditions.

- Static magnetic field of 3.0 Tesla or less
- Maximum spatial magnetic gradient field of 1600 Gauss/cm or less
- Normal operating mode: Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning or less (i.e., per scanning sequence)

STATIC MAGNETIC FIELD

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-RELATED HEATING

In non-clinical testing, the Instinct Endoscopic Hemoclip produced maximum temperature rise of 1.9°C during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a 3.0 Tesla MR System (Excite, HDx, GE Healthcare, Milwaukee, WI, Software 14X.M5) at an MR system reported whole-body-averaged SAR of 2.9 W/kg (associated with a calorimetry measured whole body averaged value of 2.7 W/kg).

IMAGE ARTIFACTS

MR image quality may be compromised if the area of interest is within approximately 30 mm of the position of the Instinct Endoscopic Hemoclip as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, HDx, GE Healthcare, Milwaukee, WI, Software 14X.M5). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic clip.

There is a potential risk of rebleeding of small arteries due to torque placed on the Instinct Endoscopic Hemoclip during MRI.

FOR US PATIENTS ONLY

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International, 2323 Colorado Avenue, Turlock, CA 95382

Phone: 888-633-4298 (toll free) / 209-668-3333 from outside the US

Fax: 209-669-2450 **Web:** www.medicalert.org

WARNINGS

Failure to deploy clip or improper deployment of a clip may require surgical intervention. Wire cutters may be used to cut wire and coil spring.

Electrosurgical devices used after clipping may cause patient injury.

SYSTEM PREPARATION

1. Uncoil device. Verify smooth handle operation and clip action. Open clip by gently moving handle spool distally (away from handle thumb ring). Once clip is fully open, do not continue advancing handle spool as clip may prematurely detach from catheter. Close clip by moving handle spool proximally until clip is fully closed. **Caution:** Do not continue to pull handle spool beyond tactile resistance as this may prematurely deploy clip.

2. Ensure clip is in closed position.

INSTRUCTIONS FOR USE

1. Visually determine desired tissue site.
2. With clip closed and without holding handle spool, advance device in small increments into accessory channel of gastroscope or colonoscope. **Caution:** Holding handle spool during clip advancement may prematurely deploy clip (See Fig. 1).
3. Endoscopically confirm that device has exited endoscope.
4. Position distal tip of device toward targeted site.
5. Open clip and advance device into contact with the targeted site. Clip can be rotated by turning handle until clip is in proper position. **Note:** Rotation may be limited by clinical circumstances and patient anatomy, among other factors (See Fig. 2).
6. When satisfied with clip position, close clip onto tissue by using slight pressure on handle spool, until tactile resistance is felt. **Note:** Ensure clip is pressed firmly against target site for maximum tissue capture. **Caution:** Do not continue to pull handle spool beyond tactile resistance until ready to deploy clip; otherwise, clip may not reopen.
7. Clip position may now be assessed prior to deployment.
8. If clip is not in desired position, clip may be reopened and repositioned up to 5 times. **Note:** Reopening and closing capability may be limited by clinical circumstances and patient anatomy, among other factors.
9. To permanently deploy clip, pull handle spool toward handle thumb ring until clip detaches. **Note:** If separation of clip is not immediate, gently move catheter back and forth or use other endoscopic maneuvers to separate catheter from clip. **Caution:** Failure to completely release clip may cause clip to pull away from closure site.
10. After clip deployment, continue to apply slight pressure on handle spool as device is removed from endoscope (See Fig. 3).

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

Note: A patient information card is provided with each device; please complete and provide to patient to ensure safety for future MR procedures.

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.