






EN 4	Hemospray® Endoscopic Hemostat Instructions for Use
CS 6	Endoskopické antihemoragikum Hemospray® Návod k použití
DA 9	Hemospray® endoskopisk hæmostat Brugsanvisning
DE 12	Hemospray® endoskopisches Hämostat Gebrauchsanweisung
EL 15	Ενδοσκοπική συσκευή αιμόστασης Hemospray® Οδηγίες χρήσης
ES 18	Dispositivo hemostático endoscópico Hemospray® Instrucciones de uso
ET 21	Endoskoopiline hemostaat Hemospray® Kasutusjuhend
FI 24	Endoskooppinen Hemospray®-hemostaatti Käyttöohjeet
FR 26	Dispositif d'hémostase endoscopique Hemospray® Mode d'emploi
HR 29	Hemospray® endoskopski hemostat Upute za uporabu
HU 32	Hemospray® endoszkópos vérzéscsillapító Használati utasítás
IT 35	Dispositivo endoscopico per emostasi Hemospray® Istruzioni per l'uso
NL 38	Hemospray® endoscopisch hemostaticum Gebruiksaanwijzing
NO 41	Hemospray® endoskopisk hemostatikum Bruksanvisning
PL 44	Endoskopowy przyrząd hemostatyczny Hemospray® Instrukcja używania
PT 47	Dispositivo hemostático endoscópico Hemospray® Instruções de utilização
RO 50	Hemostat endoscopic Hemospray® Instrucțiuni de utilizare
SL 53	Hemospray® endoskopski hemostat Navodila za uporabo
SV 55	Hemospray® endoskopisk blodstillare Bruksanvisning

* 1 2 8 9 8 / 0 3 2 5 *

Hemospray® Endoscopic Hemostat

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Construction of Device

The Hemospray introducer consists of a handle, powder chamber and delivery catheter. The handle houses a CO₂ cartridge, a pressure regulator, and valves to control movement of gas which is controlled definitively by the on/off switch on top of the distal end of the handle. The powder chamber is contained within the handle and the handle is joined to hold the chamber in place. This chamber houses a minimum of 20 g of Hemospray powder. The powder delivery catheter connects to the distal, front side of the handle. Two catheters (220 cm) are provided with each Hemospray assembly.

Performance Characteristics

Hemospray is an inert, bentonite powder developed for endoscopic hemostasis. The powder is delivered by use of a carbon dioxide powered delivery system and through a catheter inserted through the working channel of an endoscope which provides access to the site of the bleed. Each device contains a minimum of 20 g of 100% Bentonite powder.

Device Compatibility

Refer to the package label for the minimum channel size required for this device.

Patient Population

This device is indicated for adult use only.

Intended User

This device is restricted to a trained healthcare professional.

Contact with Body Tissue

Tissue contact for this device is in line with the intended use.

Operating Principle

An endoscope is introduced to the body via one of the body's natural orifices, such as the mouth or anus. After the endoscope has reached the bleeding site, the delivery catheter is advanced through the endoscope and then the handle is attached to the proximal end of the catheter. The CO₂ cartridge in the delivery system is activated by twisting the red activation knob on the handle. The user

presses the deploy button on the handle to release CO₂ from the cartridge. The CO₂ pressure drives the powder from the powder chamber through the catheter. The powder is then locally dispersed to the bleeding site within the GI Tract.

INTENDED USE

This device is used for hemostasis of nonvariceal gastrointestinal bleeding.

INDICATIONS FOR USE

The Hemospray Endoscopic Hemostat is used in adult patients for the hemostasis of nonvariceal gastrointestinal bleeding due to endoscopic procedures and medical conditions such as inflammatory lesions, GI vascular abnormalities, tumors, and ulcers.

CLINICAL BENEFITS

To stop nonvariceal gastrointestinal bleeding.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired target site. Also contraindicated in patients who have gastrointestinal fistulas, are suspected of having a gastrointestinal perforation, or are at high risk of gastrointestinal perforation during endoscopic treatment.

WARNINGS

Keep catheter tip at least 1 cm away from bleeding site to minimize risk of embolization.

Do not use this device for any purpose other than stated intended use.

No modification to this equipment is allowed.

Hemospray has not been approved for use in pediatric populations, no safety or effectiveness data exists and would be considered off-label usage, to do so is at the professional risk of the surgeon/healthcare professional.

Patients with gastrointestinal bleeding that are on antithrombotic medication may be at an increased risk of rebleeding. Follow the relevant clinical guidelines for management of antithrombotic agents for endoscopic procedures.

Sterile and Single Use

This device is supplied sterile and designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device contamination with biological or chemical agents and/or mechanical integrity failure of device.

PRECAUTIONS

For best results, Hemospray must be delivered to the source of bleeding.

Ensure gastrointestinal lumen is not distended because Hemospray adds volume in excess of insufflation volumes during procedure. Closely monitor bowel distension and balance insufflation and Hemospray volumes as necessary.

Use of more than (3) Hemospray devices per patient may result in impaction in colon and is not recommended.

Product contains 16 g CO₂ cartridge. Contents under pressure. Do not puncture or heat above 120° F/ 49° C. Do not inhale or discharge towards face/body. Keep CO₂ dispenser and cartridge out of reach of children.

Hemospray is inert and non-toxic.

As a granular material, unintentional exposure to the powder may cause potential irritation to the skin, eyes and lungs.

In the event of unintended exposure to the powder refer to the following First Aid measures:

Skin: Wash with soap and water until clean.

Eyes: Flush with water until irritation ceases.

Inhalation: Move to area free from powder. If symptoms of irritation persist, contact physician. Inhalation may aggravate existing respiratory illness.

All endoscopic hemostatic therapies, including Hemospray, have an associated risk of rebleeding, particularly in situations where the cause of bleeding is an unresolved underlying disease. After hemostasis has been achieved, monitor patients for rebleeding per the relevant Clinical Guidelines.

Although not seen in clinical practice, there remains a theoretical risk of aspiration of Hemospray powder resulting in respiratory complications. It is prudent to restrict the use of Hemospray to 5 cm below the UES.

This device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

MRI SAFETY INFORMATION

MR The bentonite powder used in Hemospray Endoscopic Hemostat is MR Safe.

PROCEDURAL PRECAUTIONS

Like other modalities, Hemospray may not be effective for all types of bleeds. Gastrointestinal bleeding may exacerbate existing comorbidities, increasing the potential for adverse events including patient mortality.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy including: allergic reaction to medication • aspiration • cardiac arrhythmia or arrest • fever • hemorrhage • hypotension • infection • perforation • respiratory depression or arrest.

Use of Hemospray in the presence of bowel obstruction and/or an anastomosis may pose a risk of injury due to over-distention.

Hemospray may occlude ducts and orifices which communicate with the main bowel lumen. Use caution when using Hemospray in the vicinity of these orifices.

Others including: powder impaction in colon or embolization.

When spraying in retroflexed position, Hemospray powder may adhere to the outside of the endoscope. This may result in difficulty repositioning/removing the endoscope, particularly if passing through a strictured area.

CLINICAL STUDIES

Clinical data summary information that was, in part, the basis for granting the de novo can be found on the Cook Medical website at CookMedical.com/HemosprayData

HOW SUPPLIED

The device is supplied Gamma sterilized in a tray with a peel open top. This device is intended for single use only.

STORAGE

Store in a dark, dry location, away from temperature extremes.

INSPECTION OF DEVICE

Do not use the device if the sterile packaging is damaged or unintentionally opened before use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working conditions, do not use. Please notify Cook for return authorization.

DEVICE PREPARATION

1. Remove device from package.

NOTES

Do not over rotate knob as this could damage the device (See Fig. 1).

Do not test device prior to insertion into endoscope accessory channel as this may increase risk of catheter occlusion.

2. Activate CO₂ cartridge by turning red activation knob until it stops.

WARNING

Do not connect this device to alternate CO₂ sources.

INSTRUCTIONS FOR USE

CAUTION

Ensure gastrointestinal lumen is not distended because Hemospray adds volume during procedure.

1. Before inserting catheter into accessory channel, identify bleeding site, remove as much blood as possible, then flush accessory channel with air.

PRECAUTION

To avoid catheter occlusion, **do not** place catheter directly in contact with blood and/or mucosa, including any pooled blood and **do not** aspirate blood while catheter is in accessory channel.

NOTE

To limit fluid from entering the working channel of the catheter, temporarily occlude the proximal end of the catheter by placing a thumb over the red catheter hub, while advancing the catheter down the accessory channel (See Fig. 2).

2. Slowly advance catheter through accessory channel in short increments until catheter tip is visualized endoscopically (See Fig. 3).
3. To ensure proper visibility, catheter tip should be 1-2 cm away from bleeding site at all times (See Fig. 4).
4. Attach the catheter to the handle, ensuring the connection is secure (See Fig. 5).

NOTE

Device is now active and ready for use. **Do not** press trigger button until powder deployment is desired.

5. To allow powder deployment, turn red valve to open position such that the valve is parallel with the direction of spray (See Fig. 6).

PRECAUTION

The first press of the trigger may not deploy an adequate amount of powder to achieve hemostasis; press again to continue the movement of powder through the system.

NOTES

Allow time for hemostasis before aspirating. Hemostasis may not be immediate because the powder takes time to absorb. Multiple applications may be required; however no more than (3) Hemospray devices should be applied per patient. Irrigation of treated site is not required after deployment. If irrigation is desired, allow time for hemostasis before irrigating with water.

To avoid potential endoscope occlusion, do not aspirate powder into endoscope channel. Scope suction can be turned off or disconnected temporarily to avoid accidental aspiration of powder into the endoscope channel, which may cause endoscope occlusion.

6. To deploy powder, hold handle upright and depress red trigger button for 1-2 seconds and release (See Fig. 7). Continue applying in short bursts until site is completely covered with powder and no active bleeding is visualized.

NOTE

If catheter becomes occluded, turn red valve to closed position, remove catheter from endoscope and replace with extra catheter provided in package. Then repeat steps 1-7 of "Instructions for Use". If more powder is needed and spray canister and/or CO₂ cartridge is empty, a new device is required.

7. Prior to removing Hemospray device from patient, turn red valve to closed position (See Fig. 8). To remove device from endoscope, grasp catheter and pull slowly.

DISPOSAL OF DEVICE

Upon completion of procedure, depressurize device by rotating activation knob counterclockwise until CO₂ cartridge depressurizes completely.

This device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contradiction, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, the competent authority, and/or regulatory authority of the country where the device was used.

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If symbol appears on product label, X = maximum insertion portion width

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If symbol appears on product label, X = minimum accessory channel

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
The symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class 1B). However, the device contains cobalt as a stainless steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

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RxOnly

cookmedical.com/symbol-glossary

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina
27105
USA

Cook Medical Europe LTD
O'Halloran Road
National Technology Park
Limerick
Ireland

William A. Cook Australia Pty. Ltd.
95 Brandl Street,
Eight Mile Plains
QLD 4113
Australia

2025-03
12898_0325
© COOK 2025
cookmedical.com