
EN 3	Universal Active Cord Instructions for Use
DA 5	Universel strømførende ledning Brugsanvisning
DE 8	Universelles HF-Kabel Gebrauchsanweisung
ES 11	Cable de conexión universal Instrucciones de uso
FI 14	Yleiskäyttöinen aktivointijohto Käyttöohjeet
FR 17	Cordon actif universel Mode d'emploi
IT 20	Cavo diatermico universale Istruzioni per l'uso
NL 24	Universeel actief snoer Gebruiksaanwijzing
NO 27	Universell strømførende ledning Bruksanvisning
PL 30	Uniwersalny przewód zasilania Instrukcja użycia
PT 33	Cabo ativo universal Instruções de utilização
SV 36	Strömförande universalkabel Bruksanvisning

* 1 4 5 3 5 / 0 2 2 5 *


1 Monopolar Device Plug

DA Stik til monopolarøst udstyr • **DE** Stecker für monopolar Instrument • **ES** Enchufe monopolar para dispositivo • **FI** Monopolaarisen laitteen pistoke • **FR** Fiche de raccordement au dispositif monopolaire • **IT** Spinotto per dispositivi monopolari • **NL** Monopolaire apparaatstekker • **NO** Plugg for monopolar anordning • **PL** Wtyczka do urządzenia jednobiegunowego • **PT** Ficha para o dispositivo monopolar • **SV** Plugg till monopolar produkt

2 Generator Plug

DA Generatorstik • **DE** Generatorstecker • **ES** Enchufe del generador • **FI** Generaattorin pistoke • **FR** Fiche du générateur • **IT** Spinotto per generatore • **NL** Generatorstekker • **NO** Generatorplugg • **PL** Wtyczka generatora • **PT** Ficha do gerador • **SV** Generatorplugg

1 Monopolar Device Plug

DA Stik til monopolarøst udstyr • **DE** Stecker für monopolar Instrument • **ES** Enchufe monopolar para dispositivo • **FI** Monopolaarisen laitteen pistoke • **FR** Fiche de raccordement au dispositif monopolaire • **IT** Spinotto per dispositivi monopolari • **NL** Monopolaire apparaatstekker • **NO** Plugg for monopolar anordning • **PL** Wtyczka do urządzenia jednobiegunowego • **PT** Ficha para o dispositivo monopolar • **SV** Plugg till monopolar produkt

2 Generator Plug

DA Generatorstik • **DE** Generatorstecker • **ES** Enchufe del generador • **FI** Generaattorin pistoke • **FR** Fiche du générateur • **IT** Spinotto per generatore • **NL** Generatorstekker • **NO** Generatorplugg • **PL** Wtyczka generatora • **PT** Ficha do gerador • **SV** Generatorplugg

Universal Active Cord

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Performance Characteristics

Active Cords are used for the transfer of electrical current to electrosurgical devices from a high-frequency (HF) electrosurgical generator. These devices connect Cook electrosurgical instruments to an electrosurgical generator. The maximum rated input voltage for this device is 5 kVp-p (2.5 kVp).

Patient Population

The device is indicated for adult use only.

Intended User

This device should only be used by or under the supervision of physicians who have been trained on its use.

Contact with Body Tissue

This device is not patient contacting. This device is not intended to come in contact with patient body tissue.

Operating Principle

Active cords are reusable devices that conduct HF electrical current from the electrosurgical generator to the electrosurgical accessory. The basic mechanism of action is conduction of current through the metal wire which is insulated from the external environment. The connectors on either end of the active cord allow proper attachment to the electrosurgical generator and to the electrosurgical instrument while maintaining electrical contact and ensuring conduction of the electrical current.

INTENDED USE

This device is used to connect Cook monopolar electrosurgical accessories to compatible electrosurgical generators.

INDICATIONS FOR USE

The active cords are used with the Cook electrosurgical accessories for the completion of endoscopic procedures such as polypectomy, sphincterotomy, and papillotomy.

CLINICAL BENEFITS

Delivers monopolar current from an electrosurgical unit to electrosurgical accessories for the completion of endoscopic procedures such as polypectomy, sphincterotomy, and papillotomy.

CONTRAINDICATIONS

None known.

WARNINGS

- Do not use this device for any purpose other than stated intended use.
- This device is not intended to be used within a sterile field.
- Any electrosurgical accessory constitutes a potential electrical hazard to patient and operator. Possible adverse effects include: burns • cardiac arrhythmia • fulguration • nerve and/or muscle stimulation.
- The electrosurgical effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the device.
- Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.
- Do not use this device with an output higher than the rated voltage 5 kVp-p (2.5 kVp). This could cause thermal injury to the patient, operator, or assistant, and could also damage the endoscope, device, and/or active cord.
- Do not use this device with any generator setting which may output a voltage exceeding the maximum input voltage for this device: 5 kVp-p (2.5 kVp). This could cause thermal injury to the patient, operator, or assistant, and could also damage the endoscope, device, and/or active cord.
- If cutting is not observed at normal operating settings, this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.
- The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.
- Skin-to-skin contact should be avoided (for example between the patient's arms and body) by using dry gauze.
- Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.
- Avoid contact between this device and the patient's body, or any other electrodes.
- This device must only be used in conjunction with a type BF or CF generator (generators that provide protection

from electric shock in compliance with IEC 60601-1). Use of a generator other than type BF or CF could cause electric shock to the patient, operator, or assistant.

- For recommendations on current settings, refer to the monopolar generator's recommended power settings.
- This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere, or in the presence of explosive gases. Doing so may result in explosion.
- In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy, or electrosurgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using this device on these patients.
- For patients with electrically conductive implants, monopolar diathermy, or electrosurgical cautery can result in localized thermal injury at the implant site.
- Additional warnings may be listed in procedural steps below.

PRECAUTIONS

Switch the electrosurgical unit to the "off" position when it is not in use.

Do not wrap or coil the cord tightly, as this may damage or cause the wire to break.

Do not wrap or coil cord around metal or other wiring/cords, as this may cause electromagnetic or other interference.

Do not remove the active cord from the device or electrosurgical generator by pulling on the wire, as this may damage or cause the wire to break.

POTENTIAL COMPLICATIONS

Provided by electrosurgical accessory manufacturer.

TECHNICAL DESCRIPTION

See fig. 1 for ACU-1 Active Cord and fig. 2 for ACU-1-VL Active Cord.

Refer to the package label for device dimensions.

Maximum rated input voltage for this device is 5 kVp-p (2.5 kVp).

This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6°F).

WARNING

No modification of this equipment is allowed.

Device Compatibility

This device is recommended for use with Cook monopolar electrosurgical accessories.

This device must be used with an electrosurgical accessory with a 3mm diameter connector pin.

Electrosurgical generator - type BF or CF (generators that provide protection from electric shock in compliance with IEC60601-1).

HOW SUPPLIED

This device is supplied non-sterile and is reusable.

STORAGE

Store in a dry location. Recommended storage temperature 22°C (71.6°F).

Device may be exposed to temperatures up to 50°C (122°F) during transportation.

INSPECTION OF DEVICE

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

DEVICE PREPARATION

PRECAUTION

Before using this device, inspect device for damage such as, damaged wire insulation, loose connections, cracked plug, metal corrosion, etc. If damaged in any way, do not use. Dispose of the device per institutional guidelines for biohazardous medical waste.

WARNINGS

For recommendations on current settings, refer to the monopolar generator's recommended power settings.

Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

1. With the electrosurgical unit off, prepare the equipment.

INSTRUCTIONS FOR USE

PRECAUTION

Before using this device, ensure the device has been cleaned and dried per the reprocessing instructions.

WARNINGS

Ensure the electrosurgical unit is in the "off" position before connecting the active cord to the electrosurgical accessory or the generator. Failure to do so could result in injury to the patient or user.

The active cords must be used with an electrosurgical accessory with a 3mm diameter connector pin. Active Cord fittings should fit snugly into both the device handle and the electrosurgical unit. If either plug/fitting does not fit snugly, do not use. Dispose of the device per institutional guidelines for biohazardous medical waste.

1. Connect active cord to the Cook monopolar electrosurgical device by placing the accessory plug/ fitting onto the device pin until fully seated. The accessory plug/fitting should fit snugly into the device handle.
2. Connect the active cord to the monopolar receptacle of the electrosurgical unit. The active cord should fit snugly into the electrosurgical unit.

- Follow the recommendations provided by the electrosurgical unit manufacturer and the electrosurgical accessory manufacturer to perform the procedure. The maximum rated input voltage for this device is 5 kVp-p (2.5 kVp).
- Upon completion of the procedure, turn the electrosurgical unit off.
- Disconnect the active cord from the electrosurgical unit and accessory by gripping and pulling the plug/fitting. Never disconnect any plug/fitting by pulling the active cord wire, as this may damage or cause the wire to break.
- Clean the active cord as soon as possible after each use following the reprocessing instructions.

REPROCESSING

WARNINGS

This device is not intended to be sterilized.

This device is not intended to be cleaned or disinfected by a machine.

This device is not intended to be submerged in water.

Frequent reprocessing can damage this product. Observe system preparation precautions for checking damage to the device. If the device is damaged in any way, do not continue cleaning or use. Dispose of the device per institutional guidelines for biohazardous medical waste.

Equipment/Supplies

- Disposable lint-free cloth
- Soft-bristled brush
- Cold tap water
- pH neutral detergent (e.g., Steris Prolystica etc.)
- Clean, filtered compressed air

Cleaning

- Prepare a cleaning solution using a pH neutral detergent according to the detergent manufacturer's instructions.
- Using a clean, lint-free cloth saturated in cold tap water, remove any visible liquids or contaminants from the device.
- Rinse the device of visible liquids or contaminants using cold running water for thirty (30) seconds.
- Using the prepared cleaning solution and a soft bristle brush, brush all surfaces of the device for one (1) minute ensuring particular attention is paid to small crevices in and around plugs/fitings.
- Rinse the device under cold running water for thirty (30) seconds.
- Thoroughly dry the device using a soft, lint-free cloth. Clean, filtered compressed air may be used to dry small crevices in and around plugs/fitings.
- Inspect the device for damaged wire insulation, loose connections, cracked plug, metal corrosion, etc. If device is damaged in any way do not continue use. Dispose of the device per institutional guidelines for biohazardous medical waste.

- Loosely coil the active cord and place in the original package or similar container.

DISPOSAL OF DEVICE

After the procedure, inspect the device for damaged wire insulation, loose connections, cracked plug, metal corrosion, etc. If device is damaged in any way do not continue use. Dispose of device in accordance with institutional guidelines.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, the competent authority, and/or regulatory authority of the country where the device was used.

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© COOK 2025
Wilson-Cook Medical, Inc.
Cook Medical Europe LTD
William A. Cook Australia Pty. Ltd.