**EN 4** **Esophageal • Pyloric • Colonic Balloon Dilators**  
Instructions for Use

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

1a.

**Eclipse TTC Wire Guided Balloon Dilator****Quantum TTC® Colonic Balloon Dilator****Quantum TTC® Esophageal Balloon Dilator**

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

**Please read all instructions before using this device.**

**DEVICE DESCRIPTION**

Balloon length, balloon diameter, and recommended maximum inflation pressure are provided in **Table 1**, and on the product label and labeling on the device.

**Table 1**

Product	Balloon Length	Recommended Maximum Inflation Pressure	Balloon Diameter
ECL-6X5.5	5.5 cm	120 psi / 8 atm	18 Fr / 6 mm
QDC-6X5.5			
QD-6X8	8 cm		
ECL-8X5.5	5.5 cm	120 psi / 8 atm	24 Fr / 8 mm
QDC-8X5.5			
QD-8X8	8 cm		
ECL-10X5.5	5.5 cm	90 psi / 6 atm	30 Fr / 10 mm
QDC-10X5.5			
QD-10X8	8 cm		
ECL-12X5.5	5.5 cm	90 psi / 6 atm	36 Fr / 12 mm
QDC-12X5.5			
QD-12X8	8 cm		
ECL-14X5.5	5.5 cm	60 psi / 4 atm	42 Fr / 14 mm
QDC-14X5.5			
QD-14X8	8 cm		
ECL-16X5.5	5.5 cm	60 psi / 4 atm	48 Fr / 16 mm
QDC-16X5.5			
QD-16X8	8 cm		
ECL-18X5.5	5.5 cm	50 psi / 3.5 atm	54 Fr / 18 mm
QDC-18X5.5			
QD-18X8	8 cm		
ECL-20X5.5	5.5 cm	50 psi / 3.5 atm	60 Fr / 20 mm
QDC-20X5.5			
QD-20X8	8 cm	90 psi / 6 atm	

**Patient Population**

The device is indicated for adult use only.

**Intended Users**

Use of this device is restricted to a trained healthcare professional.

**INTENDED USE**

This device is used to dilate strictures of the gastrointestinal tract, including strictures of the esophagus, pylorus, duodenum and colon.

**INDICATIONS FOR USE**

These devices are accessory endoscopic devices used to treat adult patients with esophageal strictures, gastric strictures, and colonic strictures of the gastrointestinal tract. **This does not include balloons used for dilation of the biliary tree or in the treatment of achalasia.**

**CLINICAL BENEFITS**

This device relieves dysphagia and lumen narrowing/obstruction in the gastrointestinal (GI) tract.

**CONTRAINDICATIONS**

Those specific to the primary endoscopic procedure to be performed in gaining access to the desired dilation site.

Those specific to dilation, including: uncooperative patient • asymptomatic rings, webs, or strictures • inability to advance the balloon dilator through the strictured area • coagulopathy • known or suspected perforation • severe inflammation or scarring near the dilation site.

**WARNINGS**

This device is designed for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Recommended 100% balloon inflation pressure can be found on the labeling of the balloon dilator.

Do not advance the balloon dilator if resistance is encountered. Assess the cause of the resistance to determine if dilation should be reattempted.

The Eclipse TTC Wire Guided Balloon Dilator (ECL) contains nickel which may cause an allergic reaction and other potential systemic complications related to nickel allergy.

**PRECAUTIONS**

- Do not use this device for any purpose other than stated intended use.
- Refer to the package label for the minimum channel size required for this device.
- Do not pre-inflate the balloon. Prior to insertion of the

balloon dilator, negative pressure is mandatory to maintain balloon deflation.

- The endoscope should remain as straight as possible when advancing or withdrawing the balloon dilator.
- The entire balloon should be extended beyond the tip of endoscope, and be completely visualized and positioned, before inflation.
- This device is used in conjunction with an inflation device or manometer and fluid-filled syringe. Do not use air or gaseous substances to inflate the balloon, as this will result in reduced balloon effectiveness.
- During withdrawal of the balloon dilator from the endoscope, negative pressure is mandatory to maintain balloon deflation.
- Apply a water-soluble lubricant to the balloon to allow easier passage through the accessory channel.

## POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy, including: perforation • hemorrhage • aspiration • fever • infection • allergic reaction to medication • hypotension • respiratory depression or arrest • cardiac arrhythmia or arrest.

## HOW SUPPLIED

This device is supplied non-sterile and may not be sterilized by any method. Store in a dark, dry, cool place. Avoid extended exposure to light and extremes of temperature and humidity.

## SYSTEM PREPARATION

1. Prepare the inflation device per the manufacturer's instructions.
2. Attach the inflation device connecting tube to the sideport of the y-connecting cap of the balloon dilator for the Eclipse (*See fig. 1*) or to the hub of the balloon dilator for the Quantum. (*See fig. 1a*) Secure by rotating the fittings. **Note:** If use of a stopcock is desired, attach the appropriate fittings to the balloon dilator and the inflation device connecting tube. To secure the attachment, rotate the fittings. The arm of the stopcock indicates which port is closed.
3. Aspirate all residual air from the balloon.
4. Create and maintain negative pressure on the balloon dilator using the inflation device. The balloon dilator is now ready for placement through the accessory channel of endoscope.

## INSTRUCTIONS FOR USE

### Balloon Introduction

1. Remove the protective sheath from the balloon.
2. Apply a water-soluble lubricant to the balloon to allow easier passage through the accessory channel. **Note:** Maintain negative pressure on the balloon dilator during introduction.
3. Introduce the deflated balloon into the accessory channel

and advance in short increments until the dilator is completely visualized endoscopically.

4. Using fluoroscopic or endoscopic monitoring, continue advancement until the desired position is achieved.
5. Before inflation, ensure the balloon fully exits the endoscope and is in the desired position within the stricture. **Note (if applicable):** To use with a wire guide, remove the stylet and flush the catheter lumen. Insert the wire guide.
6. Refer to the inflation device Instructions for Use for detailed instructions on using the inflation device to inflate and deflate balloon dilators.
7. Proceed with dilation. Refer to 100% inflation pressure found on the balloon dilator labeling.

### Balloon Withdrawal

8. Once the balloon is fully deflated, create and maintain negative pressure on the balloon dilator using the inflation device, prior to balloon withdrawal.
- Caution:** A compromised balloon may prohibit removal from the endoscope accessory channel. Removal of the endoscope along with the compromised balloon may be required.
9. To remove the deflated balloon from the accessory channel, straighten the endoscope tip and withdraw the balloon using a continuous twisting motion. **Caution:** Balloon must be completely deflated with all fluid removed before withdrawal. **Note (if applicable):** Previously positioned wire guide may be left in position to facilitate introduction of other wire-guided devices.

## EQUIPMENT REQUIRED

- Cook Inflation Device (sold separately)
- Wire guide (not required and sold separately) - wire guide compatibility is provided on the product label where applicable.
- Endoscope with compatible accessory channel - minimum accessory channel is provided on the product label.
- Water-soluble lubricant

## DISPOSAL OF DEVICE

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

## PATIENT COUNSELING

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

## SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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© COOK 2021
**William A. Cook Australia Pty. Ltd.**
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**