COOK®

MEDICAL

Patient Care Manual

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This manual provides information on the use and care of the Passport low profile gastrostomy device. These are only manufacturer guidelines and are not intended to replace physician's instructions.

Patient Information

As the caregiver, ask the patient's physician or health care provider to help complete the following section in order to note all of the physician's instructions for the patient. The information will be helpful when you begin feedings and will also give any substitute caregiver the necessary information to give the feedings correctly and safely.

Your Patient's Specific Needs

Name _____ has been fitted with a Cook Passport.

Manufacturer's Lot# _____

French Size _____ cm Length _____

Date Passport was inserted _____

Formula

Formula Name _____

Concentration _____

Ready-To-Use Formula

Total amount to be fed per day : _____ formula and _____ water.

Powdered Formula

Total amount to be fed per day : _____ formula.

Mixed _____ powder to _____ water.

Total additional water to be fed per day : _____ water.

Checking Residual (optional per physician's instructions)

- Check residual before feeding or every _____ hour(s). Call the patient's physician or health care provider if amount is greater than _____ cc.

Type of Feeding

Bolus

Feeding schedule: Every _____ hour(s).

For each feeding: _____ cc of formula.

RATE: Deliver over _____ minutes or _____ cc/minute.

Give _____ cc of water after each feeding.

Continuous

RATE: Feed _____ cc of formula/hour for _____ hour(s) or

continuously. (Follow manufacturer's instructions for setting pump rate.)

Give _____ cc of water every _____ hour(s).

Decompression (if prescribed)

Bolus

Decompress for _____ minutes before/after feeding.

Continuous

Decompress for _____ minutes before/after feeding every _____ hour(s).

Special Instructions

For Help with Problems, Call:

Physician _____ Phone # _____

Health Care

Provider _____ Phone # _____

Basics of Tube Feeding

The Passport low profile gastrostomy device is used to replace a gastrostomy tube through an existing stoma tract. This device provides continued access to the stomach for the delivery of nutrition, medication, gastric drainage or decompression.

The physician has prescribed the amount and type of formula to be fed each day, the amount of water to be given with the feedings, and the particular feeding method. The formula contains the proper mix of nutrients (proteins, carbohydrates, fats, vitamins, and minerals) that the body needs every day to maintain itself.

The patient's physician has prescribed when the feedings should be given and the method, either bolus or continuous feeding method. With the bolus method, the patient is fed a large amount quickly, several times a day. This is accomplished using a bulb or plunger syringe. The continuous feeding method provides the same amount of formula but at a slower, constant flow rate over a period of 8-24 hours. Often this method is accompanied with the use of a pump.

To prevent vomiting, the patient should be in an erect or semi-erect position during the feeding and for 30-60 minutes after. Patients in bed should be propped up to at least a 30° angle. You can use pillows to prop up older patients. Patients able to move around can sit up in a chair during feeding or, if the physician allows, stand and walk around during feeding.

Feeding the Patient

Preparation for Feeding

Thoroughly wash hands with soap and water before preparing formula or having any contact with the Passport.

The Passport should be checked before each feeding to be sure that it is in place and not clogged by following the steps below.

1. Observe the Passport exit site for signs of infection. Signs include drainage, unusual tenderness, swelling, or redness. Contact your physician if any of these symptoms are observed.

Note: The physician may or may not recommend step #2, described in the next 2 paragraphs below, prior to feeding.

2. Draw an appropriate amount of air into a 60 cc syringe then attach the tip of the syringe to the soft part of the Direct Kinect Adapter. Place the hard plastic tip of the Direct Kinect adapter into the circular valve until it is firmly engaged. Place a stethoscope on the left side of the abdomen just above the waist and depress the plunger. Listen for a bubbling or gushing sound with the stethoscope. **If you cannot hear the sound, do not attempt the feeding. Contact the patient's physician or healthcare provider and report the problem.** After you have heard the bubbling or gushing sound, gently draw back on the syringe plunger to withdraw the stomach contents or residual.

Note: You may not get any residual if the stomach is empty.

If the amount of the residual is less than the amount prescribed by the patient's physician in the "Patient Information" section, inject it back into the tube, draw up an appropriate amount of water into the syringe and rinse the Passport by injecting the water into it. If the residual is more than the physician prescribes, inject the contents into the tube and wait 30-60 minutes. Check the residual again. If the residual is still high, do not attempt the feeding and report the problem to the patient's physician or healthcare provider. Be sure to reinject the residual.

Flush the Passport with 30-50 cc of water every 8 hours and before and after each feeding unless otherwise directed by the patient's physician.

You may now begin giving the formula according to the method (bolus or continuous) the physician prescribed and noted in the "Patient Information" section.

Bolus Method:

1. Fill a 60 cc irrigation syringe as prescribed and attach to the soft part of the Direct Kinect Adapter.
2. Place the hard plastic tip of the Direct Kinect Adapter into the circular valve of the Passport until it is firmly engaged.
3. Inject the formula until prescribed amount of formula has been delivered to the patient.
4. Upon completion of injecting the prescribed amount of formula, flush the tube with 30-50 cc of water to prevent clogging.
5. Rinse the feeding syringe well with warm water and allow it to air dry. Store the feeding syringe until the next feeding.

Continuous Method:

Note: Steps may vary according to feeding pump manufacturer's guidelines for use.

1. Follow the physician's instructions for type and amount of formula. Insert the tip of the pump set directly into the Passport circular valve.

Note: The Passport circular valve fits most pump set tips. If your pump set tip does not firmly engage the Passport, use the right angle universal tube supplied with the Passport set. Clamp the right angle universal tube and place the "L" shaped tip into the circular valve until it engages. Place the pump set tip into the soft "Bolus" end of the right angle universal tube, then release the clamp.

2. Check that the tubing is connected to the pump properly. Open the clamp on the pump tubing.
3. Follow the manufacturer's instructions to set the pump at the rate prescribed by the patient's physician noted in the "Patient Information" section.
4. **Caution:** Do not allow the entire formula to empty before adding water. Air in the tube may enter the stomach and cause the patient discomfort.
5. When the feeding container is almost empty, clamp the pump tubing closed. Place the prescribed amount of water as noted in the "Patient Information" section in the feeding container.

6. Open the clamp on the pump tubing and re-check the pump rate, setting it to the prescribed flow rate.
7. When feeding is complete, clamp the pump tubing closed, turn the pump off, then disconnect the pump tubing from the Passport.

Decompression (if prescribed)

The patient's physician may or may not want you to decompress (let the air out of the stomach) before or after feedings. To decompress the stomach, place the "L" shaped or right angle tip of the Right Angle Universal Tube into the Passport circular valve until it is firmly in place. Allow the tube to drain into a collecting cup or drainage bag for the amount of time prescribed by the physician noted in the "Patient Information" section.

Giving Medications (if prescribed)

Medications may be given through the Passport if the physician allows. Ask your physician to order prescriptions in liquid form if possible.

1. If the medication is in tablet form, dissolve the prescribed amount in 30-50 cc of warm water. Or, pour the prescribed amount of liquid medication into the medication cup.
2. Draw the medication into a syringe.
3. Use the Bolus method as outlined under "Feeding the Patient" to administer medications.
4. Inject the medication into the Passport via the Direct Kinect Adapter.
5. Flush the tube with 50 cc of water to prevent clogging unless otherwise noted by your physician and then replace the feeding adapter cap.
6. The patient should remain upright or at a 30° angle for 30 minutes after giving the medication.

General Care of the Patient and PEG

Cleaning the stoma

The stoma should be cleaned daily with mild soap and warm water and allowed to air dry. Use a cotton swab or small piece of gauze to clean in a circular motion from the stoma site outward. Rotate and clean under the bolster. Check for signs of drainage, unusual tenderness, swelling, or redness and contact your physician or healthcare provider if any of these symptoms are observed.

Note: Dressings or bandages are not recommended. It is best to keep the area exposed to air.

Caring for the Mouth

The mouth should be cared for even if the patient is not eating. Brush the teeth, gums, and tongue every day. Patients may use a mild mouthwash if the physician allows. Lips should be moistened with water or a moisturizing lotion. Report any persistent bleeding of the gums to the patient's physician or healthcare provider.

Upset Stomach

If the patient has a nauseous feeling, do not force the feeding. Wait about an hour or so before continuing to feed or until the nauseous feeling subsides. Ask the patient's physician or healthcare provider if the formula may be diluted with water or if you may give smaller, more frequent feedings. Stop feeding if vomiting occurs. After waiting about an hour or so, you may try to feed again. If the vomiting occurs again, stop feeding completely. If the upset stomach does not stop enough to resume regular feedings in 24 hours, call the patient's physician or healthcare provider.

Aspiration

Aspiration means inhaling formula into the patient's lungs. Burping up small amounts of formula or vomiting can cause aspiration. If the patient begins choking or having difficulty breathing, immediately stop the feeding. The patient should be placed on his/her side or propped up while you, the care giver, calls the patient's physician or healthcare provider.

Glossary of Terms

PEG – procedure for placing a feeding tube into the stomach through the abdominal wall

Percutaneous (per-q-TAY-nee-us) – through the skin

Endoscopic (en-doh-SKOP-ik) – method for using a lighted instrument to assist in tube placement

Gastrostomy (ga-stras-te-me) – surgical opening into the stomach

Passport – low profile gastrostomy device used to replace a gastrostomy tube through an existing stoma tract

Stoma – opening in the abdomen (the feeding tube goes through the stoma to reach the stomach)

Nutrition – the proteins, carbohydrates, fats, vitamins, minerals, and water the body needs every day to maintain itself

Flush – wash out the feeding tube

Residual – stomach contents

Stethoscope – an instrument used for listening to sounds in the body

Decompression – let the air out of the stomach

Aspiration – inhaling formula into the lungs

Tube Migration – movement of the feeding tube

Problem Solving

Problem	Possible Causes	Actions
Tube clogged	<ul style="list-style-type: none"> Inadequate flushing after each use. 	<p>Prevention: Flush the tube with warm water before and after each feeding.</p> <p>To Unclog: Flush the tube using a 60 cc syringe filled with warm water. The tube may also be massaged to aid in unclogging. Care must be taken to avoid excessive pressure on the PEG tube. If both methods fail in unclogging the tube, call the patient's physician or healthcare provider.</p>
Accidental removal of the tube	<ul style="list-style-type: none"> Confused patients may try to remove the tube. 	<p>Remain calm and contact the patient's physician or healthcare provider immediately and save the Passport.</p>
Leakage around the base of the Passport	<ul style="list-style-type: none"> Passport may be too small 	<p>Contact your physician to determine if the Passport should be changed to a larger size.</p>
Skin around the stoma site seems to be infected (redness, soreness, pain, swelling)	<ul style="list-style-type: none"> Placement of the Passport Gastric leakage 	<p>Clean the stoma several times a day following the steps listed in the "Cleaning the Stoma" section. If skin irritation (redness, soreness, pain, swelling, unusual drainage) continues, contact the patient's physician or healthcare provider.</p>

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Passport is a registered trademark of Cook Endoscopy.
Cook is a registered trademark of Cook Incorporated.
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.