EN
2

Pancreatic Dilation Catheters

CS
2

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DA
3

Katetre til pancreasdilatation

NL
4

Pancreasdilatatiekatheters

FR
5

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INTENDED USE

This device is used to dilate stenotic minor papilla.

This device is indicated for adult use only.

NOTES

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to ERCP.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Do not use this device for any purpose other than stated intended use.

Use of this device restricted to a trained healthcare professional.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

If unexpected resistance is met during the procedure, do not continue advancement of catheter without first determining the cause of the unexpected resistance and taking remedial action.

If bleeding occurs during the procedure, stop the procedure immediately and take remedial action.

INSTRUCTIONS FOR USE

1. Introduce device into accessory channel of endoscope over pre-positioned wire guide. Advance in small increments until device is endoscopically visualized exiting scope.
2. Proceed with cannulation.
3. Under fluoroscopic monitoring, advance catheter to appropriate position to dilate stricture.
4. Proceed with dilation.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.