EN Zilver® Biliary Stent

CS Biliární stent Zilver®

DA Zilver® galdestent

NL Zilver® biliaire stent

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DE Zilver® Gallenstent

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ENGLISH

INTENDED USE

This device is used in palliation of malignant neoplasms in the biliary tree.

STENT DESCRIPTION

This stent is an open-ended cylinder of nitinol. The design of the stent resists radial compression and easily conforms to bends in the bile duct anatomy without collapsing. There are radiopaque markers on each end of the stent to assist in fluoroscopic visualization of the stent position.

DELIVERY SYSTEM DESCRIPTION

The stent is compressed onto a flexible delivery system. The delivery system accepts a .035 inch wire guide.

NOTES

Do not use this device for any other purpose than the stated intended use.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedure to be performed in conjunction with stent placement.

Additional contraindications include, but are not limited to: inability to pass the wire guide or stent through the obstructed area, biliary duct strictures of benign etiology, biliary obstruction preventing either endoscopic or percutaneous cholangiography, concurrent perforated bile duct, those patients for whom endoscopic or percutaneous procedures are contraindicated, patients with coagulopathy, concurrent bile duct stricture, strictures greater than 8 cm in length, and any use other than those specifically outlined under Intended Use.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Additional complications that can occur in conjunction with biliary stent placement include, but are not limited to: trauma to the biliary tract or duodenum; obstruction of the pancreatic duct; stent migration; tumor ingrowth, tumor overgrowth of stent ends, or excessive hyperplastic tissue ingrowth.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

This stent must be placed under fluoroscopic monitoring.

This stent must only be placed using the delivery system provided with the set.

This stent is intended for palliative treatment only. Alternate methods of therapy should be investigated prior to placement.

After stent placement, alternate methods of treatment such as chemotherapy and irradiation **should not be** administered as this may increase the risk of stent migration due to tumor shrinkage, stent erosion, and/or mucosal bleeding.

Long-term patency with this stent has not been established. Periodic evaluation of the stent is advised.

If the wire guide or stent cannot advance through the obstructed area, do not attempt to place the stent.

Assessment must be made to determine the necessity of sphincterotomy or balloon dilation prior to stent placement. In the event sphincterotomy or balloon dilation is required, all appropriate cautions, warnings and contraindications must be observed.

MR CONDITIONAL

Non-clinical testing has demonstrated that this stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient field of 720 Gauss/cm or less
- Whole body averaged specific absorption rate (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla), respectively.

In non-clinical testing, this stent produced maximum temperature rises of 0.1, 3.8, 0.8, and 0.1 degrees C (for a single stent at 1.5 Tesla, a pair of overlapping stents at 1.5 Tesla, a single stent at 3 Tesla, and a pair of overlapping stents at 3 Tesla, respectively) at whole body averaged specific absorption rates (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) of MR scanning in a 1.5 Tesla/64 MHz General Electric MR scanner, a 1.5 Tesla Magnetom Siemens Medical Solutions MR Scanner (to evaluate a pair of overlapping stents), and a 3 Tesla Excite General Electric MR scanner.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of this stent. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

Heating in the MRI environment for stents with fractured struts is not known.

WARNINGS

This stent is **not intended to be removed**. Attempts to remove the stent after placement may cause damage to the surrounding mucosa. This stent is considered to be a permanent implant.

This device is **not** intended to be deployed through the wall of a previously placed or existing metal stent. Doing so could result in difficulty or inability to remove introducer.

A complete diagnostic evaluation should be performed prior to placement to measure the stricture length and determine the proper stent length.

The safety and effectiveness of this device for use in the vascular system has not been established.

The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

EQUIPMENT REQUIRED

- Zilver Biliary Stent Delivery System with appropriate length stent
- .035 inch wire guide (recommended size)
- ERCP cannula
- Sphincterotome (if necessary)
- Biliary dilation balloon with inflation device (if necessary)

IMPORTANT: This stent system is intended for single use only. Do **not** attempt to reload a deployed or partially deployed stent.

PROCEDURE PREPARATION

1. With endoscope in place, introduce a wire guide, floppy tip first, and advance until it is fluoroscopically visualized in position through stricture.
2. Prepare previously selected stent delivery system by irrigating inner lumen and stent with sterile water.

INSTRUCTIONS FOR USE

SHORT WIRE EXCHANGE

I. If using Intra Ductal Exchange (IDE) port and a pre-positioned Ultra Short Wire (USW).

1. Unlock USW from Wire Guide Locking Device and advance stent delivery system over pre-positioned wire guide ensuring that wire guide exits guiding catheter at 3 mm IDE port.
2. Introduce stent delivery system into endoscope accessory channel and relock wire guide.
3. With elevator open, advance device in short increments until it is endoscopically visualized exiting scope.
4. Under fluoroscopic guidance, pass stent delivery system through papilla and into common bile duct. Advance until stent is fluoroscopically visualized through stricture. Fluoroscopically visualize radiopaque markers on either end of stent and position stent so that it bridges stricture completely. **Note:** Stents bridging papilla should extend beyond papilla and into duodenum approximately 0.5 cm after deployment.
5. Once stent delivery system is in correct position for deployment, loosen cap on proximal end of Y-body. (See Fig. 1) Stent is now ready for deployment.
6. Unlock wire guide from Wire Guide Locking Device **and completely remove wire guide from device**.
7. Begin deployment of stent by holding wire guide port hub stationary and **slowly** pulling back on Y-body, while monitoring position of stent fluoroscopically. **Warning:** Do **not** push delivery system up into bile duct after stent deployment has been initiated.
8. Continue pulling back on Y-body until it is fluoroscopically confirmed that stent is completely deployed.
9. Carefully remove delivery system from within expanded stent, and from endoscope. Dispose of delivery system per institutional guidelines for biohazardous medical waste. If necessary, complete expansion of stent using a biliary dilation balloon with a diameter less than stent nominal diameter. **Warning:** Never use a rigid dilator, as axial force may dislodge stent. **Warning:** These metal biliary stents are not intended to be repositioned or removed after deployment in bile duct. In case of accidental deployment or improper placement (immediately following deployment), stent should be left in place and placement of a second stent should be attempted to achieve desired result.

LONG WIRE EXCHANGE

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide.

Note: For best results, wire guide should be kept wet.

1. Remove stylet and advance introducer with preloaded stent over pre-positioned wire guide and into accessory channel of endoscope until proximal end of wire guide exits PWP. (See Fig. 1)

Note: Ensure that wire guide does not exit at 3 mm IDE port.

REFER TO STEPS 3-5 AND STEPS 7-9 IN “SECTION I.”

Note: It is recommended to re-advance the outer sheath to its pre-deployment position prior to removing the system.

Upon completion of procedure, remove wire guide and endoscope from patient. Dispose of wire guide per institutional guidelines for biohazardous medical waste.

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