COOK®

MEDICAL

CE
0123

FS-PRECU

EN
3

Fusion® Needle Knife

CS
6

Fusion® Jehlový papilotom

DA
9

Fusion® papillotomkanylekniv

NL
12

Fusion®-naaldmespapillotoom

FR
15

Papillotome-couteau-aiguille Fusion®

DE
18

Fusion® Nadelmesser

EL
22

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
ENGLISH

INTENDED USE

This device is used for accessing the common bile duct when standard methods of cannulation have been exhausted.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with papillotomy.

Contraindications to papillotomy include, but are not limited to: coagulopathy.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP.

Potential adverse events associated with papillotomy include, but are not limited to: pancreatitis, bleeding, perforation, and infection.

USER QUALIFICATIONS

This device should only be used by or under the supervision of physicians who have been trained in ERCP procedures and in precut papillotomy.

WARNINGS

Do not use this device for any purpose other than the stated intended use.

This device is supplied sterile and intended for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Any electrosurgical device constitutes a potential electrical hazard to the patient and the operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

The electrosurgical effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.

Before using this device, follow the recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through the proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

For patients with electrically conductive implants, monopolar diathermy or electrosurgical cautery can result in localized thermal injury at the implant site.

Do not use this device with an output higher than the rated voltage: 1.67 kVp-p (835 Vp). This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, device and/or active cord.

Do not use an active cord with a maximum rated voltage less than 1.67 kVp-p (835 Vp). This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, device and/or active cord.

If cutting is not observed at normal operating settings, this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.

The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.

Skin-to-skin contact should be avoided (for example between the patient's arms and body) by using dry gauze.

Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.

Avoid contact between the active cords and the patient's body, or any other electrodes.

This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electric shock in compliance with IEC 60601-1). Use of a generator other than type BF or CF could cause electric shock to the patient, operator, or assistant.

For recommendations on current settings for papillotomy, refer to the monopolar generator's recommended power settings for papillotomy.

This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere, or in the presence of explosive gases. Doing so may result in an explosion.

In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electrosurgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the device in these patients.

This device must be used with an endoscope that is intended for use with high frequency (HF) electrosurgical devices. If the device is used with an endoscope not compatible with high frequency (HF) electrosurgery, a burn to the patient or the operator may occur.

To avoid neuromuscular stimulation and minimize electric arc, activate the cutting wire only when it is in contact with the tissue to be cauterized.

Additional warnings may be listed in procedural steps below.

PRECAUTIONS

Switch the electrosurgical unit to the "off" position when it is not in use.

The elevator should remain open/down when advancing or retracting the device.

The wire guide diameter and the inner lumen of the wire-guided device must be compatible.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

Safe and effective electrosurgery is dependent not only on equipment design but, to a large extent, on factors under the control of the operator. Any electrosurgical accessory constitutes a potential electrical hazard to the patient and the operator.

In order to enhance safety and effectiveness, it is important that these cautions as well as the instructions supplied with the electrosurgical unit be read, understood and followed when using this device.

TECHNICAL DESCRIPTION

Figure 1 summarizes the device's characteristics. Prior to starting the procedure in which the device will be used, become familiarized with the distal mark location with respect to the cutting wire. Handle reference markings may be used for approximation of needle extension.

Refer to the package label for device dimensions, minimum endoscope compatibility, and wire guide compatibility.

This device is recommended for use with Cook electrosurgical compatible wire guides.

This device must only be used with an active cord compatible with a 3 mm diameter connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

The maximum rated input voltage for this device is 1.67 kVp-p (835 Vp).

This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6°F).

Store in a dry location. Recommended storage temperature 22°C (71.6°F). Device may be exposed to temperatures up to 50°C (122°F) during transportation.

Components highlighted with an asterisk (*) denote applied parts (the part of the device that comes in contact with the patient). Prior to starting the procedure in which the device will be used, become familiarized with the distal mark location with respect to the cutting wire.

WARNING

No modification of this equipment is allowed.

SYSTEM PREPARATION

WARNING

For recommendations on current settings for papillotomy, refer to the monopolar generator's recommended power settings for papillotomy.

Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return

electrode to the electrosurgical unit is maintained throughout the procedure.

Do not extend or retract the cutting wire while the device is coiled or the stylet wire is in place, as this may cause damage to the papillotome and render it inoperable. Extend and retract the cutting wire only when the papillotome is straight and the stylet wire is removed.

1. Upon removing the device from the package, uncoil the papillotome and carefully remove the precurved stylet wire from the distal tip of the device.
2. With the electrosurgical unit off, prepare the equipment.

INSTRUCTIONS FOR USE

1. With the needle knife fully retracted into the sheath and the finger ring fully retracted to the most proximal position and handle lock locked in place with the thumbscrew, advance the tip of the device into the endoscope accessory channel and continue to advance in short increments until the device is endoscopically visible.

Note: If using the Intra Ductal Exchange (IDE) port and a wire guide (see fig.1) Attach the Wire Guide Locking Device to the endoscope accessory port. Remove the Wire Stop handle from its retaining clip. Ensure that the Wire Stop is disengaged from the Proximal Wire Port (PWP) hub and ensure that the distal tip of Wire Stop is proximal to IDE port. Insert the distal tip of the wire guide into the IDE port and advance until flush with the distal tip of the papillotome. Advance the Wire Stop handle until it reaches the PWP hub and connect the Luer lock securely to the PWP hub. With the cutting wire fully retracted into the catheter, advance the tip of papillotome through the cap of the Wire Guide Locking Device and continue advancing until endoscopically visible. Disengage the Wire Stop by releasing the Wire Stop handle from the PWP hub and remove the Wire Stop from the catheter.

Note: If using the Proximal Wire Port (PWP) and a pre-positioned long wire guide (see fig.1). For the best results, the wire guide should be kept wet with sterile water or saline. Remove the Wire Stop. With the cutting wire fully retracted into the catheter, advance the wire guide through the PWP until it reaches the distal end of the catheter. Advance the tip of the papillotome through the accessory cap and continue advancing until it is endoscopically visible.

WARNING

The device should never be connected to the active cord prior to inserting the device through the endoscope to avoid injury to the patient and/or equipment resulting from improper electrical circuit grounding.

It is essential to tighten the thumbscrew after adjusting the cutting wire length and before proceeding with papillotomy.

Do not extend the needle knife length beyond what is minimally necessary to complete the procedure. The device is intended to be adjustable to accommodate a wide range of patient anatomy.

2. Once the endoscope position has been set and the device tip is visualized existing the endoscope, position the catheter and extend the needle knife to the desired length

(see Fig 1.) Needle knife extension may be adjusted by loosening the thumbscrew and advancing or retracting the adjustable ring on the handle. If the endoscope position is adjusted after extending the needle, re-examine the needle length and make adjustments as necessary. Once final adjustments are made, lock the thumbscrew immediately prior to cutting.

PAPILLOTOMY

WARNING

Using the endoscope monitor, verify that the distal tip of the catheter has exited the endoscope by visualizing the distal end of the catheter is past the elevator. Failure to confirm the distal tip of the catheter has fully exited the endoscope may result in contact between the cutting wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

Examine all components to be used during the procedure. Do not use a wire guide or device that has been cut, burned, or damaged. Do not use a device with a handle that has been cracked, broken or has a damaged metal active post. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

If a non-insulated wire guide is used in the papillotomy, it must be removed prior to applying electrosurgical current. The device should not be energized prior to performing papillotomy. Energizing the cutting wire prior to use will cause premature cutting wire fatigue and will compromise the cutting wire's integrity.

To avoid unintended disconnection during the procedure and damage to the device, the device should only be used with compatible active cords (Reference Technical Description for compatibility).

1. Attach the active cord to the connector on the device handle. The active cord fittings should fit snugly in both the device handle and the electrosurgical unit.
2. Following the electrosurgical unit manufacturer's instructions, verify the desired settings and proceed with the papillotomy.

WARNING

It is essential to move the cutting wire while applying current. Maintaining the cutting wire in one position may cause excessive focal coagulation, charring of tissue and/or damage to the cutting wire.

The user must maintain contact with tissue when applying electrocautery current. Failure to do so may result in increased current resulting in a broken cutting wire, damage to the endoscope, and/or patient injury.

When applying current, ensure the cutting wire is completely out of the endoscope. Contact of the cutting wire with the endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

3. When performing papillotomy, use a continuous motion while activating the device.

- Upon the completion of the papillotomy, turn the electrosurgical unit off and retract the cutting wire into the sheath.
- Disconnect the active cord from the device handle and from the electrosurgical unit. Clean and store the active cord per the manufacturer's instructions..

CANNULATION

- After papillotomy, the device may be used to cannulate the ductal system.
- Following cannulation, contrast may be injected through the injection port to fluoroscopically confirm the position of the device.
- Once access to the ductal system is obtained and contrast is injected (if applicable), a wire guide may be advanced to facilitate the introduction of compatible devices (See fig. 1).

- A previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices. If a wire guide is to remain in place while the device is withdrawn, utilize the following steps:

Note: If using the Intra Ductal Exchange (IDE) port and a wire guide prior to removing the device, utilize the reference marks on the catheter to ensure that the IDE port is within the ductal system. Fluoroscopically visualize the radiopaque band at the IDE port. Retract the wire guide until the radiopaque distal tip of wire guide passes the band; a disengagement from the wire guide lumen will occur. Advance the disengaged wire guide to maintain ductal access.

DEVICE REMOVAL FROM ENDOSCOPE

- Remove the needle knife, leaving the wire guide in place (if applicable).

Note: If using the Intra Ductal Exchange (IDE) port and a wire guide lock the wire guide into the Wire Guide Locking Device and remove the papillotome from the endoscope accessory channel.

Note: If using the Proximal Wire Port (PWP) and a pre-positioned long wire guide remove the device using standard long wire exchange technique.

DEVICE DISPOSAL

- Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.