EN

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ENGLISH

INTENDED USE

This device is used to dilate strictures of the biliary tree. The Biliary Inflation Device is used for inflation, deflation and pressure monitoring of the Biliary Balloon Dilator.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Sphincteroplasty is associated with increased incidence of pancreatitis.

CONTRAINDICATIONS

Those specific to ERCP and any procedures to be performed in conjunction with balloon dilation, including, but not limited to: asymptomatic strictures; inability to advance the balloon dilator through the strictured area; known or suspected perforation; severe inflammation or scarring near the dilation site.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Do not inflate the balloon prior to introduction into the scope, as this may cause damage to the scope.

The duodenoscope must remain as straight as possible with the elevator in the open position when advancing or withdrawing the balloon dilator.

The entire balloon must be extended beyond the tip of the duodenoscope and be fluoroscopically visualized and positioned before inflation.

The balloon dilator is used in conjunction with the inflation device or a manometer and fluid-filled syringe. **DO NOT** use air or gaseous substances to inflate the balloon, as this will result in reduced balloon effectiveness.

During withdrawal of the balloon dilator from the duodenoscope, negative pressure is mandatory to maintain balloon deflation.

A compromised balloon may prohibit removal from the endoscope accessory channel. Removal of the endoscope along with the compromised balloon may be required.

WARNINGS

Do not exceed the recommended balloon inflation pressure as listed on the catheter tag of the balloon dilator.

Do not advance balloon dilator if resistance is encountered. Assess the cause of resistance to determine if dilation should be re-attempted.

SYSTEM PREPARATION

1. Ensure the plunger handle of the inflation device is rotated completely clockwise and is pushed all the way in. To move the plunger, squeeze and hold the plunger lever located to the side of the syringe cylinder. Release the lever to lock the plunger in position.
2. Rotate the stopcock and the connecting tube fittings until tight. The stopcock is provided to facilitate purging of air.
3. Verify the sideport of the stopcock is closed by turning the "off" arm of the stopcock to the sideport.
4. While holding the inflation device, submerge the connecting tube in a minimum of 30 ml/cc of sterile water, saline or dilute contrast. Fill the device with a minimum of 15 ml/cc by compressing the plunger lever and pulling the plunger handle. Release the plunger lever to lock the plunger in place. (See fig. 1)
5. To completely purge the inflation device of air, hold the inflation device upright and compress the plunger lever while simultaneously grasping and pushing the plunger handle until all air is expelled. (See fig. 2) While maintaining the position of the plunger handle, release the plunger lever to lock the plunger in place. Release the plunger handle. To add fluid to the syringe, repeat steps 4 & 5.
6. Attach the male Luer lock fitting of the stopcock to the hub of the balloon, and secure by rotating the fittings. (See fig. 3)
7. While compressing the plunger lever, pull the plunger handle until vacuum is indicated on pressure gauge. (See fig. 4) Maintain a vacuum with the handle, then release the plunger lever. The balloon dilator is now ready for placement through the accessory channel of the duodenoscope.

INSTRUCTIONS FOR USE

Introduction

1. Remove the protective sheath from the balloon. Apply a water soluble lubricant to the balloon to allow easier passage through the accessory channel.
2. Introduce the deflated balloon into the accessory channel of the duodenoscope, then advance in short increments until the dilator is completely visualized endoscopically. **Note:** The balloon may be introduced over a pre-positioned .035 inch or smaller wire guide. Prior to advancing the device over a pre-positioned wire guide, flush the wire guide lumen with sterile water or saline to facilitate advancement.
3. Following cannulation, contrast may be injected to fluoroscopically confirm the position of the device. If a wire guide is in place, remove the wire guide and inject contrast through the wire guide port.
4. Using fluoroscopic monitoring, continue advancement into the duct until the desired position is achieved. **Note:** The position may be verified using the metal bands located at either end of the balloon as fluoroscopic markers. If the wire guide was removed for contrast injection, flush the contrast from the catheter with sterile water or saline before reintroducing the wire guide, observing the necessary cautions.

Balloon Inflation

Note: The balloon dilator must be inflated with **LIQUID ONLY**.

5. For rapid inflation, compress the plunger lever, then push the plunger handle forward until significant resistance is met. (See fig. 5) Release the plunger lever while maintaining plunger handle position. (See fig. 6)
6. Rotate the plunger handle clockwise to increase to the specified balloon pressure for the balloon dilator in use. Recommended pressures are listed on the catheter tag of the balloon dilator. Pressure may be decreased by rotating the plunger handle counterclockwise.
7. Monitor the balloon fluoroscopically during the procedure.

Balloon Deflation

8. To decrease balloon pressure and deflate, rotate the plunger handle counterclockwise until the pressure on the inflation device decreases below 20 psi. Compress the plunger lever, then pull plunger handle until negative pressure, as seen on the pressure gauge, is created. (See fig. 7) Release the plunger lever to lock in negative pressure. Release the plunger handle. **Note:** The indicator on the pressure gauge should be in the vacuum portion of the gauge during deflation.
9. Maintain negative pressure and observe the balloon for complete deflation. Negative pressure may be increased by continued counterclockwise rotation of the plunger handle.
10. Observe balloon size. If the balloon is not completely deflated, open the stopcock by rotating the "off" arm to close the balloon catheter port. Purge the inflation device of air. To continue balloon deflation, turn the stopcock to close the open sideport. Create negative pressure again as indicated in steps 8 & 9. Repeat these steps as necessary until the balloon is completely deflated.

Balloon Withdrawal

Caution: The balloon must be thoroughly deflated, with all fluid removed, before withdrawal.

11. To remove the deflated balloon from the accessory channel, straighten the duodenoscope tip and open the elevator. Maintain negative pressure while withdrawing the balloon using a continuous twisting motion. The previously placed wire guide may be left in position, in order to facilitate introduction of other wire guided devices.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.