EN Soehendra® Lithotriptor

DA Soehendra® litotriptor

NL Soehendra® Lithotriptor

FR Lithotripteur Soehendra®

DE Soehendra® Lithotriptor

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ENGLISH

INTENDED USE

This device is used to mechanically crush stones in the biliary duct when other methods of endoscopic removal have failed.

NOTES

This device is supplied non-sterile and reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Handle and cable are sold separately.

CONTRAINDICATIONS

Those specific to ERCP and mechanical lithotripsy.

Those associated with biliary stone extraction include, but are not limited to: an ampullary opening inadequate to allow for the unimpeded passage of the stone and basket.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to contrast or medication, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with basket extraction include, but are not limited to: impaction of the object, localized inflammation, pressure necrosis.

PRECAUTIONS

Wire cutters must be available when using this device.

Assessment of stone size and ampullary orifice must be made to determine necessity of sphincterotomy prior to this procedure. The length of intramural segment must be evaluated before determining advisability of mechanical lithotripsy. In event sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

Reusability of device depends in large part on care of device by user. Factors involved in prolonging life of this device include, but are not limited to: thorough cleaning following instructions included in this booklet.

WARNINGS

Due to mechanical pressure generated with this device, basket fragmentation in common bile duct is a possibility that may require surgical intervention.

Only select Cook biliary Soft Wire Baskets are recommended for use with this device.

Do not remove sheath from basket wire. Doing so may result in basket fragmentation and consequently require surgical intervention.

Due to varying compositions of biliary stones, stone fracture may not be possible. If stone cannot be fractured, continued rotation of handle may cause basket wire to break, requiring surgical intervention.

INSTRUCTIONS FOR USE

1. Examine Luer lock fittings on cable and handle for a proper fit. Ratchet on handle must attach securely to allow rotation in one direction only.
2. With basket fully retracted into sheath, introduce a soft wire multifilament basket, minimum of 220 cm in length, into endoscope accessory channel. Advance device in short increments until basket sheath exits endoscope.
3. Fluoroscopically visualize extraction target, then advance sheath of basket beyond target.
4. After confirming desired position of basket sheath relative to stone, advance basket out of sheath by pushing forward on handle. Place basket around stone and slowly retract handle until entrapment is achieved. **Note:** It may be necessary to move basket back and forth along target in order to capture object.
5. Once stone is completely captured in basket, firmly tighten basket wires around stone.
6. Slightly kink basket wire and sheath below handle and reinforcement sheath. (See fig. 1)
7. Using wire cutters, cut handle off basket. Do not remove sheath from basket wire. (See fig. 2)
8. Fluoroscopically monitor basket and stone to ensure continued entrapment while slowly removing endoscope.
9. Gently feed lithotripsy cable over basket sheath and advance cable, under fluoroscopic monitoring, until it reaches basket in duct.
10. Insert basket wire and sheath through male Luer lock fitting on lithotripsy handle, then through hole in rotating rod of handle.
11. Wrap wire and sheath several times around rod to secure them to rod. (See fig. 3)
12. To fracture stone, apply pressure by **slowly** rotating handle. **Note:** Slow, constant pressure is often most effective method of fracturing a stone. Sudden release of pressure may indicate fracture of stone, fracture of basket, or loss of stone from basket. **Note:** Handle will only move in one direction unless ratchet is released.
13. After stone has been fractured, assess need to remove fractured stone from duct.
14. Release ratchet on handle of lithotriptor to remove basket. Examine basket to confirm that it is not broken or damaged in any way, then dispose of basket per institutional guidelines for biohazardous medical waste.
15. Upon completion of procedure, clean and prepare handle and cable for future use following cleaning instructions that follow. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, bends or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.

CLEANING INSTRUCTIONS

1. Immediately after use, submerge handle and cable in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
2. Wipe both handle and cable with a textured cloth soaked in enzymatic solution.
3. Gently brush handle and cable with a soft nylon bristle brush and enzymatic solution to remove all foreign matter.
4. Remove handle and cable from enzymatic solution, then rinse both thoroughly with clean, running water.
5. Ultrasonically clean handle and cable for a minimum of 15 minutes at 25 °C (77 °F) or per ultrasonic cleaner manufacturer's instructions.
6. Remove handle and cable from ultrasonic cleaner, then rinse thoroughly with clean, running water.
7. Gently wipe handle and cable dry with a lint-free cloth.
8. Package handle and cable for sterilization following AAMI recommended practices.
9. Sterilize handle and cable using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4 °C (130 °F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4 °C (130 °F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilization, it should be taken out of service immediately.

AUTOClave CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	132.2 °C (270 °F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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


"Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers", *AAMI Technical Information Report*, No. 12 – 1994.

"Sterilization, Part 1 – Good Hospital Practices", *AAMI Standards and Recommended Practices*, Volume 1.1 – 1995.

 **Wilson-Cook Medical, Inc.**
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

© 2011 Cook Medical

RxOnly

EC	REP
-----------	------------

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland