**EN**  
7

**Saeed® Multi-Band Ligator**

**CS**  
10

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**DA**  
13

**Saeed® ligator med flere bånd**

**NL**  
16

**Saeed® multiband ligator**

**FR**  
19

**Ligateur multiple de Saeed®**

**DE**  
22

**Saeed® Multiband-Ligator**

**EL**  
26

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















12.

**INTENDED USE**

This device is used to endoscopically ligate esophageal varices at or above the gastroesophageal junction or to ligate internal hemorrhoids.

**NOTES**

This device is supplied non-sterile.

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

**CONTRAINDICATIONS**

Those specific to primary endoscopic procedure to be performed in gaining access to desired banding site.

Those specific to esophageal banding include, but are not limited to: cricopharyngeal or esophageal narrowing or stricture, tortuous esophagus, diverticula, known or suspected esophageal perforation, asymptomatic rings or webs, coagulopathy.

Those specific to hemorrhoidal banding include, but are not limited to: severe inflammatory bowel disease, portal hypertension.

Use of ligation bands is contraindicated in patients with a known hypersensitivity to latex.

**WARNINGS**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

**POTENTIAL COMPLICATIONS**

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those which can occur with esophageal banding include, but are not limited to: retrosternal pain, nausea, laryngeal laceration, esophageal perforation, stricture formation, obstruction.

Those which can occur with hemorrhoidal banding include, but are not limited to: severe pain, bleeding, urinary symptoms, swelling and edema, tissue ulceration, bands slipping off.

Hemorrhoidal banding may result in severe pain if the procedure is performed below the dentate line.

**PRECAUTIONS**

Refer to package label for minimum channel size required for this device.

Do not use with an overtube. Use with an overtube may result in mucosal pinching/entrapment and/or barrel detachment.

Band ligation may not be effective when applied to small varices.

Esophageal ligation devices are not intended for ligation of varices below gastroesophageal junction.

Current literature addresses management of acutely bleeding esophageal varices and does not address prophylactic use of banding.

Passing endoscope over a previously placed band may dislodge band.

Prior to assembling device, routine endoscopic examination is recommended to confirm diagnosis requiring treatment of esophageal varices or internal hemorrhoids.

It is vital that the integrity of the working channel is intact as grooves or other obstructions in the working channel can potentially cause the string to catch, resulting in band deployment difficulty.

Use of an endoscope in a sound state of repair is a prerequisite for a successful multi-band ligation procedure.

### COMPONENT PARTS

Barrel with preloaded bands and attached trigger cord, Multi-Band Ligator handle, Loading catheter, Irrigation adapter.

### SYSTEM PREPARATION

1. Examine features of handle. It has two positions which control rotation. **Firing** position (*fig. 1*) allows handle to be rotated in forward direction only. **Two-way** position (*fig. 2*) allows handle to rotate in both directions. Prior to introducing endoscope, keep handle in two-way position.
2. Insert handle into endoscope accessory channel following instructions below for appropriate endoscope.

Olympus - With rubber cap on accessory channel, lift plug on cap and insert stem of handle. (*fig. 3a*)

Pentax - Remove rubber cap from accessory channel. Place stem of handle through cap. Replace cap and inserted handle, as a unit, onto Luer lock fitting of accessory channel. (*fig. 3b*)

Fujinon - Place stem of handle through rubber cap attached to endoscope accessory channel. (*fig. 3c*)

**Note:** The irrigation adapter may be used to puncture the white self-sealing valve prior to introducing the loading catheter. (*fig. 3d*)

3. Introduce loading catheter through white seal in handle and advance in short increments until it exits tip of endoscope. (*fig. 4*) **Note:** Loading catheter is identical on both ends and may be introduced either way.
4. Attach trigger cord to hook on end of loading catheter, leaving approximately 2 cm of trigger cord between knot and hook. (*fig. 5*) Withdraw loading catheter and trigger cord up through endoscope and out through handle. Dispose of loading catheter per institutional guidelines for biohazardous medical waste.
5. Attach barrel to tip of endoscope, ensuring barrel is advanced onto tip as far as possible. (*fig. 6*) **Note:** When placing the barrel onto the distal end of the endoscope, ensure that the trigger cord does not become pinched between the barrel and the endoscope.
6. With endoscope tip straight, place trigger cord into slot on spool of handle (*fig. 7a*) and pull down until knot is seated in hole of slot. (*fig. 7b*) **Note:** Knot must be seated into hole or handle will not function properly.
7. With handle in two-way position, slowly rotate handle clockwise to wind trigger cord onto handle spool until it is taut. (*fig. 8*) **Note:** Care must be taken to avoid deploying a band while winding trigger cord.
8. Check endoscopic view. To maximize visualization, position of trigger cord may be altered by rotating barrel. **Note:** Endoscopic view broadens after each band deployment.
9. Lubricate endoscope and exterior portion of barrel. **Caution:** Do not place lubricant inside barrel. **Caution:** Do not apply alcohol to device.

Device is now ready for ligation of esophageal varices or internal hemorrhoids.

### INSTRUCTIONS FOR LIGATION OF ESOPHAGEAL VARICES

1. With handle in **two-way** position, introduce endoscope into esophagus. (*fig. 9*) After intubation, place handle in **firing** position. (*fig. 10*)
2. Visualize selected varix and aspirate it into barrel. **Caution:** Prior to band deployment, ensure endoscopist's hand is positioned on handle of device rather than endoscope controls.

3. Maintain suction and deploy band by rotating handle clockwise until band release is **felt**, indicating deployment. (fig. 10) **Note:** If band will not deploy, place handle in two-way position and loosen trigger cord slightly. Place handle in firing position and continue with procedure.
4. Release suction button of endoscope, insufflate air, then withdraw scope slightly to release ligated varix. **Note:** An irrigation adapter is provided with each device. If irrigation of accessory channel is desired to clear viewing field, attach adapter to a syringe filled with sterile water and insert into white seal of handle. Irrigate as necessary. (fig. 10)
5. Repeat ligation process as needed. **Note:** More than one ligation band for each varix may be required to control acute bleeding.
6. If more bands are required, remove endoscope and attach a new device. **Note:** An average of 3 to 4 ligation sessions may be required to obliterate varices.

#### **INSTRUCTIONS FOR LIGATION OF INTERNAL HEMORRHOIDS**

**Note:** A flexible sigmoidoscopic screening exam of hemorrhoidal site must be performed prior to assembling and use of device.

1. With handle in **two-way** position, introduce endoscope into rectum. (fig. 9)
2. Retroflex endoscope to visualize selected internal hemorrhoid. (fig. 11) After retroflexion, place handle in **firing** position (fig. 10), then aspirate internal hemorrhoid into barrel. (fig. 12) Prior to band deployment, ensure endoscopist's hand is positioned on handle of device rather than endoscope controls.
3. Maintain suction of internal hemorrhoid and deploy band by rotating handle clockwise until band release is **felt**, indicating deployment. (fig. 10) **Note:** If band will not deploy, place handle in two-way position and loosen trigger cord slightly. Place handle in firing position and continue with procedure.
4. Release suction button of endoscope, insufflate air, then advance scope slightly to release ligated hemorrhoid. **Note:** An irrigation adapter is provided with each device. If irrigation of accessory channel is desired to clear viewing field, attach adapter to a syringe filled with sterile water and insert into white seal of handle. Irrigate as necessary. (fig. 10)
5. Repeat ligation process as needed.
6. Place handle in two-way position when ligation process is completed, **then** straighten scope.
7. If more bands are required, remove endoscope and attach a new device.

#### **REMOVING THE MULTI-BAND LIGATOR**

1. Upon completion of ligation procedure, remove endoscope from patient.
2. Dismantle device as follows:

*If all the bands have been fired:*

- Remove handle and attached trigger cord from accessory channel.
- Remove barrel from endoscope tip.

*If any unfired bands remain on barrel:*

- Place handle in two-way position.
- Loosen trigger cord from spool, then remove handle from accessory channel cap.
- Detach trigger cord from handle slot.
- Remove barrel from endoscope tip, then pull trigger cord through channel and out endoscope tip.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**