EN

PEG Balloon Replacement Tube

DA

PEG ballonerstatningssonde

NL

Vervangende PEG-ballonsonde

FR

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DE

PEG-Ballon-Ersatzsonde

EL

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



ENGLISH

INTENDED USE

This device is intended to be used as a replacement gastrostomy tube through an established gastrostomy tract.

CONTENTS

- 1 Balloon Replacement Gastrostomy Tube
- 1 Twist Lock
- 1 Prefilled 10 cc Luer Syringe (14 FR)
- 2 Prefilled 10 cc Luer Syringes (18 FR, 24 FR)
- 2 pks. Water-Soluble Lubricant
- 1 Instruction Booklet
- 1 Patient Care Manual

NOTES

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those associated with placement and use of a balloon replacement gastrostomy tube include, but are not limited to: patient without an established gastrostomy tract, signs of infection such as redness, edema or purulent drainage present, sepsis, severe gastroesophageal reflux, large ascite, or diffuse inflammatory, infectious, or neoplastic disease involving the walls of the abdomen or anterior stomach, gastrointestinal obstruction or proximal small bowel fistulae.

POTENTIAL COMPLICATIONS

Those associated with placement and use of a balloon replacement gastrostomy tube include, but are not limited to: bronchopulmonary aspiration and pneumonia, respiratory distress or airway obstruction, peritonitis or septic shock, colocutaneous, gastrocolocutaneous or small bowel fistula, gastric dilatation, sigmoid intra-abdominal herniation and volvulus, persistent fistula following tube removal, esophageal injury, necrotizing fasciitis, candida cellulitis, improper placement or inability to place tube, tube dislodgment or migration, hemorrhage, and tumor metastasis.

Additional complications include, but are not limited to: pneumoperitoneum, peristomal wound infection and purulent drainage, stomal leakage, bowel obstruction, gastroesophageal reflux (GERD), and blockage or deterioration of balloon replacement tube.

PRECAUTIONS

Benefits of a balloon replacement tube to patient must be weighed against risks associated with any indwelling gastrostomy feeding tube.

During placement and use, care must be taken to avoid cutting, crimping, or damaging components.

Do not modify balloon replacement tube in any way.

Follow instructions and Patient Care Manual supplied with each kit. It is essential for Patient Care Manual to accompany patient and be explained to all people responsible for care of patient.

All PEG Systems, including Balloon Replacement Gastrostomy Tube, are radiopaque. Proper location and integrity of any internal component can be visualized by x-ray.

Balloon must be inflated with **sterile water only**. Do **not** use air or saline, as they may cause premature deflation.

Manufacturer guidelines are not intended to replace physician's recommendations.

When placing balloon replacement tube in obese patients, all anatomical structures must be identified prior to placement.

Replacement of balloon replacement tube is recommended every 29 days or at discretion of healthcare provider.

WARNINGS

Do not over inflate balloon. Feeding into over inflated balloon may result in damage to gastric mucosa, tube migration, tube failure, sepsis, infection or death.

Bolster should rest gently on skin surface. Excessive traction on feeding tube may cause premature removal, fatigue or failure of device.

INSTRUCTIONS FOR USE

Tube Placement

1. Evaluate existing stoma tract to ensure a well-established gastrostomy tract into stomach.
2. Push pivotal bolster above 10 cm mark.
3. Generously lubricate replacement tube and stoma site with **water-soluble lubricant**.
4. Gently advance tube into stoma site until balloon is completely in stomach. (See fig. 1)
5. Securely attach pre-filled syringe to access port on balloon replacement tube labeled "BALLOON" and inflate. (See fig. 2)
6. Withdraw tube gently until tension is felt as balloon rests against stomach wall. (See fig. 3)
7. Slide pivotal bolster down shaft of tube until bolster rests gently against abdominal wall.
8. Secure twist lock around bolster collar, being careful not to crimp it. (See fig. 4)
Important: Use twist lock to secure bolster to tube. This will help to prevent future migration of tube and reduce need to constantly reposition or pull on tube.
9. Dispose of any unused items per institutional guidelines for biohazardous medical waste.
10. Note centimeter marking on tube that is closest to bolster and record it on patient's chart and on patient information sheet in Patient Care Manual.
11. Record physician's instructions for feeding and administration of prescribed medications in Patient Care Manual.
12. Patient should remain NPO for 24 hours unless otherwise directed by a physician.

BALLOON SIZE CHART	
French Size	Inflated Balloon Volume (Sterile Water)
14	5 cc
18	20 cc
24	20 cc

Tube Removal

Caution: A replacement tube is desired, it must be placed immediately after removal. External tube removal may result in minimal trauma or bleeding that may require treatment.

Feeding tube may be removed in following manner:

1. Attach Luer lock syringe to balloon port on tube and deflate balloon completely. If balloon does not deflate, cut white tubing below balloon port. This will deflate balloon completely.
2. Gently remove tube from stoma tract.

Dispose of device per institutional guidelines for biohazardous medical waste. Stoma tract should heal and close within 24 hours.

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© 2014 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.