EN	Memory™ Soft Wire Baskets
CS	Košíky Memory™ z měkkého drátu
DA	Memory™ bløde trådkurve
NL	Memory™ Soft Wire baskets
FR	Paniers Memory™ en fil souple
DE	Memory™-Weichdraht-Körbchen
EL	Καλάθια από μαλακό σύρμα Memory™
HU	Memory™ puha drótkosarak
IT	Cestello Memory™ flessibile
NO	Memory™ myke trådkurver
PL	Koszyki Memory™ z miękkiego drutu
PT	Cestos de arame flexível Memory™
ES	Cestas de alambre blando Memory™
SV	Memory™ korgar med mjuk tråd

ENGLISH

INTENDED USE

This device is used for the endoscopic removal of biliary stones and foreign bodies.

NOTES

Select Memory Soft Wire Baskets are compatible with the Soehendra Lithotripter.

This device is reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to the primary endoscopic procedure to be performed in conjunction with the removal of the foreign body.

Perforation.

If this device is to be used for the removal of biliary stones, contraindications include an ampullary opening inadequate to allow for the unimpeded passage of the stone and basket.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy and ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction, respiratory depression or arrest, cardiac arrhythmia or arrest.

Others associated with basket extraction include, but are not limited to: impaction of the object, aspiration of foreign body, localized inflammation, pressure necrosis.

PRECAUTIONS

Refer to package label for minimum channel size required for this device and to determine basket compatibility.

If this device is to be used for removal of biliary stones, the following caution applies: Assessment of stone size and ampullary orifice must be made to determine necessity of sphincterotomy. If sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

WARNINGS

When retrieving a foreign body from upper GI tract, extreme caution must be exercised to avoid occlusion of airway.

Surgical intervention may be required if stone impaction and/or basket fragmentation occurs.

If a stone cannot be removed endoscopically with this basket, the Soehendra Lithotripter may be used with select Memory Soft Wire Baskets (see package label) to mechanically crush stone and aid in removal. Due to mechanical pressure generated with Soehendra Lithotripter, basket fragmentation and/or stone impaction in common bile duct may occur and require surgical intervention. Risk of basket fragmentation or stone impaction must be weighed against potential benefit of using lithotripter.

INSTRUCTIONS FOR USE

- With basket fully retracted into sheath, insert device into endoscope accessory channel. Advance device through channel in short increments until basket sheath exits endoscope. **Note:** When using a basket that is removable from sheath, basket may be removed to facilitate use of a wire guide. When sheath is in proper position, remove wire guide, then reintroduce basket through sheath.
- Endoscopically or fluoroscopically visualize extraction target and advance basket to appropriate position. Fluoroscopic visualization may be enhanced by injecting contrast via side arm fitting. Place sheath of basket beyond target before extending basket to entrap extraction target. **Note:** Tighten connector cap on handle before injection of contrast. (See fig. 1)
- Confirm desired position of basket sheath relative to target. With one hand holding side arm fitting, **gently** push forward on pin vise handle with other hand to advance basket out of sheath. **Caution:** Pulling on sheath while advancing or retracting basket may damage device, rendering it inoperable.
- Place basket around extraction target or foreign body and slowly retract handle until entrapment is achieved. **Note:** It may be necessary to move basket back and forth along target in order to capture object. **Note:** Connector cap on handle may be loosened in order to more freely manipulate basket. (See fig. 1)
- Withdraw basket into channel until foreign body is supported against endoscope tip, then withdraw endoscope from patient. Maintain endoscopic visualization to ensure continued entrapment of object. While withdrawing device from endoscope, wipe secretions from outer sheath. **Caution:** If difficulty is encountered during removal of basket from duct, moderate force may be applied by pulling on handle to manually fracture stone. If this proves unsuccessful, a compatible basket may be used with Soehendra Lithotripter to mechanically crush stone. If passage is still restricted, surgical intervention may be necessary.
- Prepare or dispose of extracted specimen per institutional guidelines.
- Upon completion of procedure, clean and prepare device for future use following cleaning instructions below. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, bends, or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.

CLEANING INSTRUCTIONS

Baskets Not Removable From Sheath

- Immediately after use, attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting on device, then irrigate.
- Extend basket, then submerge solution-filled device in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
- Beginning at handle and working toward tip, wipe device with a textured cloth soaked in enzymatic solution.
- Irrigate inner sheath by attaching a syringe containing a minimum of 20 cc of enzymatic solution to Luer lock fitting on handle of device. Repeat this step until all visible secretions are removed.
- Gently brush basket with enzymatic solution and a soft nylon bristle brush to remove all foreign matter.
- Remove device (or components, if applicable) from enzymatic solution, then thoroughly rinse inner sheath using a syringe filled with clean water. Rinse handle, outer sheath and extended (or separate) basket thoroughly with clean, running water.
- Ultrasonically clean device (or components, if applicable) with basket extended, for a minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.
- Remove device (or components, if applicable) from ultrasonic cleaner, thoroughly rinse inner sheath with a syringe containing a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse external portions of device and basket with clean, running water.
- Gently wipe external portions of device dry with a lint-free cloth. Hang device with basket extended to air dry lumen and basket, or use medical grade (oil-free) forced air through Luer lock fitting on handle.
- Package device for sterilization following AAMI recommended practices.
- Sterilize device using set points for desired sterilization method shown in tables.

Baskets Removable From Sheath

- Immediately after use, disconnect handle from catheter. While withdrawing basket, gently wipe secretions from wires. Attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting on catheter, then irrigate.
- Submerge solution-filled catheter and basket in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.

REFER TO STEPS 3-9 OF "BASKETS NOT REMOVABLE FROM SHEATH" ABOVE, THEN CONTINUE WITH STEP 3 BELOW:

- Gently re-assemble basket and sheath.
- Package device for sterilization following AAMI recommended practices.
- Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4°C (130°F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4°C (130°F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilization, the product should be taken out of service immediately.

AUTOClave CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	132.2°C (270°F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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
“Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers”, AAMI Technical Information Report, No. 12 - 1994.

“Sterilization, Part 1 - Good Hospital Practices”, AAMI Standards and Recommended Practices, Volume 1.1 - 1995.

RxOnly

STERILEEO

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

EC

REP

Cook Ireland Ltd.
O’Halloran Road
National Technology Park
Limerick
Ireland