
EN 14	<b>Percutaneous Endoscopic Gastrostomy Systems FLOW-J / PEG-J Gastro-Jejunal Feeding Tube</b>
CS 16	<b>Perkutánní endoskopické gastrostomické systémy Gastrojejunální výživová sonda FLOW-J / PEG-J</b>
DA 19	<b>Perkutane endoskopiske gastrostomisystemer FLOW-J-/PEG-J-gastrojejunal ernæringssonde</b>
NL 22	<b>Percutane endoscopische gastrostomiesystemen – FLOW-J/PEG-J-gastrojejunostomiesonde</b>
FR 25	<b>Systèmes pour gastrostomie endoscopique percutanée Sonde d'alimentation gastrojéjunale FLOW-J/PEG-J</b>
DE 28	<b>Perkutane endoskopische Gastrostomiesysteme Gastro-jejunale Ernährungssonde FLOW-J/PEG-J</b>
EL 32	<b>Γαστρονησιδικός σωλήνας σίτισης συστημάτων διαδερμικής ενδοσκοπικής γαστροστομίας FLOW-J/PEG-J</b>
HU 35	<b>Perkután endoszkópos gastrostomiás rendszerek FLOW-J / PEG-J gastrojejunalis táplálótubus</b>
IT 38	<b>Sistemi per gastrostomia endoscopica percutanea FLOW-J/Sondino gastrodigiuale PEG-J</b>
NO 41	<b>Perkutane endoskopiske gastrostomisystemer - FLOW-J / PEG-J gastrojejunal ernæringssonde</b>
PL 44	<b>Systemy do przezskórnej endoskopowej gastrostomii Zgłębnik do żywienia żołądkowo-czczy FLOW-J / PEG-J</b>
PT 47	<b>Sistemas de gastrostomia endoscópica percutânea FLOW-J/tubo de alimentação gastrojejunal PEG-J</b>
ES 50	<b>Sonda de alimentación gastroyeyunal FLOW-J y PEG-J para sistemas de gastrostomía endoscópica percutánea</b>
SV 54	<b>Perkutana endoskopiska gastrostomisystem FLOW-J/PEG-J gastrojejunal matnings slang</b>

**PEG ONLY**

1.

## ENGLISH

### PEG ONLY

1. Antrum
2. Feeding Adapter
3. Fundus
4. Pivotal Bolster
5. Gastric Lumen
6. Cut at X-mark
7. G-Tube

## ČESKY

### POUZE PEG

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**INTENDED USE**

This device is used for percutaneous endoscopic placement to provide enteral nutrition to patients requiring nutritional support.

**CONTENTS**

1 Gastro-Jejunal Feeding Tube, 1 - 240 cm Wire Guide, 1 Twist Lock, 1 Cable Tie, 1 Air Plug, 2 pks. Water-Soluble Lubricant (USA only), 1 Instruction Manual, 1 Patient Care Manual.

**NOTES**

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

System also allows decompression as prescribed by physician.

**CONTRAINDICATIONS**

Those associated with placement and use of a PEG-J tube include, but are not limited to: sepsis, severe gastroesophageal reflux, large ascites, or diffuse inflammatory, infectious, or neoplastic disease involving walls of abdomen or anterior stomach, gastrointestinal obstruction or proximal small bowel fistulae.

**WARNINGS**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Excessive traction on feeding tube may cause premature removal, fatigue or failure of device.

**POTENTIAL COMPLICATIONS**

Those associated with placement and use of PEG-J tube include, but are not limited to: bronchopulmonary aspiration and pneumonia, respiratory distress or airway obstruction, peritonitis or septic shock, colcutaneous, gastrocolocutaneous or small bowel fistula, gastric dilatation, sigmoid intra-abdominal herniation and volvulus, persistent fistula following PEG-J removal, esophageal injury, necrotizing fasciitis, candida cellulitis, improper placement or inability to place PEG-J tube, tube dislodgment or migration, hemorrhage, and tumor metastasis.

Others include, but are not limited to: pneumoperitoneum, peristomal wound infection and purulent drainage, stomal leakage, bowel obstruction, gastroesophageal reflux (GERD), and blockage or deterioration of PEG-J tube.

**PRECAUTIONS**

Benefit of a PEG-J tube to patient must be weighed against risks associated with any indwelling gastrostomy feeding tube.

During placement and use, care must be taken to avoid cutting, crimping, or damaging components.

Do not modify PEG, J-tube or adapters in any way.

Follow the instructions and the Patient Care Manual supplied with each kit. It is essential for the Patient Care Manual to accompany the patient and be explained to all people responsible for the care of the patient.

All PEG Systems, including PEG-J and FLOW-J, are radiopaque. Proper location and integrity of any internal component can be visualized by x-ray.

Proper positioning of gastro-jejunal feeding tube must be verified by x-ray.

Device is only compatible with appropriate Cook PEG tube. This device is not compatible with Cook Replacement PEG Tube.

Manufacturer guidelines are not intended to replace physician's recommendations.

When placing PEG-J tube in obese patients, all anatomical structures must be identified prior to placement.

When placing PEG-J tube, observe all institutional guidelines regarding gastroscopy, including removal of dentures.

PEG-J or FLOW-J replacement is recommended every **three** months or as recommended by physician.

## INSTRUCTIONS FOR USE

### Tube Placement

1. After introducing gastroscope, insufflate stomach and examine mucosa. Determine if mucosa is free of ulcerations or bleeding before proceeding. **Important:** This gastrojejunal feeding tube (J-tube) is fully radiopaque and can only be used with appropriate Cook PEG tube. Ideally, PEG tube should be placed in antrum or low in body of stomach to ease insertion.
2. Remove feeding adapter and twist lock or cable tie securing external bolster from previously placed Cook PEG tube. X mark on PEG tube should be at end. If not, tube should be cut at X mark. (See fig.1)
3. Place **non-spiked biopsy or alligator forceps** through gastroscope and push forceps through PEG tube from stomach until it is seen exiting tube.
4. **Generously lubricate** wire guide with **water-soluble lubricant**, especially at tip.
5. Place air plug over wire guide and grasp floppy tip end of wire guide with forceps. (See fig. 2)
6. Pull wire guide and forceps back into tube and place air plug securely in end of PEG tube to minimize air leakage. (See fig. 3)
7. Insufflate stomach to provide better visualization of pylorus.
8. Advance gastroscope with wire guide and forceps through pylorus and into second portion of duodenum.
9. **Straighten wire guide by holding tension at both ends.** Advance wire guide and forceps to third portion of duodenum.
10. **Generously lubricate** J-tube **with water-soluble lubricant**, especially at tip.
11. Remove air plug from G-tube, then place J-tube over wire guide. (See fig. 4)
12. While maintaining tension on wire guide, slowly advance J-tube over wire guide. Continue advancement until it dislodges forceps from wire. **Note:** J-tube should move distally when it dislodges forceps from wire.
13. Remove forceps from scope and plug feeding adapter firmly into PEG tube. (See fig. 5)
14. Endoscopically verify J-tube position.
15. **Aspirate air, then remove scope in a rotating fashion** to minimize dragging effect on J-tube.
16. Verify that J-tube has not backed through pylorus into stomach. (See fig. 6) **Note:** A cable tie may be placed around proximal gastrostomy tube at feeding adapter to avoid inadvertent removal of J-tube from PEG.
17. Carefully remove wire guide while observing endoscopically. Once again verify that J-tube has not fallen back into stomach, then remove gastroscope.
18. Using a 60 cc piston irrigation syringe, flush 50 cc of water through **yellow** J-tube port on feeding adapter. Tube should flush easily if properly placed. **Note:** Red G-tube port on feeding adapter is used to access stomach.
19. Replace twist lock or cable tie on external bolster of PEG tube, being careful not to crimp tube. **Important:** Use twist lock or cable tie to secure pivotal bolster to tube. This will help prevent future migration of tube and reduce need to constantly reposition or pull on tube. **Warning:** Excessive traction on gastric feeding tube may cause premature removal, fatigue or failure of device.

- Record physician's instructions for feeding and administration of prescribed medications in Gastro-Jejunal Patient Care Manual.
- Perform abdominal x-ray following placement of J-tube to verify proper tube position.

### Tube Removal

**Note:** Prior to removal, we recommend purchase of a new G-J tube, or universal adapter to plug into Cook PEG tube if only PEG tube feeding is desired.

**Caution:** If replacement tube is desired, it must be placed immediately after removal. External tube removal may result in minimal trauma or bleeding that may require treatment.

#### G-J Tube Removal

- Unplug G-J adapter from PEG tube, then gently pull G-J tube to remove it completely from PEG tube.
- Dispose of G-J tube per institutional guidelines for biohazardous medical waste.
- To resume PEG feeding, place universal adapter or bolus adapter into PEG and resume feeding as prescribed.

#### PEG Tube Removal (External/Traction Method)

The PEG tube is designed for removal using the external/traction method shown below. If this method is not possible, another method of removal such as endoscopic or surgical should be utilized.

Feeding tube may be removed without need for endoscopy in following manner:

- Remove twist lock or cable tie and bolster from PEG tube.
- Grasp gastric feeding tube near stoma site. While slowly rotating tube, gently push 1-2 cm of tube into stomach to separate tube from stoma tract. **Warning:** If tube does not rotate freely within tract, do not attempt to use traction as a method of removal.
- Hold gastric tube near stoma site and apply counter pressure by placing fingers of other hand around base of tube.
- Loosely cover stoma site with either a towel or a drape.
- Holding gastric feeding tube straight, apply steady traction to tube until internal dome emerges through abdominal wall. **Caution:** Tube must be pulled straight out of stoma tract.
- Stoma tract should heal and close within 24 hours.

**Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.**

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**