EN

Howell Biliary Introducer Needle / Brush / Forceps

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ENGLISH

INTENDED USE

The Howell Biliary Introducer is used with the Howell Biliary Introducer Needle for aspiration biopsy in the biliary system, and if desired to introduce the Howell Biliary Introducer Brush to collect cells from the biliary system or the Howell Biliary Introducer Forceps to obtain endoscopic mucosal tissue biopsies and/or foreign body retrieval.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Forceps are supplied non-sterile and reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

Coagulopathy.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

If introducer is used in conjunction with needle, needle must be fully retracted into catheter prior to introduction, advancement, and withdrawal from introducer. Failure to retract needle may result in damage to introducer.

Reusability of forceps depends, in large part, on care of device by user. Factors involved in prolonging life of forceps include, but are not limited to: advancing through accessory channel in short increments, maintaining handle alignment with accessory channel, gently withdrawing device from accessory channel, avoiding loops and/or kinks in cable and use of excessive force, coiling sheath in a minimum of 8-inch (20 cm) diameter, thorough cleaning following instructions included in this booklet.

INSTRUCTIONS FOR USE

Howell Biliary Introducer

1. With endoscope elevator open, introduce device into accessory channel of endoscope. If device is advanced over a pre-positioned .035 inch wire guide, flush wire guide lumen with sterile water prior to introduction of device to facilitate advancement.
2. Advance device in small increments until radiopaque Capture™ side ramp is endoscopically visualized exiting scope. **Note:** Capture side ramp must be completely out of endoscope in order to ensure proper operating condition of device.
3. Under fluoroscopic guidance, slowly advance introducer into duct to be biopsied. **Note:** Introducer has two radiopaque bands to assist with fluoroscopic visualization of catheter position. Selected biopsy device exits introducer at side ramp.
4. Fluoroscopically confirm position of introducer in relation to area to be biopsied. Position radiopaque marker indicating location of ramp below area to be biopsied. **Note:** Contrast may be injected through wire guide lumen. If a wire guide is in place, contrast may be injected using a side arm adapter placed over wire guide and attached to Luer lock fitting on wire guide lumen.
5. When introducer is in desired position, a compatible biopsy device may be introduced.

Note: Introducer should remain in place to facilitate the exchange of compatible devices listed below. **Note:** It is recommended that any brushing procedure be performed last in sequence as it is necessary to remove introducer and brush together.

INSTRUCTIONS FOR USE

Howell Biliary Introducer Needle

1. Ensure that needle is completely retracted into sheath, then advance needle catheter through proximal hub on introducer until tip of catheter reaches side ramp.
2. Gently advance needle catheter out of introducer until it is even with tip of introducer and targeted area.
3. Lock hub on needle handle to introducer hub by pressing into place.
4. Fluoroscopically confirm desired needle position using metal tip on outer catheter as a guide.
5. Unlock thumbscrew ring on handle of needle catheter, then slowly advance handle to advance needle out of catheter. **Note:** Depending on composition of area to be sampled, needle may not fully extend out of catheter. Gentle movement of needle will assist in full extension.
6. Remove stylet wire from Luer lock fitting on needle handle, then attach a 20 cc syringe to Luer lock fitting and aspirate. Continue aspiration while withdrawing needle from biopsy site.
7. When aspiration biopsy is complete, completely withdraw needle into catheter by pulling back on needle handle. Tighten thumbscrew ring on needle handle to lock needle in place.
8. Withdraw needle from introducer. **Note:** When removing introducer assembly from endoscope, wire guide may be left in place to facilitate advancement of compatible wire guided accessories.
9. Prepare specimen per institutional guidelines.

Upon completion of procedure, dispose of introducer and needle per institutional guidelines for biohazardous medical waste.

INSTRUCTIONS FOR USE

Howell Biliary Introducer Brush

1. Advance brush through proximal hub on introducer until it reaches side ramp.
2. Gently advance brush out of introducer until it is even with tip of introducer and area targeted for biopsy.
3. Under fluoroscopic monitoring, brush targeted area to obtain adequate cellular material.
4. Retract brush into introducer.
5. Withdraw introducer and brush as a unit from endoscope. **Note:** When removing introducer assembly from endoscope, wire guide may be left in place to facilitate advancement of compatible wire guided accessories.
6. Prepare specimen per institutional guidelines.

Upon completion of procedure, dispose of brush and introducer per institutional guidelines for biohazardous medical waste.

INSTRUCTIONS FOR USE

Howell Biliary Introducer Forceps

1. With cups closed, advance forceps through proximal hub of introducer until tip of forceps reaches side ramp.
2. Gently advance forceps out of introducer until it is even with tip of introducer and area targeted for biopsy.
3. Under fluoroscopic monitoring, open forceps cups and advance into targeted area. Using slight pressure on handle, close cups to obtain a sample. **Note:** It is not necessary to apply excessive pressure to cleanly excise tissue.
4. Maintain gentle handle pressure to keep cups closed and gently withdraw forceps from site and retract into introducer.
5. Continue to apply slight pressure on handle and withdraw forceps. **Note:** Repeat steps 1 - 4 until adequate samples are obtained. **Note:** When removing introducer assembly from endoscope, wire guide may be left in place to facilitate advancement of compatible wire guided accessories.
6. Prepare specimen per institutional guidelines.
7. **Upon completion of procedure, dispose of introducer per institutional guidelines for biohazardous medical waste.** Clean and prepare forceps for future use following cleaning instructions that follow. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, bends or breaks exist, do not use forceps. Dispose of device per institutional guidelines for biohazardous medical waste.

CLEANING INSTRUCTIONS

Howell Biliary Introducer Forceps

1. **Immediately** after use, rinse forceps with clean, running water.
2. Submerge forceps in a solution of warm water and enzymatic cleaner mixed according to enzymatic cleaner manufacturer's instructions. Soak for a minimum of 20 minutes. **Note:** When submerging forceps, insert cups into liquid first, then slowly advance cable into fluid.
3. Beginning at handle and working toward tip, wipe device with a textured cloth soaked in enzymatic solution. **Caution:** Stretching cable when wiping forceps will adversely affect performance of device.
4. With cups open, gently brush all cup surfaces and hinges with a soft nylon bristle brush to remove all foreign matter.
5. Remove device from enzymatic solution. Beginning at handle, rinse handle, cable and cups with clean, warm, running water. **Note: Thorough physical cleaning to remove all foreign matter prior to ultrasonic cleaning is essential in order to prolong performance characteristics of forceps.**
6. Ultrasonically clean device with cups opened for a minimum of 15 minutes at 25°C (77°F) per ultrasonic cleaner manufacturer's instructions. **Caution: Do not allow forceps to dry between physical cleaning and placing in ultrasonic cleaner.**
7. After removing forceps from ultrasonic cleaner, **thoroughly** rinse device with clean water to remove all residue.
8. Rinse forceps in 70% ethyl alcohol to facilitate drying of inner drive cable.
9. Gently wipe external portions of device dry with a lint-free cloth. Hang device to air dry or use medical grade (oil-free) forced air along cable and around handle and cups.
10. Lightly lubricate hinges of forceps with medical grade silicone lubricant. **Note:** Proper lubrication of forceps after each use will help maintain performance characteristics of device.
11. Loosely coil forceps in a minimum 8-inch (20 cm) diameter and package for sterilization following AAMI recommended practices.

12. Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4°C (130°F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4°C (130°F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilization, the product should be taken out of service immediately.

AUTOClave CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	132.2°C (270°F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.