EN
2

Disposable Hot Biopsy Forceps

Instructions for Use

FR
5

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Mode d'emploi

PT
8

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ENGLISH

INTENDED USE

This device is used endoscopically in conjunction with monopolar electrosurgical current to obtain gastrointestinal mucosal tissue biopsies and for removal of sessile polyps.

This device is indicated for adult use only.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired biopsy or polypectomy site.

Contraindications to gastrointestinal mucosal biopsy and polypectomy include, but are not limited to: coagulopathy and insufficiently prepped bowel.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to perforation, bleeding or hemorrhage, aspiration, fever, infection, allergic reaction to medication,

hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Potential complications with endoscopic mucosal biopsy or polypectomy include, but are not limited to: transmural burns, thermal injury to the patient, and explosion.

USER QUALIFICATIONS

This device should only be used by or under the supervision of physicians who have been trained in gastrointestinal mucosal biopsy and/or polypectomy procedures.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

The endoscope must remain as straight as possible when inserting or withdrawing the device.

The device's cups must remain closed during introduction into, advancement through, and removal from the

endoscope. If the cups are open, damage to the device and/or endoscope may occur.

Exercising the handle while the device is coiled may result in damage to the device.

Mucosal biopsy and polypectomy should be performed under direct endoscopic visualization.

Switch the electrosurgical unit to the "off" position when it is not in use.

WARNINGS

Do not use this device for any purpose other than the stated intended use.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse the device may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Any electrosurgical device constitutes a potential electrical hazard to the patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

The electrosurgical effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the device.

The device should never be connected to the active cord prior to inserting the device through the endoscope to avoid injury to the patient and or equipment resulting from improper electrical circuit grounding.

Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper selection, placement, and utilization of the patient return electrode. Ensure that a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

For patients with electrically conductive implants, monopolar diathermy or electrosurgical cautery can result in localized thermal injury at the implant site.

Do not use this device with an output higher than the rated voltage: 1.7 kVp-p (850 Vp). This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, device and/or active cord.

Do not use an active cord with a maximum rated voltage less than 1.7kVp-p (850 Vp). This could cause thermal injury to the patient, operator, or assistant, and could also damage the endoscope, device, and/or active cord.

This device must be used with an endoscope compatible with high frequency (HF) equipment. If the device is used with an endoscope not compatible with HF electrosurgery, a burn to the patient or the operator may occur.

If cutting is not observed at normal operating settings, this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.

The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.

Skin-to-skin contact should be avoided (for example between the arms and body of the patient) using dry gauze.

Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.

Avoid contact between the active cords and the patient's body, or any other electrodes.

This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electric shock in compliance with IEC 60601-1). Use of generator other than type BF or CF could cause electric shock to the patient, operator, or assistant.

For recommendations on current settings for gastrointestinal mucosal biopsy and/or polypectomy, refer to the monopolar generator's recommended power settings for gastrointestinal mucosal biopsy and/or polypectomy.

This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere or in the presence of explosive gases. Doing so may result in an explosion.

In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electrosurgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the device in these patients.

Using accessories others than those recommended by the electrosurgical generator manufacturer may result in decreased immunity or increased emissions.

This device must be used with an endoscope that is intended for use with high frequency (HF) electrosurgical devices. If the device is used with an endoscope not

compatible with high frequency (HF) electrosurgery, a burn to the patient or the operator may occur.

To avoid neuromuscular stimulation and minimize electric arc, activate the device only when the cups are in contact with the tissue to be cauterized.

Failure to isolate tissue to be biopsied or polyp to be removed by pulling it away from mucosal wall may result in fulguration of normal mucosa and/or perforation.

These single-use biopsy forceps should only be used to biopsy tissue where possible bleeding or hemorrhage will not present a danger for patients. Adequate plans for management of potential bleeding or hemorrhage and appropriate airway management should be in place.

Additional warnings may be listed in procedural steps below.

TECHNICAL DESCRIPTION

Figure 1 and Figure 2 summarize the device's characteristics.

Refer to the package label for device dimensions, and minimum endoscope compatibility.

This device must only be used with an active cord compatible with a 3 mm connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

Maximum rated input voltage for this device is 1.7 kVp-p (850 Vp).

Device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6° F).

Store in a dry location. Recommended storage temperature is 22°C (71.6°F). Device may be exposed to temperatures up to 50°C (122°F) during transportation.

Components marked with an asterisk (*) denote applied parts (the part of the device that comes into contact with the patient, as defined by IEC60601-1), see figs. 1 and 2.

SYSTEM PREPARATION

WARNING

No modification of this equipment is allowed

For recommendations on current settings for endoscopic mucosal biopsy and/or polypectomy, refer to the monopolar generator's recommended power settings for endoscopic mucosal biopsy and/or polypectomy.

Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

1. Upon removing the device from the package, uncoil the device. Remove the protective cap from the device's cups. Open and close the cups, verifying smooth handle operation and cup action.
2. With the electrosurgical unit off, prepare the equipment.

INSTRUCTIONS FOR USE

1. Endoscopically visualize the area to be biopsied or the polyp to be removed.

WARNING

Do not use a device that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

The device should never be connected to the active cord prior to inserting the device through the endoscope to avoid injury to the patient and or equipment resulting from improper electrical circuit grounding.

The device should not be energized prior to performing gastrointestinal mucosal biopsy and/or polypectomy. Energizing the device prior to use will cause premature fatigue of the device and will compromise the device's integrity.

To avoid unintended disconnection during procedure and damage to the device, device should only be used with compatible active cords (Reference Technical Description for compatibility).

2. With the cups closed, insert the device's cups into the accessory channel of the endoscope.
3. The active cord fittings should fit snugly into both the device handle and the electrosurgical unit. Following the electrosurgical unit manufacturer's instructions for settings, verify the desired settings and then activate electrosurgical unit. **Note:** Maximum rated input voltage for this device is 1.7 kVp-p (850 Vp) in coagulation mode only.
4. Advance the device in 1-2 cm increments until the device's cups are visualized exiting endoscope. Note: Keep the end of the device that is extending from accessory channel straight at all times. Allowing the device to hang from the accessory channel may cause damage to the device.

WARNING

Verify that the device's cups have exited the endoscope by visualizing them on the endoscope monitor. Failure to confirm the device's cups have fully exited endoscope may result in contact between the device and the endoscope while current is applied. Contact of the device's cups with the endoscope may result in grounding, injury to patient and/or operator, as well as damage to the endoscope and/or device.

5. Once the device's cups are visualized exiting the endoscope, advance the device to the desired biopsy

site or polyp, then open the cups and advance them into the tissue to be biopsied or the polyp to be removed.

WARNING

User must maintain contact with tissue when applying electrocautery current. Failure to do so may result in increased current resulting in a broken device, damage to the endoscope, and/or patient injury.

When applying current, ensure the device's cups do not come in contact with the endoscope. Contact of the device's cups with the endoscope may result in grounding, injury to the patient and/or operator, as well as damage to the endoscope and/or device.

6. Using slight pressure on the handle, close the device' cups around the tissue or polyp. Isolate the tissue by gently pulling away from the mucosal wall.
7. Maintain gentle handle pressure to keep the cups closed and gently withdraw the device from the site.
8. Continue to apply slight pressure on the handle and withdraw the device from the channel. While withdrawing the device from the endoscope, wipe excess secretions from cable.
9. Prepare the specimen for examination per institutional guidelines.
10. Upon completion of the procedure, turn the electrosurgical unit off. Disconnect the active cord from the handle and from the electrosurgical unit.

Upon completion of the procedure, dispose of the device(s) per institutional guidelines for biohazardous medical waste.

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A symbol glossary can be found at cookmedical.com/symbol-glossary

Pour un glossaire des symboles, consulter le site Web cookmedical.com/symbol-glossary • En cookmedical.com/symbol-glossary puede consultarse un glosario de símbolos • En symbolordlista.finns.på.cookmedical.com/symbol-glossary

If symbol appears on product label, X = maximum insertion portion width

Si le symbole est visible sur l'étiquette du produit, X = largeur maximale de la section d'insertion • Se o símbolo aparecer no rótulo do produto, X = largura máxima da parte introduzida • Si el símbolo aparece en la etiqueta del producto, X = anchura máxima de la parte que se va a insertar

If symbol appears on product label, X = minimum accessory channel

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If symbol appears on product label, X = wire guide compatibility

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If symbol appears on product label, X = quantity per box

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This device is intended
for single use only.

Wilson-Cook Medical, Inc.

4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

© COOK 2020
COOKMEDICAL.COM