EN

Howell Biliary Aspiration Needle

DA

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ENGLISH

INTENDED USE

This device is used for aspiration biopsy in the biliary system.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Needle must be fully retracted into sheath prior to introduction into, advancement through or withdrawal from endoscope. Failure to retract needle may result in damage to endoscope.

INSTRUCTIONS FOR USE

1. Uncoil and straighten catheter completely.
 2. Prior to inserting device into endoscope, retract needle completely into sheath. Ensure needle fully retracts into sheath and that no sharp tip extends beyond ball tip.
 3. Introduce device into accessory channel of endoscope and advance in short increments to desired biopsy site. Do not straighten pre-bend more than is necessary for insertion down channel. **Note:** A large-bore endoscope is preferred to minimize straightening.
 4. Position endoscope tip away from and slightly above ampulla.
 5. Place ball tip, with needle retracted, on lower lip of biliary orifice. Advance catheter into cannulating position. Proceed with cannulation of biliary system.
 6. Once cannulation is achieved, advance ball tip to appropriate position for aspiration biopsy.
 7. Carefully advance needle out of sheath. With needle in desired location, remove stylet from Luer lock fitting on handle.
 8. Attach a 20 cc syringe to Luer lock fitting, then aspirate. Continue aspiration and gently move needle in small increments back and forth within biopsy site. **Note:** Several biopsy locations may be chosen and aspiration repeated without withdrawing needle from bile duct. Between each placement of needle for biopsy, suction should be released.
 9. When aspiration is complete, withdraw needle from biopsy site and into ball tip by pulling back on Luer slip fitting, then withdraw device from endoscope.
 10. Prepare specimen for examination per institutional guidelines.
- Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.**

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.