EN	Reusable Forceps
CS	Kleště k opakovanému použití
DA	Genbrugstang
NL	Herbruikbare forceps
FR	Pince réutilisable
DE	Wiederverwendbare Zange
EL	Επαναχρησιμοποιήσιμη λαβίδα
HU	Többször használatos csipesz
IT	Pinza riutilizzabile
NO	Gjenanvendelig tang
PL	Kleszczyki wielokrotnego użytku
PT	Pinça reutilizável
ES	Pinzas reutilizables
SV	Flergångstång

18908/0411

ENGLISH

INTENDED USE

This device is used to obtain endoscopic mucosal tissue biopsies. This device may also be used for foreign body retrieval.

NOTES

This device is non-sterile and reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to the primary endoscopic procedure to be performed in gaining access to the desired biopsy or retrieval site.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Endoscope should remain as straight as possible when inserting or withdrawing forceps.

Forceps cups should be closed during introduction into, advancement through, and removal from endoscope. If cups are open, damage to forceps and endoscope may occur.

Gentle pressure should be used when operating handle of forceps. Excessive pressure will result in rigidity of forceps which may damage forceps and/or endoscope. Reusability of a device depends in large part on care of device by user. Factors involved in prolonging life of device include, but are not limited to: advancing through endoscope accessory channel in short increments, maintaining handle in alignment with endoscope accessory channel, gently withdrawing device from channel, avoiding loops and/or kinks in cable and use of excessive force, as well as thorough cleaning following instructions included in this booklet. This device should never be coiled in less than an 8 inch (20 cm) diameter.

WARNING

When retrieving a foreign body from upper GI tract, care must be taken to avoid occlusion of airway.

PRODUCT INSPECTION

Perform prior to each use.

1. Uncoil forceps beginning at handle and working toward cups, exercising care not to stretch cable. Open and close cups, verifying smooth handle operation and appropriate cup action. Become familiar with amount of handle movement required to operate cups. If any irregularities are noted, do not use. **Note:** Exercising handle while device is coiled may result in damage to forceps.

INSTRUCTIONS FOR USE

1. After endoscopically visualizing area to be biopsied or object to be retrieved, select forceps compatible with endoscope being used.
2. With cups closed, insert device into accessory channel of endoscope.
3. Advance, in small increments, until device is visualized exiting endoscope.
Note: If resistance is met while advancing forceps, straighten endoscope tip slightly. **DO NOT** force forceps through endoscope. If endoscope has an accessory elevator and resistance is met, open elevator to allow forceps to pass. Use elevator to position the forceps.
4. After advancing forceps to desired biopsy or retrieval site, open cups and advance device into tissue to be biopsied or object to be retrieved.
5. Using slight pressure on handle, close forceps around tissue or foreign body. Maintain gentle handle pressure to keep cups closed while gently withdrawing forceps from site. **Note:** It is not necessary to apply excessive pressure to cleanly excise tissue.
6. Continue to apply slight pressure on handle while withdrawing forceps from channel. If retrieving a foreign body, gently withdraw forceps until object is against tip of endoscope. Slowly withdraw endoscope, endoscopically monitoring foreign body. While withdrawing forceps from endoscope, wipe excess secretions from cable. **Note:** If endoscope has an elevator, open it before withdrawing forceps. If resistance is met while withdrawing forceps, straighten endoscope tip. **DO NOT** apply excessive force when removing forceps as damage to forceps and/or endoscope may occur. **Caution:** If forceps fail to close, slowly pull forceps cups against channel opening. Remove endoscope and forceps as a unit. Manually close forceps cups and withdraw device from endoscope.
7. Prepare specimen for examination per institutional guidelines.
8. If additional biopsies are anticipated during procedure, coil device in a minimum 8-inch (20 cm) diameter and place in a basin of clean water between samplings. **Note:** Allowing secretions to dry on inner drive wire or in coils of cable will compromise performance characteristics of device and diminish life of forceps.

CLEANING INSTRUCTIONS

Please follow carefully to ensure proper reusability of device.

1. Immediately after use, submerge device in a solution of warm water and enzymatic cleaner which has been mixed according to enzymatic cleaner manufacturer's instructions. Soak for a minimum of 20 minutes. **Note:** When submerging forceps, insert cups into liquid first and slowly advance cable into fluid.
2. Beginning at handle and working toward tip, wipe device with a textured cloth soaked in enzymatic solution. **Caution:** Stretching cable when wiping device will adversely affect performance of forceps.
3. With cups open, gently brush all cup surfaces and hinges with a soft nylon bristle brush to remove all foreign matter.
4. Remove device from enzymatic solution. Beginning at handle, rinse handle, cable and cups with copious amounts of clean, warm running water. **Note: Thorough physical cleaning to remove all foreign matter prior to ultrasonic cleaning is essential to prolong performance characteristics of forceps.**
5. Ultrasonically clean device, with cups open, for a minimum of 15 minutes at 25 °C (77 °F) per ultrasonic cleaner manufacturer's instructions. **Caution: Do not allow forceps to dry between physical cleaning and placing in ultrasonic cleaner.**
6. After removing device from ultrasonic cleaner, **thoroughly** rinse device with copious amounts of clean water to remove all residue.
7. Rinse device in 70% ethyl alcohol to facilitate drying of inner drive cable.
8. Gently wipe external portions of device dry with a lint-free cloth. Air dry forceps by hanging for a minimum of three hours, or by using medical grade (oil free) forced air along cable and around handle and cups.
9. Lightly lubricate forceps hinges with medical grade silicone lubricant. **Note:** Proper lubrication of forceps after each use will help maintain performance characteristics of device.
10. Loosely coil forceps in a minimum 8-inch (20 cm) diameter and package for sterilization following AAMI recommended practices.
11. Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4 °C (130 °F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4 °C (130 °F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilization, the product should be taken out of service immediately.

Note: AAMI stands for Association for the Advancement of Medical Instrumentation.

AUTOClave CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	134 °C (273.2 °F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber or exceed the usable chamber space.

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Obs! AAMI betyder Association for the Advancement of Medical Instrumentation.

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"Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers", *AAMI Technical Information Report*, No 12 - 1994.

"Sterilization, Part 1 - Good Hospital Practices", *AAMI Standards and Recommended Practices*, Volume 1.1 - 1995.

RxOnly

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland