EN
DA
NL
FR
DE
EL
IT
NO
PT
ES
SV

**Soehendra® Stent Retriever /
Soehendra Universal Catheter**

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ENGLISH

INTENDED USE

These devices are used in conjunction to cannulate and remove stents from the biliary and pancreatic ducts while maintaining wire guide placement.

NOTES

Soehendra Stent Retriever is supplied sterile and intended for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Soehendra Universal Catheter is supplied sterile and is reusable if the device integrity is intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Not compatible with coated wire guides.

Refer to package label for minimum channel size required for this device.

Reusability of a device depends in large part on care of device by user. Factors involved in prolonging life of Soehendra Universal Catheter include, but are not limited to: advancing through endoscope accessory channel in short increments, gently withdrawing device from channel, avoiding loops and/or kinks in catheter, thorough cleaning following instructions included in this booklet.

INSTRUCTIONS FOR USE

- Note:** If product package has been opened inadvertently and device not used, device may be resterilized **one time only**. Refer to table at end of these instructions for recommended ETO sterilization cycle parameters. Package device for sterilization following AAMI recommended practices.
- Introduce and advance catheter in small increments through endoscope accessory channel.
- Advance Cook's Standard Wire Guide through catheter in small increments until wire guide is completely through stent and well into duct.
- Remove catheter, then clean and prepare it for future use following cleaning instructions that follow. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, bends or breaks exist, do not use device. Dispose of device per institutional guidelines for biohazardous medical waste.
- Introduce and advance threaded end of stent retriever over pre-positioned wire guide until it reaches stent. *(See fig. 1)*
- Keeping end of stent retriever that is outside of scope straight, turn stent retriever handle clockwise until threaded end completely connects to stent. *(See fig. 2)* **Note:** Maximum thread engagement is recommended.
- Once stent is securely connected to retriever, open elevator of endoscope and slowly pull stent and stent retriever out of channel making sure wire guide is left in place. **Note:** Dislodgement of wire guide from duct can be prevented by simultaneously pushing wire guide and pulling stent and stent retriever out of endoscope. **Caution:** If stent detaches from stent retriever while in channel, push stent down scope and into elevator with stent retriever. Re-thread stent onto retriever, then slowly pull both out of endoscope.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

UNIVERSAL CATHETER CLEANING INSTRUCTIONS

Caution: When handling, coil catheter in a minimum of 8-inch diameter to prevent product damage.

- While withdrawing catheter from endoscope, wipe excess secretions from outer sheath.
- Immediately after use, attach syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting on catheter, then irrigate inner lumen of catheter.
- Submerge solution-filled catheter in solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
- Beginning at Luer lock fitting and working toward tip, wipe catheter with a textured cloth soaked in enzymatic solution.
- Irrigate inner lumen of catheter by attaching syringe containing a minimum of 20 cc of enzymatic solution to Luer lock fitting on catheter. Repeat this step until all visible secretions are removed.
- Remove catheter from enzymatic solution, then thoroughly rinse inner lumen of catheter using syringe filled with clean water. Rinse outer portions of catheter thoroughly with clean, running water.
- Ultrasonically clean catheter for a minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.
- Remove catheter from ultrasonic cleaner, then thoroughly rinse inner lumen with a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse outer portions of catheter with clean, running water.
- Gently wipe external portions of catheter dry with lint-free cloth. Hang device to air dry inner lumen, or use medical grade (oil free) forced air through Luer lock fitting.
- Package device for sterilization following AAMI recommended practices.
- Sterilize device using set points for desired sterilization method shown in following tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4°C (130°F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4°C (130°F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilization, the product should be taken out of service immediately.

AUTOClave CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	132.2°C (270°F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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
“Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers”, *AAMI Technical Information Report*, No. 12 - 1994.

“Sterilization, Part 1 - Good Hospital Practices”, *AAMI Standards and Recommended Practices*, Volume 1.1 - 1995.

 **Soehendra Stent Retriever Only**

RxOnly

STERILE	EO
---------	----

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

EC	REP
----	-----

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland