EN

D.A.S.H.® ERCP Catheter

CS

Katetr D.A.S.H.® pro ERCP

DA

D.A.S.H.® kateter til endoskopisk retrograd kolangiopankreatikografi

NL

D.A.S.H.® ERCP-katheter

FR

Cathéter de CPRE D.A.S.H.®

DE

D.A.S.H.® ERCP-Katheter

EL

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ENGLISH

INTENDED USE

This device is used for endoscopic cannulation of the ductal system. If preloaded, they will also aid in bridging difficult strictures during ERCP.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The Tracer Metro® is compatible with non-metal tip devices. Use of this wire guide with metal-tip ERCP devices may compromise integrity of external coating on wire guide.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling pancreatic duct may cause pancreatitis.

INSTRUCTIONS FOR USE

1. Secure wire guide or stylet in position with Tuohy-Borst adapter to avoid risk of movement.
2. Flush injection port with sterile water or saline. For best results, wire guide should be kept wet. **Note:** Tuohy-Borst adapter should be securely tightened around wire guide or stylet prior to injecting to prevent back-flow. Tuohy-Borst adapter may also be adjusted to allow wire guide movement while maintaining fluid seal.
3. Introduce catheter into endoscope accessory channel and advance in short increments until it is endoscopically visualized exiting scope.
4. Wire guide may be utilized to facilitate cannulation, if desired.
5. Fluoroscopically monitor wire guide advancement in ductal system.
6. Following cannulation, contrast may be injected through injection port to fluoroscopically confirm position of device. **Note:** When injecting, a 10 cc syringe (or smaller) allows easier injection.
7. If desired, stylet may be removed and replaced with a wire guide to maintain position and facilitate exchange of compatible wire-guided devices. To insert a compatible wire guide, loosen Tuohy-Borst adapter, remove stylet from catheter, and advance wire guide in short increments through Tuohy-Borst adapter. **Note:** Before each introduction of a wire guide, wet wire guide in a sterile water or saline bath. Flush contrast from catheter with sterile water or saline, observing necessary cautions if device is positioned in pancreatic duct. **Caution:** A contrast-filled catheter may make wire guide difficult to advance.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.