**COOK**

MEDICAL

CE  
0123

EN 2	<b>Biliary Dilation Catheters</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
CS 2	<b>Katetry pro dilataci žlučových cest</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
DA 3	<b>Katetre til galdedilatation</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
NL 4	<b>Galwegdilatatiekatheters</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
FR 5	<b>Cathéters de dilatation biliaires</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
DE 5	<b>Biliäre Dilatationskatheter</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
EL 6	<b>Καθετήρες διαστολής χοληφόρων</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
HU 7	<b>Epevezeték-tágító katéterek</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
IT 8	<b>Cateteri per dilatazione biliare</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
NO 9	<b>Galledilateringskatetere</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
PL 10	<b>Cewniki do poszerzania dróg żółciowych</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
PT 10	<b>Cateteres de dilatação biliar</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
ES 11	<b>Catéteres de dilatación biliar</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel
SV 12	<b>Biliära dilatationskatetrar</b> Soehendra®/Cotton®/Siegel-Cohen/van Andel

\* 1 8 9 2 4 / 0 6 1 9 \*

**INTENDED USE**

This device is used to dilate biliary strictures.

This device is indicated for adult use only.

**NOTES**

Store in a dry location, away from temperature extremes.

**CONTRAINDICATIONS**

Those specific to ERCP.

**WARNINGS**

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Do not use this device for any purpose other than stated intended use.

Use of this device restricted to a trained healthcare professional.

**POTENTIAL COMPLICATIONS**

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

**PRECAUTION**

Refer to package label for minimum channel size required for this device.

If unexpected resistance is met during the procedure, do not continue advancement of catheter without first determining the cause of the unexpected resistance and taking remedial action.

If bleeding occurs during the procedure, stop the procedure immediately and take remedial action.

**INSTRUCTIONS FOR USE**

1. Introduce device into accessory channel of endoscope over a pre-positioned wire guide. Advance in small increments until it is endoscopically visualized exiting endoscope and proceed with cannulation.
2. Under fluoroscopic monitoring, advance device to appropriate position to dilate stricture or papillary orifice and proceed with dilation.

**Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.**

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**