**EN**  
3

**Variable Injection Needle**

**CS**  
4

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**DA**  
6

**Variabel injektionsnål**

**NL**  
7

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**FR**  
9

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**DE**  
10

**Einstellbare Injektionskanüle**

**EL**  
12

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





## ENGLISH

### INTENDED USE

This device is used for endoscopic injection into gastrointestinal mucosa.

### NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

### CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired injection site.

Severe coagulopathy.

### POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

Those associated with mucosal injection include, but are not limited to: ulceration or necrosis of injected tissue, venous thrombosis, stricture formation.

Others to be considered if injection is performed in esophagus include, but are not limited to: chest pain, mediastinitis, pleural effusion.

## PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Needle must be fully retracted into sheath prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

## SYSTEM PREPARATION

1. Uncoil catheter and straighten completely.
2. Disengage Luer slip fitting from adjustment wheel prior to adjusting needle extension. (See fig. 1) **Caution:** Attempting to adjust needle while Luer slip fitting is engaged may result in damage to device. **Note:** A 1 cm black band at tip of catheter aids in needle visualization during injection (if applicable).
3. To ensure desired needle extension, advance Luer slip fitting and engage it firmly into adjustment wheel. (See fig. 2)
4. To adjust needle, retract it into outer sheath. (See fig. 3)
5. Hold device with Luer lock fitting facing operator and rotate adjustment wheel clockwise to lengthen needle extension, counterclockwise to shorten needle extension. (See fig. 4)
6. Verify adjustment by advancing and engaging Luer slip fitting. (See fig. 5) Needle should move in and out of outer sheath easily. Inspect needle for breaks. **Caution:** Advancing and retracting needle while catheter is coiled may result in stretching of catheter and damage to device.

## INSTRUCTIONS FOR USE

1. Securely attach a pre-filled syringe to Luer lock fitting on device, then flush air from needle. Check for patency and leaks.
2. With needle fully retracted into sheath, introduce device into accessory channel of endoscope.
3. Advance device in short increments until it is visualized exiting endoscope.
4. Endoscopically visualize desired injection site and maintain position of endoscope.
5. With sheath in endoscopic view, advance and secure needle by firmly engaging Luer slip fitting into adjustment wheel. **Note:** If desired, needle extension may be adjusted by withdrawing needle into sheath and fine tuning needle extension using adjustment wheel. **Caution:** During needle adjustment or advancement, ensure sheath of device has exited endoscope completely and is within endoscopic view.
6. Advance needle into mucosa, then using attached syringe perform injection.
7. When injection is complete, withdraw needle from injection site and retract it into sheath by pulling back on Luer slip fitting. Repeat procedure as necessary.
8. If applicable, to clear targeted area, attach a 20 cc syringe filled with sterile water to side arm fitting of device. (See fig. 6) Flush targeted area until desired effect is achieved. Repeat as necessary.
9. Upon completion of procedure, retract needle into sheath, then remove device from endoscope.

**Dispose of device(s) per institutional guidelines for biohazardous medical waste.**

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**