EN	Achalasia Balloon
CS	Balónek pro pacienty s achalazií
DA	Achalasia ballon
NL	Achalasieballon
FR	Ballonnet pour achalasia
DE	Ballon für Patienten mit Achalasia
EL	Μπαλόني διαστολής αλασίας οισοφάγου
HU	Achalasia kezelésénél asználatos ballon
IT	Palloncino per acalasia
NO	Akalasiballong
PL	Balon do achalazaji
PT	Balão para utilização em doentes com acalásia
ES	Balón para acalasia
SV	Akalasiballong

18931/0411

ENGLISH

INTENDED USE

This device is used to dilate strictures of the esophagus. [Specifically indicated for patients with Achalasia.]

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

This device is supplied non-sterile and may not be sterilized by any method.

CONTRAINDICATIONS

Those specific to upper GI endoscopy.

Those specific to dilation include, but are not limited to: uncooperative patient, asymptomatic rings or strictures, inability to advance the balloon through the strictured area, coagulopathy, known or suspect perforation, severe inflammation or scarring near the dilatation site.

POTENTIAL COMPLICATIONS

Those associated with upper GI endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Inflate balloon with air only.

Do not pre-inflate balloon.

WARNING

Do not exceed recommended balloon pressure as specified on the molded junction and package label.

EQUIPMENT REQUIRED

- Achalasia Balloon
- Three-way stopcock
- Savary-Gilliard® Wire Guide (not included)
- 60 cc Inflation Syringe (not included)
- Disposable Pressure Gauge (not included)

INSTRUCTIONS FOR USE

1. Introduce endoscope into patient and advance to proximal margin of stricture.
2. Introduce a Savary-Gilliard wire guide, floppy tip first, into accessory channel and advance in short increments until it is positioned well into stomach.
3. Remove endoscope, leaving wire guide in place.
4. Under fluoroscopic guidance, introduce achalasia balloon catheter over wire guide, and advance in short increments until dilator is in proper position.
5. Attach extension line of disposable pressure gauge to inflation hub on balloon catheter. Attach disposable pressure gauge to male Luer lock fitting on three-way stopcock, then attach a 60 ml inflation syringe to female Luer lock fitting opposite disposable pressure gauge.
6. Turn stopcock lever to "off" position (pointing toward disposable pressure gauge), then pull plunger of syringe back.
7. Fluoroscopically confirm that balloon remains in proper position, then turn stopcock lever 90° to "open" position (pointing toward side port). Under fluoroscopic visualization advance plunger of syringe to rapidly inflate balloon with air to recommended pressure (as specified on molded junction and package label). When balloon reaches appropriate pressure, turn stopcock lever to "off" position to maintain balloon pressure. **Note:** In order to reach recommended pressure, steps 6 and 7 may have to be repeated.
8. Proceed with dilation.
9. Turn stopcock lever toward inflation syringe and deflate balloon. Carefully remove deflated balloon and wire guide.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2011 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.