- EN**
3 **Wilson-Cook Achalasia Balloon**
Instructions for Use
- BG**
6 **Балон за ахалазия Wilson-Cook**
Инструкции за употреба
- CS**
9 **Balónek Wilson-Cook pro pacienty s achalazií**
Návod k použití
- DA**
12 **Wilson-Cook achalasia ballon**
Brugsanvisning
- DE**
15 **Wilson-Cook Ballon für Patienten mit Achalasie**
Gebrauchsanweisung
- EL**
19 **Μπαλόνι διαστολής αχалаσίας Wilson-Cook**
Οδηγίες χρήσης
- ES**
22 **Balón para acalasia Wilson-Cook**
Instrucciones de uso
- ET**
26 **Wilson-Cook ahhalasia balloon**
Kasutusjuhend
- FI**
28 **Wilson-Cook akalasiapallo**
Käyttöohjeet
- FR**
31 **Ballonnet pour achalasia Wilson-Cook**
Mode d'emploi
- HR**
35 **Balon za ahalaziju Wilson-Cook**
Upute za uporabu
- HU**
38 **Wilson-Cook achalasia kezelésénél
használatos ballon**
Használati utasítás
- IT**
41 **Palloncino per acalasia Wilson-Cook**
Istruzioni per l'uso

- LT**
44 **„Wilson-Cook“ achalazijos balionėlis**
Naudojimo instrukcija
- LV**
47 **Wilson-Cook ahalāzijas balons**
Lietošanas instrukcija
- NL**
50 **Wilson-Cook achalasiëballon**
Gebruiksaanwijzing
- NO**
54 **Wilson-Cook akalasiëballong**
Bruksanvisning
- PL**
57 **Balon do achalazji Wilson-Cook**
Instrukcja użycia
- PT**
60 **Balão para utilização em doentes com acalásia Wilson-Cook**
Instruções de utilização
- RO**
63 **Balon pentru acalazie Wilson-Cook**
Instrucțiuni de utilizare
- SK**
66 **Balónik na achaláziu Wilson-Cook**
Návod na použitie
- SL**
70 **Balon za ahalazijo Wilson-Cook**
Navodila za uporabo
- SV**
73 **Wilson-Cook-akalasiëballong**
Bruksanvisning
- TR**
76 **Wilson-Cook Akalazya Balonu**
Kullanma Talimatı

Wilson-Cook Achalasia Balloon

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Please read all instructions before using this device.

DEVICE DESCRIPTION

The recommended inflation pressure, balloon length and balloon diameter are provided in **Table 1**, and on the package label and the molded junction of the balloon.

Table 1

Product	Balloon Length	Recommended Inflation Pressure	Balloon Diameter
WCAB-30	8 cm	310 kPa / 3.0 atm	30 mm / 90 Fr
WCAB-35	8 cm	276 kPa / 2.7 atm	35 mm / 105 Fr

Patient Population

The device is indicated for adult use only.

Intended Users

Use of this device is restricted to a trained healthcare professional.

INTENDED USE

This device is used to dilate the lower esophageal sphincter in patients with Achalasia.

INDICATIONS FOR USE

This device is used to treat adult patients with Achalasia.

CLINICAL BENEFITS

This device relieves symptoms associated with Achalasia such as dysphagia.

CONTRAINDICATIONS

Those specific to upper GI endoscopy.

Those specific to dilation, including:

uncooperative patient • asymptomatic rings or strictures • inability to advance the balloon through the strictured area • coagulopathy • known or suspect perforation • severe inflammation or scarring near the dilatation site.

WARNINGS

- This device is supplied non-sterile, and therefore should not be used in a sterile field.

- Do not exceed the recommended balloon pressure as specified on the molded junction and package label.
- This device is designed for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.
- If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

PRECAUTIONS

- Inflate the balloon with air only. **Note:** Inflating with liquid can increase the risk of balloon damage, aspiration and/or perforation.
- Do not pre-inflate the balloon.
- Do not use this device for any purpose other than the stated intended use.

POTENTIAL COMPLICATIONS

Those associated with upper GI endoscopy, including: perforation

- hemorrhage • aspiration • fever • infection • allergic reaction to medication • hypotension • respiratory depression or arrest • cardiac arrhythmia or arrest.

HOW SUPPLIED

This device is supplied non-sterile and may not be sterilized by any method. Store in a dark, dry, cool place. Avoid extended exposure to light and extremes of temperature and humidity.

INSTRUCTIONS FOR USE

1. Introduce the endoscope into the patient and advance to the proximal margin of the stricture.
2. Introduce a wire guide, floppy tip first, into the accessory channel and advance in short increments until it is positioned well into the stomach.
3. Remove the endoscope, leaving the wire guide in place.
4. Attach a 60cc inflation syringe to the inflation device.
5. Attach the syringe extension line to the stopcock and then to the inflation hub of the balloon catheter.
6. With the plunger at the bottom of the syringe, turn the stopcock lever to the "open" position (pointing toward the side port) and then retract the plunger of the syringe back to aspirate all the residual air from the balloon.
7. While maintaining pressure, turn the stopcock lever to the "off" position (pointing toward the balloon).

8. Under fluoroscopic visualization, introduce the catheter over the wire guide, and advance in short increments until the dilator is in the proper position. **Caution:** Do not introduce the balloon without a wire guide as this may increase the risk for perforation.
9. Fluoroscopically confirm that the balloon remains in the proper position, then turn the stopcock lever 90° to the “open” position (pointing toward the side port).
10. Under fluoroscopic visualization, advance the plunger of the syringe to rapidly inflate the balloon with air up to the recommended pressure (as specified on the molded junction and the package label).
Balloon inflation may be monitored endoscopically by placing the endoscope next to the balloon catheter.
11. When the balloon reaches the appropriate pressure, turn the stopcock lever to the “off” position (pointing toward the balloon) while maintaining balloon pressure. In order to reach the recommended pressure, steps 10 & 11 may have to be repeated.
12. Proceed with dilation.
13. Turn the stopcock lever to the “open” position (pointing toward the side port) and deflate the balloon.
14. Carefully remove the deflated balloon, wire guide, and endoscope if applicable.

EQUIPMENT REQUIRED

- Achalasia Balloon
- Three-way stopcock
- Compatible Wire Guide (sold separately) - Compatible wire guide size is provided on the package label
- DS-60CC Inflation Syringe (sold separately)
- Inflation Device Handle (sold separately) - If use of an inflation handle is desired, select an inflation handle that is intended for use with a 60cc inflation syringe

DISPOSAL OF DEVICE

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken and limitations of use that the patient should be aware of.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device this should be reported to Cook Medical and the competent authority of the country where the device was used.

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William A. Cook Australia Pty. Ltd.
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.