EN	Fusion® Sphincterotome
CS	Sfinkterotom Fusion®
DA	Fusion® sphincterotom
NL	Fusion® sfincterotoom
FR	Sphinctérotome Fusion®
DE	Fusion®-Sfinkterotom
EL	Σφιγκτηροτόμος Fusion®
HU	Fusion® sphincterotóm
IT	Sfinterotomo Fusion®
NO	Fusion®-sfinkterotom
PL	Sfinkterotom Fusion®
PT	Esfincterótomo Fusion®
ES	Esfincterótomo Fusion®
SV	Fusion® sfinkterotom

ENGLISH

INTENDED USE

This device is used for cannulation of the ductal system and for sphincterotomy. If preloaded, also aids in bridging difficult strictures during ERCP.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with sphincterotomy. Contraindications to sphincterotomy include but are not limited to: coagulopathy and inability to properly position the sphincterotome cutting wire.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Maximum rated input voltage for this device is 1.5 kVp-p.

Any electrosurgical accessory constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

Before using this device, follow recommendations provided by electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of patient return electrode.

Ensure a proper path from patient return electrode to electrosurgical unit is maintained throughout procedure.

Switch electrosurgical unit to "off" position when it is not in use.

When applying current, ensure cutting wire is completely out of endoscope. Contact of cutting wire with endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope.

If a non-protected wire guide is used in sphincterotome, it must be removed prior to applying electrosurgical current.

Do not over flex or bow tip beyond 90°, as this may damage or cause cutting wire to break.

Elevator should remain open/down when advancing or retracting sphincterotome.

If preloaded, use of wire guide with metal tip ERCP devices may result in damage to external coating and/or tip of wire guide.

SYSTEM PREPARATION

- Upon removing device from package, uncoil and straighten sphincterotome. Carefully **remove precurved stylet from** cannulating tip. **Note:** Do not apply manual pressure to tip or cutting wire of sphincterotome in an attempt to influence orientation, as this may result in damage to device. **Note:** Do not exercise handle while device is coiled or precurved stylet is in place, as this may cause damage to sphincterotome and render it inoperable.
- With electrosurgical unit off, prepare equipment. Active cord fittings should fit snugly into both device handle and electrosurgical unit.
- Attach Wire Guide Locking Device to endoscope accessory port (if applicable).

INSTRUCTIONS FOR USE

- I. If using Intra Ductal Exchange (IDE) port and short wire.** (See fig. 1)
 - Remove Wire Stop from its retaining clip. Ensure that Wire Stop is disengaged from Proximal Wire Port (PWP) hub and that distal tip of Wire Stop is proximal to IDE port.
 - Insert distal tip of wire guide into IDE port and advance until it is flush with distal tip of sphincterotome.
 - Advance Wire Stop handle until it reaches PWP hub and connect Luer lock securely to PWP hub.
 - Advance tip of sphincterotome through cap of Wire Guide Locking Device and continue advancing until it is endoscopically visible.
 - Disengage Wire Stop by releasing handle from PWP hub and remove Wire Stop from catheter.
 - Following cannulation, contrast may be injected through injection port to fluoroscopically confirm position of device.
 - Following electrosurgical unit manufacturer's instructions, verify desired settings and proceed with sphincterotomy.
 - Upon completion of sphincterotomy, turn electrosurgical unit off.
 - Disconnect active cord from device handle and from electrosurgical unit. Wipe active cord with a damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping active cord tightly may damage device.

Note: Previously placed wire guide may be left in position, in order to facilitate introduction of other wire-guided devices. If wire guide is to remain in place while device is withdrawn, utilize following steps:

- Prior to removing device, utilize reference marks on catheter to ensure that IDE port is within ductal system.
- Fluoroscopically visualize radiopaque band at IDE port. Retract wire guide until radiopaque distal tip of wire guide passes band, a disengagement from wire guide lumen will occur.
- Advance disengaged wire guide to maintain ductal access.
- Lock wire guide into Wire Guide Locking Device and remove sphincterotome from endoscope accessory channel.

II. If using Proximal Wire Port (PWP) and a pre-positioned long wire guide.

Note: For best results, wire guide should be kept wet.

- Remove Wire Stop (if applicable).
- Advance sphincterotome over pre-positioned wire guide ensuring that wire exits catheter at PWP.
- Continue advancing device until it is endoscopically visible.

REFER TO STEPS 6-9 IN "SECTION I", THEN RESUME WITH STEP 4 BELOW:

- Remove device using standard long wire exchange technique.

III. If using Wire Control Port (WCP) and a short wire or long wire guide.

- Advance tip of sphincterotome through cap of Wire Guide Locking Device and continue advancing until it is endoscopically visible.

REFER TO STEPS 6-9 IN "SECTION I", THEN RESUME WITH STEP 2 BELOW:

- Introduce a wire guide into WCP and advance into duct of choice, separating wire guide from catheter. (See figs. 3 & 4)
- To withdraw sphincterotome from endoscope, ensure wire guide is locked in place, and then pull back on catheter, allowing wire guide lumen to separate from wire guide until metal band is visible at Wire Guide Locking Device and resistance is felt. (See fig. 5)
- Unlock wire guide from Wire Guide Locking Device, completely remove sphincterotome from wire guide, and re-lock wire guide.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.