EN	Memory™ Hard Wire Baskets
CS	Košíky z tvrdého drátu Memory™
DA	Memory™ hårde trådkurve
NL	Memory™ Hard Wire baskets
FR	Paniers Memory™ en fil rigide
DE	Memory™ -Hartdraht-Körbchen
EL	Καλάθια από σκληρό σύρμα Memory™
HU	Memory™ kemény drótkosarak
IT	Cestello Memory™ rigido
NO	Memory™ -kurver med stiv tråd
PL	Koszyki Memory™ z twardego drutu
PT	Cesto de arame rígido Memory™
ES	Cestas de alambre duro Memory™
SV	Memory™ korgar med styv tråd

ENGLISH

INTENDED USE

This device is used for the endoscopic removal of biliary stones and foreign bodies.

NOTES

This device is reusable if integrity remains intact.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in conjunction with removal of foreign body.

Perforation.

If device is to be used for biliary stone extraction, contraindications include an ampullary opening inadequate to allow for unimpeded passage of stone and basket.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy and ERCP include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction, respiratory depression or arrest, cardiac arrhythmia or arrest.

Others associated with basket extraction include, but are not limited to: impaction of object, aspiration of foreign body, localized inflammation, pressure necrosis.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

If device is to be used for removal of biliary stones, assessment of stone size and ampullary orifice must be made to determine necessity of sphincterotomy. If sphincterotomy is required, all appropriate cautions, warnings and contraindications must be observed.

This device should never be coiled in less than an 8-inch (20 cm) diameter.

WARNINGS

This device is **not** compatible with the Soehendra® Lithotripter or any other mechanical lithotripter.

When retrieving a foreign body from the upper GI tract, extreme caution must be exercised to avoid occlusion of airway.

Surgical intervention may be required if impaction occurs.

INSTRUCTIONS FOR USE

1. With basket fully retracted into sheath, insert device into endoscope accessory channel. Advance device through channel in short increments until basket sheath exits endoscope.
2. Endoscopically or fluoroscopically visualize extraction target and advance basket to appropriate position. Fluoroscopic visualization may be enhanced by injecting contrast via side arm fitting. Place sheath of basket beyond target before extending basket to entrap extraction target.
3. After confirming desired position of basket sheath relative to target, advance basket out of sheath by pushing forward on handle. Place basket around extraction target and slowly retract handle until entrapment is achieved. **Note:** It may be necessary to move basket back and forth along target in order to capture object.
4. Withdraw basket into channel while maintaining entrapment. Support foreign body or biliary stone against endoscope tip, and withdraw endoscope from patient. Maintain endoscopic visualization to ensure continued entrapment of object. While withdrawing device from endoscope, wipe secretions from outer sheath. **Caution:** If difficulty is encountered when removing basket from duct, moderate force may be applied by pulling on handle to manually fracture stone. If passage is still restricted, surgical intervention may be necessary.
5. Prepare or dispose of extracted specimen per institutional guidelines.
6. Upon completion of procedure, clean and prepare device for future use following cleaning instructions below. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, bends, or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.

CLEANING INSTRUCTIONS

1. Immediately after use, attach a syringe containing a minimum of 20 cc of warm water and enzymatic cleaner solution to Luer lock fitting on device, then irrigate.
2. Extend basket, then submerge solution-filled device in a solution of warm water and enzymatic cleaner. Soak for a minimum of 20 minutes.
3. Beginning at handle and working toward tip, wipe device with a textured cloth soaked in enzymatic solution.
4. Irrigate inner sheath by attaching a syringe containing a minimum of 20 cc of enzymatic solution to Luer lock fitting on handle of device. Repeat this step until all visible secretions are removed.
5. Gently brush basket with enzymatic solution and a soft nylon bristle brush to remove all foreign matter.
6. Remove device from enzymatic solution, then thoroughly rinse inner sheath using a syringe filled with clean water. Rinse handle, outer sheath and extended basket thoroughly with clean, running water.
7. Ultrasonically clean device, with basket extended, for a minimum of 15 minutes at 25°C (77°F) or per ultrasonic cleaner manufacturer's instructions.
8. Remove device from ultrasonic cleaner, then thoroughly rinse inner sheath with a syringe containing a minimum of 20 cc of clean water to remove all residue. Thoroughly rinse external portions of device and basket with clean, running water.
9. Gently wipe external portions of device dry with a lint-free cloth. Air dry lumen and basket by hanging device, with basket extended, for three hours, or by using medical grade (oil free) forced air through Luer lock fitting on handle.
10. Package device for sterilization following AAMI recommended practices.
11. Sterilize device using set points for desired sterilization method shown in tables.

ETO CYCLE PARAMETERS

PARAMETER	SET POINT
Prevacuum Pressure	2.0 psia (25.9" Hg Vac)
Relative Humidity	60% RH
Humidity Dwell	20 minutes
Chamber Temperature	54.4°C (130°F)
Exposure Pressure	To Be Determined ¹
Exposure Dwell	2 hours
ETO Concentration	600 mg/L
Post Vacuums	3
Post Vacuum Pressure	2.0 psia (25.9" Hg Vac)
Air Wash ²	5 minutes
Aeration Temperature	54.4°C (130°F)
Aeration Dwell	12 hours
Air Changes During Aeration	90 per hour

¹ To be calculated by the Institution to achieve an ETO concentration of 600 mg/L.

² Repetitive vacuums from ambient pressure to 13.7 psia over 5 minutes.

Note: ETO sterilizer chamber loading should not exceed 11 products per cubic foot of chamber space or exceed the usable chamber space. If any damage to a product is noted prior to sterilizing, it should be taken out of service immediately.

AUTOClave CYCLE PARAMETERS

PARAMETER	SET POINT
Cycle Type	Prevacuum
Prevacuum Pressure	2.5 psia (24.8" Hg Vac)
Chamber Exposure Temperature	132.2°C (270°F)
Chamber Exposure Pressure	41.8 psia (27.1 psig)
Exposure Dwell	5 minutes

Note: Autoclave chamber loading should not exceed 10 products per cubic foot of chamber space or exceed the usable chamber space.

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No compatible con litotricia

Ej litotripsikompatibel

“Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers”, *AAMI Technical Information Report*, No 12 - 1994.

“Sterilization, Part 1 - Good Hospital Practices”, *AAMI Standards and Recommended Practices*, Volume 1.1 - 1995.

Wilson-Cook Medical, Inc.
4900 Bethania Station Road
Winston-Salem, North Carolina 27105
USA

© 2013 Cook Medical

RxOnly

STERILE	EO
---------	----

EC	REP
----	-----

Cook Ireland Ltd.
O'Halloran Road
National Technology Park
Limerick
Ireland