EN

Lehman Sphincter of Oddi

Manometry Catheter

DA

Lehman manometrikateter til

sphincter Oddi

NL

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FR

Cathéter Lehman de manométrie

du sphincter d'Oddi

DE

Lehman Sphinkter

Oddi Manometrie-Katheter

EL

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ENGLISH

INTENDED USE

This device is used for motility studies in the common bile duct, pancreatic duct, Sphincter of ODDI and duodenum.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP.

POTENTIAL COMPLICATIONS

Those associated with ERCP and Sphincter of Oddi manometry include, but are not limited to: pancreatitis, cholangitis, sepsis, perforation, hemorrhage, aspiration, fever, infection, hypotension, allergic reaction to contrast or medication, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Set up and standardization may vary. Follow recording and perfusion system manufacturer's instructions.

Only Cook aspiration handles and wire guides are recommended for use with this device.

INSTRUCTIONS FOR USE

1. Attach red and yellow perfusion lumens to transducers and fill with sterile water.
2. Flush wire guide/aspiration lumen using a syringe filled with sterile water. Remove syringe, then place distal tip of manometry catheter over pre-positioned wire guide.
3. Advance manometry catheter until wire guide exits Luer lock fitting. Continue to advance manometry catheter over wire guide and into endoscope accessory channel until it is in desired position. **Note:** If opting to aspirate through wire guide/aspiration lumen, wire guide must be removed.
4. Perform Sphincter of Oddi manometry per institutional guidelines.

Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.