|                        |                                                         |
|------------------------|---------------------------------------------------------|
| <b>EN</b><br><b>2</b>  | <b>Dilation Syringe</b><br>Instructions for Use         |
| <b>DA</b><br><b>3</b>  | <b>Dilatationssprøjte</b><br>Brugsanvisning             |
| <b>DE</b><br><b>4</b>  | <b>Dilatationsspritze</b><br>Gebrauchsanweisung         |
| <b>EL</b><br><b>5</b>  | <b>Σύριγγα διαστολής</b><br>Οδηγίες χρήσης              |
| <b>ES</b><br><b>6</b>  | <b>Jeringa de dilatación</b><br>Instrucciones de uso    |
| <b>FI</b><br><b>7</b>  | <b>Laajennusruisku</b><br>Käyttöohjeet                  |
| <b>FR</b><br><b>8</b>  | <b>Seringue de dilatation</b><br>Mode d'emploi          |
| <b>HU</b><br><b>9</b>  | <b>Dilatációs fecskendő</b><br>Használati utasítás      |
| <b>IT</b><br><b>11</b> | <b>Siringa di dilatazione</b><br>Istruzioni per l'uso   |
| <b>NL</b><br><b>12</b> | <b>Dilatatiespuit</b><br>Gebruiksaanwijzing             |
| <b>NO</b><br><b>13</b> | <b>Dilatasjonssprøyte</b><br>Bruksanvisning             |
| <b>PL</b><br><b>14</b> | <b>Strzykawka do poszerzania</b><br>Instrukcja użycia   |
| <b>PT</b><br><b>15</b> | <b>Seringa de dilatação</b><br>Instruções de utilização |
| <b>SV</b><br><b>16</b> | <b>Dilatationsspruta</b><br>Bruksanvisning              |

## Dilation Syringe

**Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.**

**CAUTION: US federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

## DEVICE DESCRIPTION

### Performance Characteristics

Dilation syringes inflate, deflate, and monitor the pressure of Cook GI dilation balloons. The pressure gauge on the syringe features a fluorescent face for visibility in darkened endoscopy suites. Dilation syringes come packaged with a stopcock for added pressure control. Refer to the table below for the pressure gauge accuracy.

| Accuracy | Pressure Range |                 |
|----------|----------------|-----------------|
| 3%       | Vac - 50 psi   | Vac - 3.25 atm  |
| 2%       | 50 - 150 psi   | 3.25 - 9.75 atm |
| 3%       | 150 - 200 psi  | 9.75 - 13.0 atm |

**Note:** Vacuum is abbreviated as 'Vac'.

### Device Compatibility

This device is compatible with 60 cc syringe inflation handles and Cook GI dilation balloons. Refer to GI dilation balloon labeling for the maximum recommended inflation pressure.

### Patient Population

This device is indicated for adult use only.

### Intended User

The use of this device is restricted to a trained healthcare professional.

### Contact with Body Tissue

This device is not patient contacting.

### Operating Principle

Once a GI dilation balloon is in place, the dilation syringe is used to inflate the balloon to the recommended pressure and maintain that pressure for the desired amount of time. Once the dilation procedure is finished, the dilation syringe is used to deflate the balloon until a vacuum is created and the balloon can be removed from the scope.

## INTENDED USE

This device is used to inflate, deflate, and monitor pressure of the Cook line of balloon dilators during GI balloon dilation.

## INDICATIONS FOR USE

The dilation inflation syringe is used with the Cook GI dilation balloons for endoscopic treatment of strictures throughout the primary GI tract and achalasia.

## CLINICAL BENEFITS

This device supports the Cook line of balloon dilators during GI balloon dilation.

## CONTRAINDICATIONS

Contraindications include those specific to dilation.

## WARNINGS

Do not use this device for any purpose other than stated intended use.

### Single Use

This device is supplied non-sterile and designed for single use only. Attempts to reprocess and/or reuse may lead to contamination with biological or chemical agents and/or mechanical integrity failure of device.

## POTENTIAL COMPLICATIONS

This device is used with Cook GI dilation balloons. Read the appropriate Instructions for Use for complications associated with the procedure and device being used.

## HOW SUPPLIED

This device is supplied non-sterile.

## STORAGE

Store in a dry location, away from temperature extremes.

## INSPECTION OF DEVICE

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends, and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

## INSTRUCTIONS FOR USE

1. Assemble and use the Dilation Syringe according to inflation handle manufacturer's instructions.

## DISPOSAL OF DEVICE

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

## PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of.

## SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to this device, this should be reported to Cook Medical, the competent authority, and/or regulatory authority of the country where the device was used.

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## RxOnly

This device is intended  
for single use only.

[cookmedical.com/symbol-glossary](https://cookmedical.com/symbol-glossary)

**Wilson-Cook Medical, Inc.**  
4900 Bethania Station Road  
Winston-Salem, North Carolina  
27105  
USA

**Cook Medical Europe LTD**  
O'Halloran Road  
National Technology Park  
Limerick  
Ireland

**William A. Cook Australia Pty. Ltd.**  
95 Brandl Street,  
Eight Mile Plains  
QLD 4113  
Australia

2025-01  
18945\_0125  
© COOK 2025  
[cookmedical.com](https://cookmedical.com)