COOK®


MEDICAL

CE
0123


EN	Zimmon Needle Knife Papillotomes
3	
CS	Jehlové papilotomy Zimmon
4	
DA	Zimmon papillotomkanylekniv
6	
NL	Zimmon-naaldmespapillotomen
8	
FR	Papillotomes-couteau-aiguille Zimmon
11	
DE	Zimmon Nadelmesser-Papillotome
13	
EL	Θηλοτόμοι – Βελονοτόμοι Zimmon
15	
HU	Zimmon tűkések papillotómok
18	
IT	Papillotomi ad ago Zimmon
20	
NO	Zimmon nålekniv-papillotomer
22	
PL	Papilotomy igłowe Zimmon
24	
PT	Papilótomos de agulha cortante Zimmon
26	
ES	Papilótomos de aguja/precorte Zimmon
28	
SV	Zimmon nålknivspapillotomer
30	

* 1 8 9 5 3 / 0 7 1 9 *

1	Injection Port Injekční port Injektionsport Injectiepoort Orifice d'injection Injektionsport Θύρα έγχυσης Befecskendezési nyílás Raccordo di iniezione Injektionsport Port do wstrzykiwania Orificio de injeção Orificio de inyección Injektionsport	2	Active Cord Connection Připojení aktivního vodiče Tilslutning til strømførende ledning Aansluiting actief snoer Connexion du cordon actif HF-Kabelanschluss Σύνδεση ενεργού καλωδίου Aktív vezeték csatlakozója Connessione per cavo diatermico Tilkobling til strømførende ledning Złącze przewodu zasilania Ligação do cabo ativo Conexión del cable de conexión Anslutning strömförande kabel	3*	Catheter Katetr Kateter Katheter Cathéter Katheter Καθετήρας Katéter Catetere Kateter Cewnik Cateter Catéter Kateter	4*	Cutting Wire Diskéční drát Skæretråd Snijdraad Fil de coupe Schneidedraht Χορή τομή Vágódrót Filo di taglio Skjæretråd Cięciwa tnąca Fio de corte Alambre de corte Skärtråd	5	Distal Markings Distální značky Distale markeringer Distale markeringen Repères distaux Distale Markierungen Περιφερικές σημάνσεις Disztális jelzések Marker distali Distale merker Oznaczenia dystalne Marcas distales Marcas distales Distala markeringar
---	---	---	---	----	---	----	--	---	--

INTENDED USE

These devices are used for accessing the common bile duct when standard methods of cannulation have been exhausted.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with papillotomy.

Contraindications to papillotomy include, but are not limited to: coagulopathy.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP.

Potential adverse events associated with papillotomy include, but are not limited to: pancreatitis, bleeding, perforation, and infection.

USER QUALIFICATIONS

This device should only be used by or under the supervision of physicians who have been trained in ERCP procedures.

WARNINGS

Do not use this device for any purpose other than stated intended use.

This device is designed for single use only. Attempts to reprocess, sterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Any electrosurgical device constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

The electrosurgical effect is greatly influenced by the size and configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.

Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper selection, placement, and utilization of patient return electrode. Ensure that a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

For patients with electrically conductive implants, monopolar diathermy or electrosurgical cautery can result in localized thermal injury at the implant site.

Do not use this device with an output higher than the rated voltage: 2 kVp-p (1 kVp). This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, device and/or active cord.

Do not use an active cord with a maximum rated voltage less than 1 kVp. This could cause thermal injury to the patient, operator, or assistant and could also damage the endoscope, device and/or active cord.

If cutting is not observed at normal operating settings this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.

The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.

Skin-to-skin contact should be avoided (for example between the patient's arms and body) by using dry gauze.

Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.

Avoid contact between the active cords and the patient's body, or any other electrodes.

This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electric shock in compliance with IEC 60601-1). Use of generator other than type BF or CF could cause electric shock to the patient, operator, or assistant.

For recommendations on current settings for papillotomy, refer to the monopolar generator's recommended power settings for papillotomy.

This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere, or in the presence of explosive gases. Doing so may result in an explosion.

In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electrosurgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the device in these patients.

Using accessories others than those recommended by the electrosurgical generator manufacturer may result in decreased immunity or increased emissions.

This device must be used with an endoscope that is intended for use with high frequency (HF) electrosurgical devices. If the device is used with an endoscope not compatible with high frequency (HF) electrosurgery, a burn to the patient or the operator may occur.

To avoid neuromuscular stimulation and minimize electric arc, activate the cutting wire only when it is in contact with the tissue to be cauterized.

Additional warnings may be listed in procedural steps below.

PRECAUTIONS

Switch the electrosurgical unit to the "off" position when it is not in use.

The elevator should remain open/down when advancing or retracting the device.

TECHNICAL DESCRIPTION

Figure 1 and Figure 2 summarize the device's characteristics.

Refer to the package label for device dimensions, minimum endoscope compatibility, and wire guide compatibility.

This device must only be used with an active cord compatible with a 3 mm diameter connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (supplied non-sterile).

The maximum rated input voltage for this device is 2 kVp-p (1 kVp).

Device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6° F).

Store in a dry location. Recommended storage temperature 22°C (71.6° F). Device may be exposed to temperatures up to 50°C (122° F) during transportation.

Components highlighted with an asterisk (*) denote applied parts (part of the device that comes into contact with the patient). Prior to starting the procedure in which the device will be used, become familiarized with distal mark location with respect to the cutting wire.

WARNING

-No modification of this equipment is allowed.

SYSTEM PREPARATION

WARNINGS

- For recommendations on current settings for papillotomy, refer to the monopolar generator's recommended power settings for papillotomy.

-Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode.

Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

1. Upon removing the device from the package uncoil the device.
2. With the electrosurgical unit off, prepare the equipment.

INSTRUCTIONS FOR USE

1. Advance the tip of the device into the endoscope accessory channel and continue to advance in short increments until the device is endoscopically visible.

WARNING

-Using the endoscope monitor, verify that the cutting wire has exited the endoscope by visualizing the distal end of the catheter is past the elevator. Failure to confirm the cutting wire has fully exited the endoscope may result in contact between the cutting wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

2. Once the tip is visualized exiting the endoscope, position the device at the desired site.

PAPILLOTOMY

WARNINGS

- Do not use a wire guide or device that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

- The device should never be connected to the active cord prior to inserting the device through the endoscope to avoid injury to patient and or equipment resulting from improper electrical circuit grounding.

- The device should not be energized prior to performing papillotomy. Energizing the cutting wire prior to use will cause premature cutting wire fatigue and will compromise the cutting wire's integrity.

- To avoid unintended disconnection during procedure and damage to device, the device should only be used with compatible active cords (Reference Technical Description for compatibility).

1. Attach the active cord to the connector on the device handle. The active cord fittings should fit snugly into both the device handle and the electrosurgical unit.
2. Following the electrosurgical unit manufacturer's instructions, verify the desired settings and proceed with the papillotomy.

PRECAUTION

-It is essential to move the needle knife while applying current.

WARNINGS

- The user must maintain contact with tissue when applying electrocautery current. Failure to do so may result in increased current resulting in a broken cutting wire, damage to the endoscope, and/or patient injury.

- When applying current, ensure the cutting wire is completely out of endoscope. Contact of the cutting wire with the endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope.

3. Insert the needle knife into the papillary orifice and, using a continuous motion, activate the device.
4. Upon completion of papillotomy, turn the electrosurgical unit off.
5. Disconnect the active cord from the device handle and from the electrosurgical unit. Clean and store the active cord per the manufacturer's instructions.

CANNULATION

1. After papillotomy, the device may be used to cannulate the ductal system.
2. Completely remove the cutting wire from the catheter (see fig. 2).
3. Following cannulation, contrast may be injected through the injection port or the Luer lock fitting on the handle to fluoroscopically confirm the position of the device.
4. Once access to the ductal system is obtained and contrast is injected (if applicable), a wire guide may be advanced to facilitate the introduction of compatible devices (see fig. 2). Prior to advancing the wire guide, inject water through the injection port of the Luer lock fitting on the handle to flush out remaining contrast in device (if applicable).

DEVICE REMOVAL FROM ENDOSCOPE

1. If a pre-positioned wire guide is in place, remove the device using standard long wire exchange technique.

DEVICE DISPOSAL

1. Upon completion of the procedure, dispose of the device(s) per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
Wilson-Cook Medical, Inc.
Cook Medical Europe Ltd.