EN

Zimmon Needle Knife Papillotomes

CS

Jehlové papilotomy Zimmon

DA

Zimmon kanylekniv-papillotomer

NL

Zimmon naaldmespapillotomen

FR

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ENGLISH

INTENDED USE

These devices are used for accessing the common bile duct when standard methods of cannulation have been exhausted.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with papillotomy.

Contraindications to papillotomy include, but are not limited to: coagulopathy.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The maximum rated input voltage for this device is 2 kVp-p.

Any electrosurgical accessory constitutes a potential electrical hazard to the patient and the operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

Before using this device, follow the recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through the proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

Switch the electrosurgical unit to the "off" position when it is not in use.

When applying current, ensure the needle knife is completely out of the endoscope. Contact of the needle knife with the endoscope may cause grounding, which can result in patient injury, operator injury, a broken needle knife, and/or damage to the endoscope.

If a non-protected wire guide is used in the papillotomy, it must be removed prior to applying electrosurgical current.

The elevator should remain open/down when advancing or retracting the papillotomy.

SYSTEM PREPARATION

1. Upon removing the device from the package, uncoil the papillotomy.
2. With the electrosurgical unit off, prepare the equipment. The active cord fittings should fit snugly into both the device handle and the electrosurgical unit.

INSTRUCTIONS FOR USE

1. Advance the tip of the papillotomy into the endoscope accessory channel and continue to advance in short increments until the device is endoscopically visible.
2. Once the tip is visualized exiting the endoscope, position the papillotomy at the desired site.
3. Following the electrosurgical unit manufacturer's instructions, verify the desired settings.
4. Insert the needle knife into the papillary orifice and, using a continuous motion, activate the needle knife. **Caution:** It is essential to move the needle knife while applying current.
5. Upon completion of papillotomy, turn the electrosurgical unit off by disconnecting the active cord, and completely remove the cutting wire from the catheter.
6. The papillotomy may then be used to cannulate the ductal system.

Note: Once access to the ductal system is obtained, a wire guide may be advanced to facilitate the introduction of compatible devices. (See fig. 1)

7. Following cannulation, contrast may be injected through the injection port (if applicable) or the Luer lock fitting on the handle to fluoroscopically confirm the position of the device.
8. Wipe the active cord with a damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping the active cord tightly may damage the device.

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.