18954/1013

ENGLISH

INTENDED USE

This device is used for accessing the common bile duct when standard methods of cannulation have been exhausted and for papillotomy. This device is supplied sterile and is intended for single use only.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to ERCP and any procedures to be performed in conjunction with papillotomy.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with ERCP.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The maximum rated input voltage for this device is 1.67 kVp-p.

Any electrosurgical accessory constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through the proper placement and utilization of patient return electrode. Ensure a proper path from patient return electrode to electrosurgical unit is maintained throughout procedure.

Switch electrosurgical unit to "off" position when it is not in use.

When applying current, ensure needle knife is completely out of endoscope. Contact of needle knife with endoscope may cause grounding, which can result in patient injury, operator injury, a broken needle knife, and/or damage to endoscope.

If a non-protected wire guide is used in papillotomy, it must be removed prior to applying electrosurgical current.

Elevator should remain open/down when advancing or retracting papillotomy.

SYSTEM PREPARATION

- Upon removing device from package, uncoil papillotomy and carefully remove tip protector from distal tip of device. **Note:** Do not extend or retract needle knife while the device is coiled, as this may cause damage to papillotomy and render it inoperable.
- With electrosurgical unit off, prepare the equipment. Active cord fittings should fit snugly into both device handle and electrosurgical unit.

INSTRUCTIONS FOR USE

- With needle knife fully retracted into sheath, advance tip of papillotomy into endoscope accessory channel and continue to advance in short increments until device is endoscopically visible.
- Once tip is visualized exiting endoscope, position catheter and extend needle knife to desired length. **Note:** Extension may be adjusted by loosening thumbscrew and advancing or retracting adjustable ring on handle (if applicable). (See fig. 1) **Caution:** It is essential to tighten thumbscrew (if applicable) after adjusting needle knife length and before proceeding with papillotomy.
- Following electrosurgical unit manufacturer's instructions, verify desired settings.
- Insert needle knife into papillary orifice and, using a continuous motion, activate needle knife. **Caution:** It is essential to move needle knife while applying current.
- Upon completion of papillotomy, turn electrosurgical unit off, disconnect active cord, and completely remove cutting wire from catheter or retract needle knife into sheath (if applicable).
- Papillotomy may then be used to cannulate ductal system.

Note: Wire guide may be placed in order to facilitate introduction of compatible devices (if applicable). (See fig. 1 or 2)

- Following cannulation, contrast may be injected through injection port (if applicable) or through Luer lock fitting on the handle (if applicable) to fluoroscopically confirm position of device.
- Wipe active cord with a damp cloth to remove all foreign matter. Store in a loose coil.

Note: Wrapping active cord tightly may damage device.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.