**COOK®**

**MEDICAL**


**CE**  
**0123**

**PT, PTG**

<b>EN</b> <b>3</b>	<b>Billroth II Sphincterotomes</b>
<b>CS</b> <b>5</b>	<b>Sfinkterotomy Billroth II</b>
<b>DA</b> <b>7</b>	<b>Billroth II sphincterotomer</b>
<b>NL</b> <b>10</b>	<b>Billroth II-sfincterotomen</b>
<b>FR</b> <b>12</b>	<b>Sphinctérotomes Billroth II</b>
<b>DE</b> <b>15</b>	<b>Billroth-II-Sfinkterotome</b>
<b>EL</b> <b>18</b>	<b>Σφιγκτηροτόμοι Billroth II</b>
<b>HU</b> <b>21</b>	<b>Billroth II sphincterotómok</b>
<b>IT</b> <b>23</b>	<b>Sfinterotomi per Billroth II</b>
<b>NO</b> <b>26</b>	<b>Billroth II-sfinkterotomer</b>
<b>PL</b> <b>28</b>	<b>Sfinkterotomy Billroth II</b>
<b>PT</b> <b>31</b>	<b>Esfincterótomos Billroth II</b>
<b>ES</b> <b>33</b>	<b>Esfinterótomos Billroth II</b>
<b>SV</b> <b>36</b>	<b>Billroth II sfinkterotomer</b>

\* 1 8 9 5 5 / 0 7 1 9 \*

- 1 Injection Port  
 Injekčni port  
 Injektionsport  
 Injectiepoort  
 Orifice d'injection  
 Injektionsport  
 Θύρα έγχυσης  
 Befecskendezési nyílás  
 Raccordo di iniezione  
 Injeksjonsport  
 Port do wstrzykiwań  
 Ofício de injeção  
 Orificio de inyección  
 Injektionsport

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## ENGLISH

### INTENDED USE

This device is used for the cannulation of the ductal system and for sphincterotomy.

### CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with sphincterotomy.

Contraindications to sphincterotomy include, but are not limited to: coagulopathy and inability to properly position the sphincterotome cutting wire.

### POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

### USER QUALIFICATIONS

The sphincterotome should only be used by or under the supervision of physicians who have been trained in ERCP procedures.

### WARNINGS

Do not use this device for any purpose other than stated intended use.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify cook for return authorization.

Any electrosurgical accessory constitutes a potential electrical hazard to the patient and the operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

The electrosurgical effect is greatly influenced by the size and the configuration of the active electrode. Therefore, it is impossible to determine the exact effect achieved in a given control setting. If the proper generator setting is not known, one should set the unit at a power setting lower than the recommended range and cautiously increase the power until the desired effect is achieved. Excess power may lead to patient injury, or may damage the integrity of the cutting wire.

Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through the proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

For patients with electrically conductive implants, monopolar diathermy or electrosurgical cautery can result in localized thermal injury at the implant site.

Do not use this device with an output higher than the rated voltage 1.5 kVp-p (750 Vp). This could cause thermal injury to

the patient, operator, or assistant and could also damage the endoscope, sphincterotome and/or active cord.

Do not use this device with an active cord which has a maximum voltage rating less than 1.5 kVp-p (750 Vp). This could cause thermal injury to the patient, operator, or an assistant.

If cutting is not observed at the normal operating settings, this may indicate faulty application of the neutral electrode or poor contact in its connections. Confirm the neutral electrode and its connectors are attached before selecting a higher output power.

The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth.

Skin-to-skin contact should be avoided (for example between the patient's arms and body) by using a dry gauge.

Monitoring electrodes should be placed as far as possible from the surgical area. Needle monitoring electrodes are not recommended.

Avoid contact between the active cords and the patient's body, or any other electrodes.

This device must only be used in conjunction with a type BF or CF generator (generators that provide protection from electric shock in compliance with IEC60601-1). Use of generator other than type BF or CF could cause electric shock to the patient, operator, or assistant.

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

This device is not intended to be used in the presence of flammable liquid, in an oxygen enriched atmosphere or in the presence of explosive gases. Doing so may result in explosion.

In patients with pacemakers or implantable cardiac defibrillators, monopolar diathermy or electrosurgical cautery can result in electrical reset of the cardiac device, inappropriate sensing and/or therapy, tissue damage around the implanted electrodes, or permanent damage to the pulse generator. A cardiologist should be consulted prior to using the sphincterotome in these patients.

Using accessories other than those recommended by the electrosurgical generator manufacturer may result in decreased immunity or increased emissions.

This device must be used with an endoscope that is intended for use with high frequency (HF) electrosurgical devices. If the device is used with an endoscope not compatible with HF electrosurgery, a burn to the patient or the operator may occur.

To avoid neuromuscular stimulation and minimize electric arc, activate the cutting wire only when it is in contact with the tissue to be cauterized.

Additional warnings may be listed in procedural steps below.

### PRECAUTIONS

Switch the electrosurgical unit to the "off" position when it is not in use.

Do not over flex or bow tip beyond 90°, as this may damage or cause the cutting wire to break.

The elevator should remain open/down when advancing or retracting the sphincterotome.

### TECHNICAL DESCRIPTION

See fig 1 for PTG-20-6-BII-NG and fig 2 for PT-5-5-BII-SOEHENDRA sphincterotomes. Components marked with an asterisk (\*) denote applied parts (the part of the device that comes into contact with the patient, as defined by IEC60601-1), see figs. 1 and 2. Prior to starting the procedure in which the device will be used, become familiarized with distal mark location with respect to the cutting wire.

Refer to the package label for device dimensions, minimum endoscope compatibility and wire guide compatibility.

This device is recommended for use with Cook electrosurgical compatible wire guides.

This device must only be used with an active cord compatible with a 3 mm diameter connector. This device has only been verified to be compatible with the following Cook active cords: ACU-1 and ACU-1-VL (Supplied Non-Sterile).

Maximum rated input voltage for this device is 1.5 kVp-p (750 Vp).

This device should be used in a professional healthcare facility in an environment with a temperature less than 37°C (98.6° F).

Store in a dry location. Recommended storage temperature 22°C (71.6° F). Device may be exposed to temperatures up to 50°C (122° F) during transport.

### WARNING

No modification of this equipment is allowed.

### SYSTEM PREPARATION

#### PRECAUTIONS

Do not apply manual pressure to the tip or the cutting wire of the sphincterotome in an attempt to influence orientation, as this may result in damage to the device.

Do not exercise the handle while the device is coiled, as this may cause damage to the sphincterotome and render it inoperable.

1. Upon removing the device from package, uncoil and straighten the sphincterotome.

### WARNINGS

For recommendations on current settings for sphincterotomy, refer to the monopolar generator's recommended power settings for sphincterotomy.

Before using this device, follow the recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to the electrosurgical unit is maintained throughout the procedure.

2. With the electrosurgical unit off, prepare the equipment. Active cord fittings should fit snugly into both the device handle and the electrosurgical unit.

## INSTRUCTIONS FOR USE

### INSERTION

#### WARNING

The sphincterotome should never be connected to the active cord prior to inserting the sphincterotome through the endoscope. Doing so may injure the patient and/or damage the equipment due to improper electrical circuit grounding.

1. Advance the tip of sphincterotome into the endoscope accessory channel and continue to advance in short increments until the device is endoscopically visible. **Note:** If using a double lumen device, it may be placed over a pre-positioned wire guide. Prior to advancing the device over a pre-positioned wire guide, flush the wire guide lumen.

### CANNULATION

#### WARNING

For the PTG-20-6-BII-NG, using the endoscope monitor, verify that the cutting wire has fully exited the endoscope. For the PT-5-5-BII-SOEHENDRA, using the endoscope monitor, verify that the cutting wire has exited the endoscope by visualizing the distal end of the catheter is past the elevator. Failure to confirm that the cutting wire has fully exited the endoscope may result in contact between the cutting wire and the endoscope while current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

2. Using endoscopic visualization, confirm the cutting wire has exited the endoscope.
3. Cannulate using the sphincterotome.
4. Following cannulation, confirm the position fluoroscopically via contrast injection through the injection port or visualization of the radiopaque marker.

### SPHINCTEROTOMY

5. Attach the active cord to the connector on the sphincterotome handle. Active cord fittings should fit snugly into both the device handle and the electrosurgical unit.

#### WARNING

To avoid unintended disconnection during procedure, and damage to the device, device should only be used with compatible active cords (Reference Technical Description for compatibility).

6. If using a wire guide that is not compatible with electrosurgical procedures, remove the wire guides from the sphincterotome before activating the electrosurgical current. If the wire is not removed, there is a risk of thermal injury to the patient or the operator.

#### WARNINGS

Do not use a wire guide (or sphincterotome) that has been cut, burned, or damaged. Damaged insulation may cause unsafe currents in either the patient or the operator. Leakage current to the patient or user may increase at the site of the damaged insulation.

The sphincterotome should not be energized prior to performing sphincterotomy. Energizing the cutting wire prior

to use will cause premature cutting wire fatigue and will compromise the cutting wire's integrity.

Verify that the cutting wire has exited the endoscope by visualizing it on the endoscope monitor. Failure to do so may result in contact between the cutting wire and the endoscope while the electrical current is applied. This may cause grounding, which can result in patient injury, operator injury, a broken cutting wire and/or damage to the endoscope.

User must maintain constant contact with the tissue when applying electrosurgical current. Failure to do so may result in increased current resulting in broken cutting wire, damage to the endoscope and/or patient injury.

When applying current, ensure the cutting wire is completely out of the endoscope. Contact of cutting wire with the endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to the endoscope.

If a wire guide that is not compatible with electrosurgical procedures is used in the sphincterotome, it must be removed prior to applying the electrosurgical current.

Cook wire guides packaged with this device may be left in place.

7. Following electrosurgical unit manufacturer's instructions, verify desired settings and proceed with sphincterotomy.
8. Upon completion of sphincterotomy, turn the electrosurgical unit off.
9. Disconnect the active cord from the device handle and from the electrosurgical unit. Clean and store the active cord per the active cord manufacturer's instructions.

#### DEVICE REMOVAL FROM ENDOSCOPE

1. A previously placed wire guide may be left in position, in order to facilitate introduction of compatible devices (if applicable).

#### DEVICE DISPOSAL

1. Upon completion of the procedure, dispose of the device(s) per institutional guidelines for biohazardous medical waste.

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© 2019 Cook Medical
**Wilson-Cook Medical, Inc.**
**Cook Medical Europe Ltd.**