EN	PreCurved Double Lumen Sphincterotome
CS	Preformovaný dvojlumenný sfinkterotom
DA	Buede sphincterotom med dobbelt lumen
NL	Voorgevormde sfincterotoom met dubbel lumen
FR	Sphinctérotome à double lumière précourbée
DE	Vorgebogenes doppelumiges Sphinkterotom
EL	Προκεκαμμένος σφιγκτηροτόμος διπλού αυλού
HU	Előhajlított kétlumenű sphincterotóm
IT	Sfinterotomo precurvato a lume doppio
NO	Forhåndsbøyd sfinkterotom med dobbeltlumen
PL	Zakrzywiony sfinkterotom dwukanałowy
PT	Esfinterótomo pré-curvado de duplo lúmen
ES	Esfinterótomo de doble luz precurvado
SV	Förböjd sfinkterotom med dubbellumen

ENGLISH

INTENDED USE

This device is used for cannulation of the ductal system and for sphincterotomy.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP and any procedures to be performed in conjunction with sphincterotomy.

Contraindications to sphincterotomy include, but are not limited to: coagulopathy and inability to properly position sphincterotome cutting wire.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Maximum rated input voltage for this device is 1.5 kVp-p.

Any electrosurgical accessory constitutes a potential electrical hazard to patient and operator. Possible adverse effects include, but are not limited to: fulguration, burns, nerve and/or muscle stimulation and cardiac arrhythmia.

Before using this device, follow recommendations provided by electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of patient return electrode.

Ensure a proper path from patient return electrode to electrosurgical unit is maintained throughout procedure.

Switch electrosurgical unit to "off" position when it is not in use.

When applying current, ensure cutting wire is completely out of endoscope. Contact of cutting wire with endoscope may cause grounding, which can result in patient injury, operator injury, a broken cutting wire, and/or damage to endoscope.

If a non-protected wire guide is used in sphincterotome, it must be removed prior to applying electrosurgical current.

Do not over flex or bow tip beyond 90°, as this may damage or cause cutting wire to break.

Elevator should remain open/down when advancing or retracting sphincterotome.

SYSTEM PREPARATION

- Upon removing device from package, uncoil and straighten sphincterotome. Carefully remove precurved stylet from cannulating tip. **Note:** Do not apply manual pressure to tip or cutting wire of sphincterotome to influence orientation, as this may result in damage to device. **Note:** Do not exercise handle while device is coiled or precurved stylet is in place, as this may cause damage to sphincterotome and render it inoperable.

- With electrosurgical unit off, prepare equipment. Active cord fittings should fit snugly into both device handle and electrosurgical unit.

INSTRUCTIONS FOR USE

- Advance tip of sphincterotome into endoscope accessory channel and continue to advance in short increments until device is endoscopically visible. **Note:** This sphincterotome may be placed over a pre-positioned wire guide. Prior to advancing device over a pre-positioned wire guide, flush wire guide lumen.

- Following cannulation, contrast may be injected to fluoroscopically confirm position of device.

Note: If applicable, move slide to "off" position prior to injecting over pre-positioned wire guide. To manipulate wire guide, move slide to "on" position. (See fig. 1)

- Following electrosurgical unit manufacturer's instructions, verify desired settings and proceed with sphincterotomy.

- Upon completion of sphincterotomy, turn electrosurgical unit off.

Note: Previously placed wire guide may be left in position, in order to facilitate introduction of compatible devices (if applicable).

- Disconnect active cord from device handle and from electrosurgical unit. Wipe active cord with a damp cloth to remove all foreign matter. Store in a loose coil. **Note:** Wrapping active cord tightly may damage device.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

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© 2013 Cook Medical
Wilson-Cook Medical, Inc.
Cook Ireland Ltd.