EN Haber RAMP™ Catheter

CS Katetr Haber RAMP™

DA Haber RAMP™ kateter

NL Haber RAMP™ katheter

FR Cathéter Haber RAMP™

DE Haber RAMP™ Katheter

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ENGLISH

INTENDED USE

This device is used for endoscopic cannulation of the ductal system.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL COMPLICATIONS

Potential complications associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to the package label for the minimum channel size required for this device.

Cook Endoscopy wire guides are recommended for use with this device.

Injection of contrast during ERCP should be monitored fluoroscopically. Overfilling of the pancreatic duct may cause pancreatitis.

INSTRUCTIONS FOR USE

1. Introduce the catheter into the endoscope and advance in short increments until the device is endoscopically visualized exiting the scope.
2. Following cannulation, contrast may be injected to fluoroscopically confirm the position of the device by attaching a syringe to the appropriate injection port of the catheter fitting. **Note:** Remove cap from injection port prior to flushing. This catheter has a radiopaque tip to assist with fluoroscopic visualization.
3. A wire guide may be placed into either wire guide port and into the duct to maintain position and facilitate the exchange of compatible wire guided devices. **Note:** Before each introduction of a wire guide into a device, wet the wire guide in a sterile water or saline bath. A contrast-filled catheter may make the wire guide difficult to advance. Flush the contrast from the catheter with sterile water or saline before each introduction of the wire guide, observing the necessary cautions if the device is positioned in the pancreatic duct.
4. Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

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Wilson-Cook Medical, Inc.
Cook Ireland Ltd.